



M-Tech/Monoclonal Antibodies, Inc.

Product Catalogue – 2022



Legionella monoclonal antibodies, antigens and control slides.

Table of Contents:

	Page
Legionella pneumophila “All in One Conjugate”	3
Legionella Species b - m “All in One Conjugate”	4
Legionella pneumophila monospecific conjugates and antigens	5
Legionella b-m monospecific conjugates and antigens	6
Monospecific conjugates (NOT IN PANVALENT)	7
DFA Control Slides	8-9
Contact Information	10

ALL FITC CONJUGATES are in Evans blue counter stain and contain merthiolate as preservative.

ALL ANTIGEN SUSPENSIONS are heat killed antigens suspended in PBS plus a clear mordant to assist in attachment of cells to slide.

Reference: HEW publication no.[CDC] 79-8375 – “ Legionnaires’ “the disease, the bacterium and methodology. Detection of Legionnaires’ Disease Bacteria in Clinical Specimens by Direct Immunofluorescence

by William B. Cherry & Roger M. McKinney

All in One Conjugates

PRODUCT # 92-148-FL (All in One Conjugate)

DFA test for Legionella pneumophila Sero-Groups 1 – 14

Panvalent FITC Conjugate produced to the following pneumophila strains:

Genus/Species: Legionella pneumophila

Serogroups: 1-14

Strains: Knoxville 1, Bellingham 1, Olda, Philadelphia 1, Togus 1, Bloomington, 2, Portland 1, Los Angeles 1, Cambridge 1, Dallas 1, Chicago 2, Chicago 8, Concord 3, IN-23-G1-C2, Leiden 1, 797-PA-H, 570-CO-H, 82A3105, 1169-MN-H

DFA Mounting Media

Product # 92-057-MF

DFA Negative Control

Product # 92-600-FL

“Why test for only one or maybe six serogroups of Legionella pneumophila when you could test for all fourteen serogroups with one conjugate?”

*Lee
Thacker –*

PRODUCT # 92-149-FL (All in One Conjugate)

DFA test for Legionella Species b - m (bozemanii - maceachernii). Panvalent FITC Conjugate produced to the following Legionella species strains:

<u>Species</u>	<u>Serogroup</u>	<u>Strain</u>
Legionella bozemanii	1	WIGA
Legionella bozemanii	2	Toronto 3
Legionella dumoffii	1	NY-23
Legionella gormanii	1	LS-13
Legionella micdadei	1	TATLOCK
Legionella longbeachae	1	Long Beach 4
Legionella longbeachae	2	Tucker 1
Legionella jordanis	1	BL-540
Legionella oakridgensis	1	Oak Ridge 10
Legionella wadsworthii	1	81-716
Legionella feeleeii	1	WO-44C-C3
Legionella feeleeii	2	691-WI-H
Legionella sainthelensi	1	Mt. St. Helens 4
Legionella anisa	1	WA-316-C3
Legionella santicrucis	1	SC-63-C7
Legionella spiritensis	1	Mt. St. Helens 9
Legionella hackeliae	1	Lansing 2
Legionella hackeliae	2	798-PA-H
Legionella maceachernii	1	PX-1-G2-E2

DFA Negative Control Conjugate

Product # 92-600-FL

DFA Mounting Media

Product # 92-057-MF

LEGIONELLA PNEUMOPHILA MONOSPECIFIC CONJUGATES AND ANTIGENS

* Designated Sero-Group reference type strains

** Suggested supplemental strain to give more complete Sero-Group coverage

<u>Serogroup</u>	<u>Strain</u>		<u>Product Number</u>	
			<u>FITC Conjugate</u>	<u>Antigen Suspension</u>
1	Knoxville 1	**	92-100-FL	92-100-H
1	Bellingham 1	**	92-101-FL	92-101-H
1	Olda	**	92-102-FL	92-102-H
1	Philadelphia 1	*	92-103-FL	92-103-H
2	Togus 1	*	92-104-FL	92-104-H
3	Bloomington 2	*	92-105-FL	92-105-H
4	Portland 1	**	92-106-FL	92-106-H
4	Los Angeles 1	*	92-107-FL	92-107-H
5	Cambridge 1	**	92-108-FL	92-108-H
5	Dallas 1	*	92-109-FL	92-109-H
6	Chicago 2	*	92-110-FL	92-110-H
7	Chicago 8	*	92-111-FL	92-111-H
8	Concord 3	*	92-112-FL	92-112-H
9	IN-23-G1-C2	*	92-113-FL	92-113-H
10	Leiden 1	*	92-114-FL	92-114-H
11	797-PA-H	*	92-115-FL	92-115-H
12	570-CO-H	*	92-116-FL	92-116-H
13	82A3105	*	92-117-FL	92-117-H
14	1169-MN-H	*	92-118-FL	92-118-H

DFA Negative Control Conjugate

92-600-FL

LEGIONELLA b-m MONOSPECIFIC CONJUGATES AND ANTIGENS

* designated Species Type Strains.

<u>Species</u>	<u>Serogroup</u>	<u>Strain</u>		<u>Product Number</u>	
				<u>FITC Conjugate #</u>	<u>Antigen Suspension #</u>
L. bozemanii	1	WIGA	*	92-120-FL	92-120-H
L. bozemanii	2	Toronto 3	*	92-121-FL	92-121-H
L. dumoffii	1	NY-23	*	92-122-FL	92-122-H
L. gormanii	1	LS-13	*	92-123-FL	92-123-H
L. micdadei	1	TATLOCK	*	92-124-FL	92-124-H
L. longbeachae	1	Long Beach 4	*	92-125-FL	92-125-H
L. longbeachae	2	Tucker 1	*	92-126-FL	92-126-H
L. jordanis	1	BL-540	*	92-127-FL	92-127-H
L. oakridgensis	1	Oak Ridge 10	*	92-128-FL	92-128-H
L. wadsworthii	1	81-716	*	92-129-FL	92-129-H
L. feeleii	1	WO-44C-C3	*	92-130-FL	92-130-H
L. feeleii	2	691-WI-H	*	92-131-FL	92-131-H
L. sainthelensi	1	Mt. St. Helens 4	*	92-132-FL	92-132-H
L. anisa	1	WA-316-C3	*	92-133-FL	92-133-H
L. santicrucis	1	SC-63-C7	*	92-134-FL	92-134-H
L. spiritensis	1	Mt. St. Helens 9	*	92-137-FL	92-137-H
L. hackeliae	1	Lansing 2	*	92-138-FL	92-138-H
L. hackeliae	2	798-PA-H	*	92-139-FL	92-139-H
L. maceachemii	1	PX-1-G2-E2	*	92-140-FL	92-140-H
DFA Negative Control Conjugate				92-600-FL	
DFA Mounting Media				92-057-MF	

DFA MONO-SPECIFIC CONJUGATES (NOT IN PANVALENT)*

<u>Species</u>	<u>Strain</u>	<u>Product Number</u>	
		<u>FITC Conjugate #</u>	<u>Antigen Suspension #</u>
L. pneumophila SG 15	Lansing 3	92-119-FL	92-119-H
L. steigerwaltii	SC-18-C9	92-135-FL	92-135-H
L. spiritensis	MSH-9	92-137-FL	92-137-H
L. jamestowniensis	JA-26-G1-E2	92-141-FL	92-141-H
L. cherrii	ORW	92-142-FL	92-142-H

DFA Mounting Media 92-057-MF

DFA Negative Control 92-600-FL

Legionella pneumophila DFA Control Slides

Catalog # 92-701-S

<u>Well</u>	<u>Serogroup</u>	<u>Strain#</u>	<u>Product #</u>
1	1	Philadelphia 1	92-103-H
2	1	Knoxville 1	92-100-H
3	1	Olda	92-102-H
4	2	Togus 1	92-104-H
5	3	Bloomington 2	92-105-H
6	4	Los Angeles 1	92-107-H
7	4	Portland 1	92-106-H
8	5	Dallas 1	92-109-H
9	5	Cambridge 1	92-108-H
10	6	Chicago 2	92-110-H
11	7	Chicago 8	92-111-H
12	8	Concord 3	92-112-H
13	9	IN-23-G1-C2	92-113-H
14	10	Leiden 1	92-114-H
15	11	797-PA-H	92-115-H
16	12	570-CO-H	92-116-H
17	13	82A3105	92-117-H
18	14	1169-MN-H	92-118-H
DFA Negative Control			92-600-FL
DFA Mounting Fluid			92-057-MF

Store at -20° Centigrade - Desiccate if possible. Prior to staining slide should be gently heat-fixed.



Legionella b-m DFA Control Slides

Catalog # 92-731-S

<u>Well</u>	<u>Serogroup</u>	<u>Description</u>	<u>Strain</u>	<u>Product #</u>
1	1	L. bozemanii	WIGA	92-120-H
2	2	L. bozemanii	Toronto 3	92-121-H
3	1	L. dumoffii	NY-23	92-122-H
4	1	L. gormanii	LS-13	92-123-H
5	1	L. micdadei	TATLOCK	92-124-H
6	1	L. longbeachae	Long Beach 4	92-125-H
7	2	L. longbeachae	Tucker 1	92-126-H
8	1	L. jordanis	BI-540	92-127-H
9	1	L. oakridgensis	Oak Ridge 10	92-128-H
10	1	L. wadsworthii	81-716	92-129-H
11	1	L. feeleii	WO-44C-C3	92-130-H
12	2	L. feeleii	691-WI-H	92-131-H
13	1	L. sainthelensi	Mt. St. Helens 4	92-132-H
14	1	L. anisa	WA-316-C3	92-133-H
15	1	L. santicrucis	SC-63-C7	92-134-H
16	1	L. hackeliae	Lansing 2	92-138-H
17	2	L. hackeliae	798-PA-H	92-139-H
18	1	L. maceachernii	PX-1-G2-E2	92-140-H

DFA Negative Control 92-600-FL

DFA Mounting Fluid 92-057-MF

Store at -20° Centigrade - Desiccate if possible. Prior to staining slide should be gently heat-fixed.

About m-Tech

M-Tech was founded in 1984 by Lee Thacker, one of the CDC's principal scientists assigned to research the 1976 Legionella outbreak. As part of his research in detecting legionella from the 1976 outbreak, he helped develop reagents for the first DFA legionella test while describing the first four serogroups of Legionella pneumophila. He also produced the first monoclonal antibodies to Legionella.

M-Tech commercialized the first DFA conjugates used in the detection of Legionella pneumophila while also developing an "all in one" conjugate that can detect up to 14 serogroups in one conjugate. m-Tech has also created a novel "all in one" conjugate for the detection of 19 Legionella Species b-m (bozemanii-maceachernii).

M-Tech is your expert in developing conjugates and antigens specific to Legionella used in research and environmental legionella testing.

m-TECH products are sold throughout the U.S, Europe and other countries around the world.
Serogroup 1-14 and Species b-p DFA tests are accurate and reliable.
Contact us by e-mail at m-tech@mindspring.com for more information on distribution in your area.

Contact Information

If you need more information we can be reached at the following:

Telephone (Toll Free U.S.)
1-888-683-2414

International
01-770-521-1960

Fax: 1-770-521-1930

Email: m-tech@mindspring.com