

INTELLIGENTLY DESIGNED
GUIDED BY
INNOVATION.




REDUCTION AND PLATE ALIGNMENT

5th Metatarsal Hook Plate System

Hook plate guide aids in proper plate alignment, fracture reduction, and homerun screw trajectory.

Nitinol drill pin allows the surgeon to pre-drill for the exact length of the plate's hooks.

Insertor threads onto the plate, allowing it to slide down a guidewire for easy, accurate placement.



PROVISIONAL PLATE PLACEMENT



UNITE[®]
FOOT & ANKLE

For more information or to schedule a case, contact your Medline UNITE representative, or visit medlineunite.com.

© 2019 Medline Industries, Inc. All rights reserved. Medline and Medline UNITE are registered trademarks of Medline Industries, Inc. / MKT19W433779 / 49