

VITRILE GLOVES

BEST PROTECTION

BEST PERFORMANCE

BEST VALUE

VITRILE GLOVES

ESSENTIAL SERIES

Vitrile, powder-free examination gloves with a more robust fit and feel

This non-sterile glove is a standard nitrile + vinyl formulation.

A more robust thickness provides a more substantial fit and feel.

An inner polymer coating makes donning easy and the absence of powder considerably reduces the risk of skin abrasions.

These premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile styles have the peace of mind to focus on what matters most: the health and Protection of themselves and their patients.

WARNING: No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.



Product Code	Size	Length	Inner Carton	Outer Carton
VWNGXS	XS	24 cm	100 pieces	10 Inner Cartons
VWNGS	S	24 cm	100 pieces	10 Inner Cartons
VWNGM	M	24 cm	100 pieces	10 Inner Cartons
VWNGL	L	25 cm	100 pieces	10 Inner Cartons
VWNGXL	XL	25 cm	100 pieces	10 Inner Cartons

Please note:

Approx cbm of each outer carton is 0.04


Approx Gross weight of each outer carton is 6 kgs

Why choose our gloves:

China is the world's largest exporter of medical gloves made from natural rubber and nitrile. These gloves are exported all over the world with USA and EU being top export Destinations. Our Gloves from China conform to standards such as the ASTM, ISO And EN.

VITRILE GLOVES

ESSENTIAL SERIES



- Industry leading gloves offer unrivalled protection, cleanliness and quality
- Nitrile construction results in products that are stronger and leaner than latex Gloves, and feature better protection against a wider range of chemicals
- Gloves are anti-static tested to protect the wearer and equipment.
- Textured fingertips enhance better grip and tactile sensitivity for safer and More efficient processes
- Beaded cuffs add strength to the gloves, reducing the risk of tearing and Increasing their durability.

VITRILE GLOVES

HOW TO USE:

- Dry hands before donning the gloves
- Take out the product after checking package
- Gloves need to be donned carefully, do not apply excessive to pull the gloves when donning
- Keep the fingernails short and change gloves frequently
- The gloves should be discarded with an appropriate waste disposal container after use according to local disposal regulations

CAUTIONS FOR USE:

- The product is for single use only
- Does not contain natural rubber latex

STORAGE CONDITION:

- The atmosphere must always be conducive to the storage of gloves
- Gloves before and after opened for use shall be stored in dark, cool And dry place.
 - Gloves need to avoid area with excessive heat.
- Gloves need to be stored at the areas with good air ventilation.
- Gloves should be shielded from direct sunlight and fluorescent lights.

VITRILE GLOVES

ESSENTIAL SERIES