

# ANTI BACTERIAL RANGE





# HAND SANITISER




**STOCK SELLING FAST!  
CALL US NOW!**



**75% alcohol**



**World Health  
Organization  
recommended  
Alcohol based  
hand rub**





# PRINT OPTIONS



**GENERIC BRANDING**  
Standard option

**CUSTOM STICKERS**  
Can be ordered as an add on



## PCA04 30ml Hand Sanitiser Gel

**Product Colour:** Clear  
**Material: Bottle:** Bottle: PE  
75% Ethyl Alcohol  
**Product Size:** 3\*8\*2 cm



## PCA05 50ml Hand Sanitiser Gel

**Product Colour:** Clear  
**Material: Bottle:** Bottle: PE  
75% Ethyl Alcohol  
**Product Size:** 4.5\*10.2\*2.5cm





## PCA06 60ml Hand Sanitiser Gel

**Product Colour:** Clear  
**Material: Bottle:** Bottle: PE  
75% Ethyl Alcohol  
**Product Size:** 9.5\*4.5\*2.7 cm



## PCA07 100ml Hand Sanitiser Gel

**Product Colour:** Clear  
**Material: Bottle:** Bottle: PE  
75% Ethyl Alcohol  
**Product Size:** 12.5\*4\*4 cm



## PCA08 300ml Hand Sanitiser Gel

**Product Colour:** Clear  
**Material: Bottle:** Bottle: PE  
75% Ethyl Alcohol  
**Product Size:** 16\*6\*6 cm



## PCA09 500ml Hand Sanitiser Gel

**Product Colour:** Clear  
**Material: Bottle:** Bottle: PE  
75% Ethyl Alcohol  
**Product Size:** 15\*8.5\*8.5 cm



## PCA010 1L Hand Sanitiser Gel

**Product Colour:** Clear  
**Material: Bottle:** Bottle: PE  
75% Ethyl Alcohol  
**Product Size:** 25.2\*8.8\*8.8 cm

# PROTECTIVE FACE MASKS

**TGA APPROVAL AWARDED**







# TGA APPROVAL ALLOWS YOU TO SUPPLY MASKS TO HEALTH CARE PROVIDERS, AND OTHERS WHO REQUIRE IT BY LAW.



Australian Government

Department of Health  
Therapeutic Goods Administration

## Australian Register of Therapeutic Goods Certificate

Issued to

**Fair Ocean Trading Australia Pty Ltd**

for approval to supply

**Fair Ocean Trading Australia Pty Ltd - Public N95 respirator**

<b>ARTG Identifier</b>	[REDACTED]
<b>ARTG Start date</b>	8/04/2020
<b>Product Category</b>	Medical Device Included Class 1
<b>GMDN</b>	[REDACTED]
<b>GMDN Term</b>	Public N95 respirator
<b>Intended Purpose</b>	Filtration Face Masks (filter efficiency 95%) designed to be placed over the nose and mouth of a member of the general public to permit normal breathing while protecting the wearer from exposure to pathogenic biological airborne particulates. The product is typically made of multiple layers of non-woven polymers, and incorporates a forming nosepiece (metal wire); it is typically secured using elastic head straps. It is a single-use device.

Manufacturer Details	Address	Certificate number(s)
[REDACTED]	[REDACTED]	[REDACTED]

### ARTG Standard Conditions

The above Medical Device Included Class 1 has been entered on the Register subject to the following conditions:

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

### Products Covered by This Entry

#### 1. Public N95 respirator

### Product Specific Conditions

No specific conditions have been recorded against this entry.

Therapeutic Goods Administration  
PO Box 100, Woden ACT 2606 Australia  
Phone: 1800 020 653  
Email: info@tga.gov.au

ARTG Identifier: [REDACTED]  
ARTG Start Date: 8/04/2020

# MARKET SITUATION

Demand has exploded around the world,  
but poor quality is causing lots of issues

 South China Morning Post

## Netherlands recalls 600,000 face masks from China due to low quality

- Top-of-the-line N95 masks did not fit well and filters did not function properly, Dutch government says
- It comes after Spain asked for Covid-19 rapid test kits imported from China to be replaced because they were unreliable



 Stuart Lau

Published: 3:45pm, 29 Mar, 2020 ▾





# RISK MANAGEMENT

To ensure the masks you are sourcing meet the needs of the clients, you need to understand the below key points:



## A. KNOW THE RELEVANT STANDARD

Different standards of masks are suitable for different groups.



TYPE OF MASK	STANDARDS	EXPLANATION	COMMENTS	PRODUCTS MEETING STANDARD
GENERAL USE	EN 149-2001 FFP1	EUROPEAN STANDARD FOR FILTERING FACE MASKS		
OFTEN USED IN HEALTH CARE	AS 4381: 2015	AUSTRALIAN STANDARD	EQUIVALENT STANDARDS	PCA021, PCA022
	CE: EN14683 Type 1 & 2	EUROPEAN STANDARD		
HIGHER GRADE OF FILTRATION	N95 - NIOSH - 42CFR84	US STANDARD	EQUIVALENT STANDARDS	PCA023, PCA024
	KN95 - GB2626 - 2006	CHINA STANDARD		
	FFP2 EN 149-2001	EUROPEAN STANDARD		
	P2 - AS/NZA 1716:2012	AUSTRALIAN/NZ STANDARD		

## B. ENSURE TGA CERTIFICATION

For use in health care, not only are medical grade face masks required, but they will need to come with the relevant TGA approvals for their factory source.

## C. CHOOSE A REPUTABLE SUPPLY PARTNER

With the amount of poor quality product flooding the market currently, you need a partner who will ensure you are getting what you paid for, not just selling you anything they can get their hands on.



# WHY CHOOSE US

When you source your protective face masks from us we give you the extra assurance by undertaking:



01

## FACTORY DUE DILIGENCE

I. RESEARCH OF FACTORY INFORMATION

II. INSPECTION OF CERTIFICATION



02

## PHYSICAL VISITS

I. INSPECTION OF FACTORY FACILITIES

II. DETERMINING CAPACITY TO SATISFY ORDERS



03

## SAMPLING OF PRODUCTION

UNDERTAKING SAMPLE BATCHES TO ENSURE THE PRODUCT MEETS THE REQUIRED STANDARD

## ISO13485 MEDICAL SUPPLY MANUFACTURER

## EN14683 LAB TESTED & CERTIFIED

**intertek**  
Intertek Testing Services

**CERTIFICATE OF REGISTRATION**

This is to certify that the management systems of:

**MEDICAL CO., LTD.**

has been registered by Intertek as conforming to the requirements of **ISO 13485:2016**

The management system is applicable to:

Manufacture and distribution of disposable sterile medical equipment covers, N95 respirator, gowns, shoe collector bag, apron, bed protection net, disposable plastic apron, gown, cap, face mask, management system, connectors, valves, coils, roller clamp, handclasp, inflators, blood smear applicator, cap with separator, syringe, blood transfer stand, urine transfer stand, urine tube and cap, negative suction system.

Intertek Testing Services, Ltd.  
1000 Eastman Avenue, Suite 101, Irvine, CA 92614, USA

**PONY 潘尼测试**  
Pony Testing International Group

**Test Report**

No. 202009170001

Sample Description: Disposable face mask for medicine

Applicant: [Redacted]

Pony Testing International Group

**PONY 潘尼测试**  
Pony Testing International Group

**Test Results**

No. 202009170001 Page 1 of 2

Sample Description	Disposable face mask for medical use	Sample Specification	Trade Mark
Applicant	[Redacted]		
Received Date	2020-09-17	Manufacturing Date of Lot No.	2020.5.9
Test Date	2020-09-17-2020-09-27	Sample Grade	
Sample Quantity	60Pcs	Test Type	Comprehending Test
Sample Status	Solid	Test Environment	To meet the requirements
Test Items	See the next page		
Reference Methods	See the next page		
Main Instruments	Bacterial Strain efficiency detector, Constant temperature incubator etc.		
Note	Limit use: EN 14683 - 2019 Sample source: self-made		
Edited by	Lu Meng		
Checked by	Sijie Shen		
Approved by	Dr. Bao		
Issued Date	2020-09-27		

**PONY 潘尼测试**  
Pony Testing International Group

**Test Results**

No. 202009170001 Page 2 of 2

Sample Description and Number	Test Items	EN136		Test Results	Reference Methods
		Typical	Typical		
R2314570A Disposable face mask for medical use	Bacterial Strain efficiency (CFU)	≤10	≤10	≤10	EN 14683 - 2019
	Differential pressure (Pa)	< -40	< -40	51	EN 14683 - 2019
	Speech resistance pressure (Pa)	≥ 1	≥ 16.8	≥ 16.8	EN 14683 - 2019
	Material characteristics (kPa)	≤ 2.0	≤ 2.0	4	EN 14683 - 2019

Sample Number and Photo: [Image of a face mask]







## OUR HIGH STANDARD MASKS RANGE



**1.2 + Million**

*reliable masks sold  
in the past week*

# KN95 Masks

( KN95-GB2626-2006 )  
( P2-AS/NZA 1716:2012 )

## Filter Performance $\geq$ 95%



## PCA024

KN95 Face Mask



Adjustable nose wire

**4** layers of protection



Therapeutic Goods Administration of Australia  
**TGA**  
CERTIFIED



## PCA024 KN95 Face Mask

**Product Colour:** White

**Product Size:** 155 x 105 x 10mm

**Material:** Outer layer: PP 60g/m<sup>2</sup>  
Second layer: Filter meltblown fabric 25g/m<sup>2</sup>

Third layer: Es hot air cotton 35g/m<sup>2</sup>

Inner layer: Hygroscopic PP fabric 35g/m<sup>2</sup>

Elastic: Polyester & Nylon

Average filtration efficiency: 95%+

**Certifications:** ISO 13485:2003, TGA Approval



# Personal Disposable Masks

## High Standard Masks

Medical Grade Masks (CE:EN14683 TYPE I,II+) (AS 4381:2015)



**COMFORTABLE**



**BREATHABLE**



### PCA022

#### Personal Disposable Mask - 50PC Box

**Product Colour:** Natural Blue

**Product Size:** Mask: 17.7 x 9.4 x 0.2 cm

Package: 19 x 10 x 8 cm

**Material:** Outer layer: PP non-woven fabric

Mid layer: Filter meltblown fabric

Inner layer: Hygroscopic PP non-woven fabric

Bridge of nose: Plastic

Elastic: Polyester & Nylon

Filtration efficiency: 95%+

**MOQ:**100

**Certifications:** Standard Grade: CE:EN149

High Grade: EN 14683/ISO 13485:2016;

TGA Approval



### PCA021

#### Personal Disposable Mask - 10PC Pack

**Product Colour:** Natural Blue

**Product Size:** Mask: 17.7 x 9.4 x 0.2 cm

Package: 24.8 x 12.5

Digital label: 6 x 10 cm

**Material:** Outer layer: PP non-woven fabric

Mid layer: Filter meltblown fabric

Inner layer: Hygroscopic PP non-woven fabric

Bridge of nose: Plastic

Elastic: Polyester & Nylon

Filtration efficiency: 95%+

**MOQ:**1000

**Certifications:** EN 14683/ISO 13485:2016

# HOW TO CORRECTLY USE A FACE MASK ?





# 75% ALCOHOL WET WIPES



**75% ALCOHOL**  
**KILLS 99% OF GERMS**

# TYPE A : 10PC SMALL PACK

**48HR**  
EXPRESS



10 Wipes Per Pack



Compact size



# TYPE B : 50PC LARGE PACK

**NEW**

**GREAT NEW 50PC LARGE  
PACK AVAILABLE!**








**75% Alcohol**



**Easy to apply plastic lid**



**Clean hands and surfaces**



# PRINT OPTIONS

## 10pc pack



Unbranded



Digital label

## 50pc pack



Unbranded



Digital label



## PCA01 75% Alcohol Wet Wipes

**Product Colour:** Blue

**Material:** Pouch: PE, Wipes: Cotton

**Product Size:**

Pack: 150(L) x 80(W)mm

Wipes: 180(L) x 140(W)mm

Digital label size: 2.5 x 5 cm

**Certifications:** CE, FDA, MSDS



## PCA02 75% Alcohol Wet Wipes - 50PC Pack

**Product Colour:** White

**Material:** Pouch: PE, Wipes: Cotton

**Product Size:**

Pack: 200(L) x 100(W) x 45(D)mm

Wipes: 150(L) x 200(W)mm

Digital label size: 4.7 x 8 cm

**Certifications:** CE, ROHS, MSDS



# MULTI-MODE THERMOMETER





# BUYERS BEWARE - TGA APPROVAL REQUIRED

Thermometers, under the Therapeutic Goods Administration's guidelines, are considered Class IIa medical devices. Sale in Australia of all medical devices requires TGA approval.

**DON'T TAKE THE LEGAL AND FINANCIAL RISK, AND BUY NON TGA APPROVED PRODUCTS. ALWAYS ASK FOR A TGA CERTIFICATE BEFORE BUYING.**

## TGA CERTIFICATION



Australian Government  
Department of Health  
Therapeutic Goods Administration

### Australian Register of Therapeutic Goods Certificate

Issued to

**Fair Ocean Trading Australia Pty Ltd**

for approval to supply

**Fair Ocean Trading Australia Pty Ltd - Thermometer, infrared, skin**

**ARTG Identifier** [REDACTED]

**ARTG Start date** 24/04/2020

**Product Category** Medical Device Included Class IIa

**GMDN** 17888

**GMDN Term** Thermometer, infrared, skin

**Intended Purpose** A hand-held, battery powered device that is intended to be used for the monitoring of body temperature. It measures the infrared energy emitted from the skin surface of the forehead so there is no need for contact with the skin.

Manufacturer Details	Address	Certificate number(s)
[REDACTED]	[REDACTED]	DV-2020-MC-12179-1

#### ARTG Standard Conditions

The above Medical Device Included Class IIa has been entered on the Register subject to the following conditions:

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

#### Products Covered by This Entry

1. Thermometer, infrared, skin

#### Product Specific Conditions

No specific conditions have been recorded against this entry.

Therapeutic Goods Administration  
PO Box 100, Woden ACT 2606 Australia  
Phone: 1800 020 653  
Email: info@tga.gov.au

ARTG Identifier: [REDACTED]  
ARTG Start Date: 24/04/2020



GET RESULTS IN **1 SECOND** WITH  $\pm 0.2^{\circ}\text{C}$  !

**MORE ACCURATE AND RELIABLE RESULTS !**

- Advanced sensor
- Save up to 20 temperature readings



**MULTI-MODE THERMOMETER**



Designed for all ages



Can take room / object temperature

**LED TEMPERATURE INDICATOR**



Red for fever

Green for regular temperature



# PRINT OPTIONS



UNBRANDED OPTION



SCREEN PRINTING

## USE THEM ANYWHERE



### PCA015

### Touchless Digital Infrared Forehead Thermometer

**Product Colour:**White

**Material:** ABS & Metal

**Certifications:** CE, FCC, TGA

**Product size:** 6 x 16 x 3cm

**Product weight:** 118g



# SAFETY PROTECTIVE GOGGLES



# High level protection

Light weight

36g

28g


Can be worn over glasses  
- PCA011

Adjustable for the perfect fit  
- PCA012



# PRINT OPTIONS

PCA011




PCA012



Unbranded

PCA011



PCA012



Screen Print

# USE THEM ANYWHERE



## PCA012 Anti Fog Retractable Goggles

**Product Colour:** Translucent and black

**Material:** PVC, PC

**Product size:** 15 x 4.5 x 6cm

**Product weight:** 28g



## PCA011 Anti Fog Safety Protective Goggles

**Product Colour:** Translucent

**Material:** PC

**Product size:** 15.5 x 5 x 4.5cm

**Product weight:** 36g





