



## Dental Unit G7 model



### Seat Size

Large Width – 54 cm  
Seat/64cm Backrest

### Package

150 \* 103 \* 116cm  
G.W. 255 KG

### Container Volume

16units/20GP  
32units/40GP

### Fiber leather

### Major Parts

1. Water heating system
2. Spittoon flushing and cup filling
3. Quiet electric lift chair (supported by (Denmark motors)
4. Fashionable control panel
5. Chair movement system with 9 memory positions
6. LED film viewer
7. Assistant control panel
8. Rotatable ceramic spittoon
9. Shadowless LED sensor dental light
10. Multifunction pedal switch
11. Pure water bottles, 600ml
12. Dentist stool with backrest
13. Imported polyether inner tubing

### Technical Specification

Supply voltage: AC110V/220V 60Hz/50Hz

Motor voltage: DC24V

Air pressure: 0.5MPa - 0.8MPa

Water pressure: 0.2MPa - 0.4MPa











Imported motor  
9 memory system  
Extra stainless steel tray  
Injected Plastic  
6LED Dual Color sensor lamp  
Laser printed touch button assistant tray

