



Portable Anaesthesia Machine LM-PH



Portable Anesthesia Machine Technical Specification

Main Frame

Overall Dimension	50*30*60CM(L*W*H)
Working Mode	Closed,Semi-Closed,Semi-Open System
Application	1~100KG Patients
Flow Meters	O2 Two Tubes(0~10L/Min)
Vaporizer Quantity	Isoflurane,Sevoflurane,Halothane(Optional One Of Three)
Tidal volume (Vt)	0 mL ~ 1500 mL
Rapid oxygen supply	35 L/min ~ 75 L/min
APL Valve	0.5~6.0KPa
CO2 Canister	1000mL
The Pressure Gauges	Two for O2,One for Airway pressure
Gas Supplies	O2 (Optional N2O)
I:E:	As the doctor operation
BPM	As the doctor operation
Monitoring	Airway pressure, respiratory mode

Ventilation Mode

Manual (Optional Portable Ventilator: A/C,IPPV)

Vaporizer

Concentration Range	0.2~5%, 0.2~8%
Flow Range	0.2 to 5L/Min
Automatic Compensation	Flow,Temperature and Pressure
Lock	Self-Lock and Interlock
Volume	250mL
Mount	Pour-fill Selectatec Mount

Working Condition

Electricity No Need, only Manual Operation

Optional

Portable Ventilator: Goods size: L20*W20*H34CM, Net weight : 5KG (Packing Size:40*27*46 cm 7KG)


Trolley: Goods size L60*W45*H120CM,Net Weight : 16KG (Packing Size: 75*47*24cm 18KG)

Aluminum Pull-rod Case: 46*34*90 cm,32KG

Packing Size

Anesthesia machine Carton Packing size	53*57*44CM(L*W*H)
G.W.	25KGS
N.W.	14KGS
CBM	0.2 m3

Remark: Above configuraions include standard and option. Please contact with sales representative for more details.



Add: Room 1901,A3 Kechuangcheng Building, No.30 Fengzhan Road,
Yuhuatai Distinct,Nanjing, Jiangsu,China
Tel: +86-25-84452678 Mobile: +86-13951660247
E-mail: Sales@healicom.com www.healicom.com