# POSTURAL AND BEHAVIOURAL HARNESSES CAREVA HARNESS SYSTEM




#### **CAREVA COMBINATION HARNESS**











# POSTURAL AND BEHAVIOURAL HARNESSES CAREVA HARNESS SYSTEM

### **CAREVA COMBINATION HARNESS**



The **Careva Harness System** is a postural and support harness that provides upper body support and prevents the person falling forward or to the side to ensure that the vehicles seatbelt works correctly.

#### What is the function of Careva Combination Harness?

This product is a full-body support, that is used to ensure the occupant is safely seated in the correct seating position for the vehicles seatbelt.



The system consists of a range of single belts that can be combined as required for users of all ages who need additional support for an upright sitting position. **Careva Combination Harness** can be used in any seat, provided there is an opening between the seat and the backrest.



Universal postural belt system that can be used in any seat, provided there is an opening between the seat and the backrest

Ouick and easy to install and transfer between vehicles

## SPECIFICATIONS

Available in: Small 4-7 years Medium 8-12 years Large 12 years - Adult

Can be used with wheelchairs with a stable backrest

Careva Harness System DOES NOT replace the vehicle safety belt


The **Careva Combination Harness** is a positioning belt designed with the help of many involved parents of children with disabilities, drivers in the community transport, Occupational Therapists and health professionals to ensure the correct postural supports are achieved.

ADABT-A-CAP

MOBILITY SOLUTIONS

Our harnesses are quick to install and can be easily transferred from across vehicles.

## BENEFITS



No modifications to your vehicle


Provides you with transport independence

# **OPTIONAL ACCESSORIES**



# THORACIC BELT

An additional option that provides additional upper body chest support



**Autoadapt** 

**REQUIREMENTS** It is a requirement to keep a *current* Doctor's Certificate in your vehicle when this Harness is in use in the vehicle

11 O'Briens Hill Road O'Connell NSW 2795

0414 736 854



www.mobilityengineering.com.au

adaptacar601@gmail.com

