

SWIVEL SEAT SYSTEMS – THE TURNY FAMILY

TURNY EVO

The Turny Evo is a programmable lowering swivel system for transferring into the vehicle.

The Turny Evo brings the seat completely out of the vehicle giving the occupant ample transfer space to get seated. Once seated the Turny Evo will lift you to a safe position inside the vehicle. It is recommended to be used with one of our low profile seats so that you can optimise the head and knee room when inside the vehicle.



It is a fully automatic, programmable swivel seat where all movements are operated by the hand held remote. This unique programmable function means that you can set the parameters of the swivelling function to optimise the amount of room generated in the vehicle and also allows it to be installed into the drivers seating position of many vehicles.

Key features of the Turny Evo include:

- Features 390mm of vertical travel
- Programmable movement path
- Full backrest adjustability
- Powered forwards and backwards movement once inside the vehicle
- Left or right side applications
- Easy to use back-up and emergency operations
- Inbuilt footrest for easier transfer
- Can be combined with the tilda to increase knee and head room
- Can be combined with carony to eliminate any lifting during the transfer
- Crash tested and CE marked



GETTING SEATED

TURNY ORBIT

Lowering swivel seat system for transferring into the vehicle.

The Turny allows the car seat to swivel out of the door opening and lower down to a height that is suitable for transfer from a wheelchair. It is recommended to be used with one of our low profile seats so that you can optimise the head and knee room when inside the vehicle.



The Turny Orbit is a fully automatic swivel seat where all movements are operated by the hand held remote.

Key features of the Turny include:

- Available in either 300mm, 395mm or 480mm vertical travel versions
- Left or Right side applications
- Easy to use back-up and emergency operations
- Can be combined with Carony to eliminate any lifting during the transfer
- Optional footrest available for easier transfer
- Crash tested and CE marked



TURNOUT SWIVEL SEAT

Seat base that allows you to swivel out of the door for easy access.

The Turnout is a swivel base that enables you to swivel the seat out of the vehicle so that you are facing outside making it easier to transfer in and out. It is recommended to be used with one of our low profile seats so that you can optimise the head and knee room when inside the vehicle. Options for the Turnout include the Caroslide and Tilda which can be used to further increase the swivelling room available in the vehicle.

Both a manual or electric version is available and Turnout can also be used in combination with the Carony transfer wheelchair system.

The key features of the Turnout include:

- Manual or electric versions available
- 2 or 4 door models to optimise swivel space
- Left or Right side applications
- Can be combined with the caroslide to increase leg room
- Can be combined with the Tilda to increase knee and head room
- Can be combined with Carony to eliminate any lifting during the transfer
- Crash tested and CE marked



GETTING SEATED

TURNOUT ACCESSORIES

TILDA

The Tilda enables the seat to be tilted to 15 degrees using 3 motions. There is an easy to use handle on the side of the unit that allows the occupant or carer to tilt the seat. When used in combination with the Carony wheelchair transfer system and the Turny Evo or Turnout this product improves the head space when entering the vehicle.

The first position of the Tilda has been crash tested allowing you to travel in a safe and comfortable position when in the vehicle and can help improve head stability for those with limited neck strength.

There is also an optional footrest that when used with a swivel base allows you to get your feet over the door frame easily and inside the vehicle.



CAROSLIDE

The Caroslide is a powered seat base that can be used to give you additional leg room when used in conjunction with the Turnout for swivelling out of the vehicle. In addition, if the occupant does not have the strength to move the seat forwards and backwards manually the Caroslide can overcome this obstacle. It can also be used with an aftermarket seat such as the Bev, Compact or Recaro seat to convert your seat to a powered seat.



Features include:

- 240mm stroke
- Operated with a remote
- Left or Right side applications
- Can be used with the Turnout for increased leg room
- Crash tested and CE marked



6 WAY BASE INTERNAL ACCESS SYSTEM (SELF DRIVE)

Swivel seat designed for internal transfers from your wheelchair

The 6-way base is a swivel seat system that enables an easy transfer from the wheelchair into a safe and secure front driver passenger position in the vehicle. Suited for vehicles with a flat floor entry into the front seating positions the 6-way base allows the driver to move forwards, backwards, rotate to 140 degrees as well as moving up and down to match the height of your wheelchair.

The 6-way base has also been crash tested giving you peace of mind while travelling in the vehicle.

The key features of the 6-way base include:

- Forward, backwards, up and down movement
- 140-degree rotation
- Easy to use toggle switches
- Ability to install the vehicle's original seat onto the base
- Complete stability inside the vehicle
- Footrest available as an option
- Crash test approved in all driving positions
- Available for both the driver and passenger seating positions



TRANSFER PLATFORMS

TIP-UP PLATE

Vehicle mounted transfer board for safe vehicle transfers.

The Tip-up plate is an aid used by people who find it difficult to transfer to the front passenger or drivers seating positions. It comes in



both a manual and electric version that incorporates a lifting function. As the transfer board is fixed onto the vehicle this reduces the safety risks involved with transferring and has a weight capacity of 120kg.

The fixed version has a transfer plate that can be folded up and out of the way or removed and stored in the vehicle. It is recommended for vehicles where the seat is lower to the ground and the user has difficulty jumping the gap between the wheelchair and the seat.

The electric version incorporates the same function as the manual however also has a lifting mechanism that allows you to transfer onto the plate and then either lift or drop to the height of the seat in the vehicle.

Key features of the Tip-up plate include:

- Manual or electric versions available
- Plate can either be removed or folded out of the way
- Electric version includes an easy to use remote handset with the option to have it wireless
- Low profile design allows this to be installed into most vehicles
- All versions have up to 120kg weight capacity
- This product has been EMC test and is CE marked



GETTING SEATED

STOW AWAY

The Stow-away is a powered transfer seat that helps you get in and out of the car seat with ease. Primarily designed for Ute's and light trucks, this powered base lifts you from the height of your car seat so that you can slide into the car safely.

When not in use the Stow-away is tucked underneath the vehicle in a sealed durable drawer, keeping the space inside of the vehicle untouched.

Key features of the Stow-away include:

- 508mm of vertical travel
- Seat pad size of 25cm (width) x 38cm (length)
- Left or Right side applications
- Easy to use back-up and emergency operations
- Safety handle for increased stability while transferring
- Optional taller seat pad for extra height while transferring
- Optional padded transfer pad to bridge the gap between the Stow-away and car seat
- 136kg weight capacity



TRANSFER SLING

KIVI PERSON LIFT



The Person lift is designed to carry someone from a wheelchair into the seat of the vehicle. By utilising a sling for the transfer the lift eliminates the strain put on the carer for transfers while safely transferring the occupant into the seat. Once you have completed the transfer the Person lift can then be stored securely in the boot of the vehicle.

It uses a pivoting arm to optimise the transfer space with easy to use attachment points for the transfer sling. For larger vehicles there is the Combi person lift that incorporates a lifting cylinder allowing the occupant to reach higher seats.

There is also an easy to use mobile vehicle bracket that allows you to use the Person lift without ever modifying the vehicle and allows the product to be easily transferred between vehicles.

Key features of the Person lift include:

- Fixed or mobile vehicle brackets
- Left or Right side applications
- 100kg weight capacity
- Combi option available for higher vehicles
- Sling sizes: Small, Medium and Large
- Easy to use push-button handset
- Can also be used inside the home either with a mobile frame or fixed wall bracket
- EMC tested and CE marked

