



HOISTING AND STOWING VSL-4400 JOEY LIFT



**MOBILITY
ENGINEERING**

VSL-4400 JOEY LIFT





VSL-4400 Joey Lift makes transporting a scooter or powerchair easy!



The **VSL-4400 Joey Lift** allows you to ride a mobility device onto the platform from either direction, press a button and the Joey lifts and stows your device into the back of your vehicle. The Joey lift is a safe, simple solution to transport your mobility device around during your travels.

Can be applied to a wide range of vehicles including, minivans, full-size vans as well as selected SUV's and station wagons.

FEATURES

-  Fully-powered, one button operation
-  Compact Design
-  Self-levelling platform
-  Platform illumination for night loading

USER REQUIREMENTS

The user or carer need to be able to do the following:



Walk from the rear of the vehicle to the driver's entry and vice versa

Stand for 2-3 minutes




Drive mobility device on/off platform

Hold down a button on hand-held control

Connect a hook to receiver and tighten strap to secure mobility device *(not required for every application)*



BENEFITS

-  Provides daily independence
-  Safe and reliable design
-  Complete solution to stow your mobility device inside the vehicle

SPECIFICATIONS

Maximum Mobility Device Wheelbase: 102cm

Maximum Mobility Device Width: 71cm

Vehicle Types: Minivans, full-size vans, select SUVs and station wagons

Weight capacity: 159kg

