

Cleönmed

-----Driving Technology to care-----



Bubble CPAP / HFNC / nCPAP

BC 70



About Us

Cleonmed was formed by a group of Biomedical Engineers in the year 2021 with an clear objective to provide innovative and cutting edge technologies assisting health care professionals in giving better patient care. With our team's in-depth knowledge over Respiratory care and Neonatal care products we are able to provide optimum solutions instead of just selling products.

We value our customer relationship and hence have best after-sale service team.




OUR NICU PRODUCTS

We provide a range of respiratory and Neonatal care products like:

Bubble CPAP
Led Phototherapy
NICU Ventilator
High Frequency Ventilator

OUR OTHER PRODUCTS

Transport Ventilator
ICU Ventilator
Portable Ventilator
OT Light
Examination Light



BC- 70

BC-70 is a high performance flow generator which can be used as both HFNC and Turbine based Bubble CPAP unit. Its ultra quiet turbine helps in maintaining low noise environment with the critical care units like ICU and NICU. Its flow range of upto 80 lpm helps it to be used for wide range of applications. It has a user friendly touchscreen display of 7.0 inch. It is compatible with any type of Humidifiers.

Technical Specifications:

Therapy Modes	Bcpap	HFNC	nCPAP
Screen Size	7.0"	7.0"	7.0"
Touch Screen	Yes	Yes	Yes
Battery Back-up	2 hrs	2 hrs	2 hrs
Flow Range	0-40 lpm	0-80 lpm	- -
Pressure Range	- -	- -	0-30 cmH2O



Company’s Vision

Cleonmed aims to keep developing unique and high performance medical equipments for Healthcare Professionals easing out their efforts in saving a life.

We understand the criticality of the situation within an ICU, NICU and OT. So we keep our equipments very easy to operate. We also use high quality components in the development of all our products so that we can delivery very accurate and precise results.

We also believe Service is the backbone of sales. Thus ensuring best After sales service support.

<i>Call Us today!</i>		Hours Of Operation	
		<i>Day</i>	<i>Hours</i>
9692819717/7008940482		Monday	9:00 am
<i>Cleonmed</i>		To	To
Peyton Sahi, Buxi Bazar Cuttack-753001, Odisha, India		Saturday	7:00 pm