

Cleonmed

-----Driving Technology to care-----



CP 20 Phototherapy

CP-20 Led Phototherapy is a specially designed phototherapy system to treat newborn jaundice by lowering the bilirubin levels in baby's blood through a process called photo-oxidation.

Calista Specifications:

No. of Blue Leds: 20

No. of white Leds: 4

Wavelength of blue Leds: 460-465nm

Irradiance: 30 – 45 $\mu\text{W}/\text{cm}^2\text{nm}$ at 35 cm

Power Supply: 12V, 5A

Led life: 45000 hrs



User friendly screen layout
With Timer function



Web: www.cleonmed.com

email: info@cleonmed.com