

CLEONMED

...Driving Technology to Care...

The model of basic health care



PM 12 / PM 10
Patient Monitor

PM10 GOLD Masimo SET

Patient Monitor

Smart, Light, Easy operating

Suitable for bedside, ICU, OR



Product Feature

- Standard parameters: ECG, NIBP, RESP, PR, SpO2, TEMP
- Colorful and Clear 12.1" & 10.4" TFT screen, Backlight buttons
- Multiple display modes optional: Standard interface, Big font, ECG standard full display, OXY, Trend table, BP trend, Viewbed
- Ambulatory blood pressure technology, anti-movement
- Special design against high frequency surgical unit, and defibrillation protection
- Support Masimo SpO2-authorization partner
- 13 kinds of arrhythmia analysis
- 15 kinds of drug dose calculation
- Various languages operating systems
- Built-in detachable rechargeable lithium battery 4 hours battery life
- To connect CMS, other bed observation and software updating with wireless and wired mode

Data Storage – stable and fast

- 8000 groups NIBP measurements
- 680 hours Trend data and Trend graphs
- 200 groups Alarm events reviewing
- 2 hours Wave forms reviewing
- Supporting for Arrhythmia analysis and drug dose calculation

Alarming – more safe and reliable

- 3 level Audible and visual alarming



Masimo SET. PM 12 GOLD Patient Monitor

Suitable for Adult, Pediatric and Neonatal.

Reliable quality gives confidence, stable system makes perfect.
Powerful storage capacity ensures data safety.

Movable handle design

Lightweight portable,
practical all-around,
convenient operation

CleomedMedical regards quality as the core, relies on creative
R&D, produces human friendly products, and doctors' reliable
partner.



NIBP list interface



Standard interface



Viewbed interface



Big font interface



Tendency diagram



OxyCRG interface



Support

HD, VGA, USB,
Internet connector output



Double invasive blood pressure (IBP)

Adopting double-channel measurement, and
realizing continuous testing blood pressure
variation of severe patients



KNOB



ONE-BUTTON NIBP



POWER

backlight button design

Optional: EtCO2, Dual-IBP, 12-leads ECG, Touch screen, Built-in thermal printer

PM 12 GOLD
PM 10 GOLD
Patient Monitor



CLEONMED
...Driving Technology to Care...

Peyton Sahi, Buxibazar, Cuttack- 753001
Odisha, India. Tel: +91 98611 41801
Email: info@cleonmed.com



Please log in for more details www.cleonmed.com

Customer service **+91 7847 857459**

- Cleonmed. All Rights Reserved. The product specifications are subject to change without further notice.
- See instructions for contraindication or precautions