



# Clinical 3DEEP<sup>®</sup> Results

Before & After Treatment Photos  
With Long Term Follow-Up



**EndyMed**<sup>™</sup>  
THE 3DEEP<sup>®</sup> COMPANY

The sleekly designed EndyMed PRO system is multi-application, upgradeable and available with multiple handpieces, customized to different applications and treatment areas. EndyMed PRO utilizes proprietary 3DEEP® RF technology and has an easy-to-use touch screen interface.



3DEEP® Large Handpiece



3DEEP® Small Handpiece



3DEEP® iFine Handpiece




EndyMed's proprietary 3DEEP® technology is a multi-source, phase-controlled radiofrequency energy source that is state-of-the-art, safe and effective.

3DEEP uses an array of multiple RF sources and a sophisticated algorithm to manage the phase of the energy flowing between sets of electrodes. The repelling forces between adjacent electromagnetic fields drive the energy vertically into the target tissue, reducing significantly the energy flowing along the skin's surface (and any potential side effects) and eliminating the need for cooling.

## Energy penetration depth of 3DEEP RF

3DEEP RF – focused and contained deep energy flow with no need for cooling. The electrode bars are shown in red and green to denote positive and negative charges, which alternate and push the energy deep into the tissue.



**EndyMed Medical** is a cutting-edge medical device company that develops energy based aesthetic treatment systems for the professional and home use markets. EndyMed's products are powered by our proprietary, FDA cleared 3DEEP® multisource radiofrequency technology, the most advanced RF treatment technology available today for high efficacy, painless wrinkle reduction and skin tightening.

**The EndyMed PRO™** is the most comprehensive high end aesthetic solution platform available today. The PRO offers impressive results with minimal discomfort, no side effects and minimal downtime after treatments. The sleekly designed PRO is a multi-application system that is available with a range of handpieces that are customized for different applications and treatment areas.

**The PURE™** by EndyMed offers professional grade superior skin treatments with minimal downtime in all in one easy to use, compact affordable platform. The PURE provides customized treatment parameters, an extensive range of handpieces, and the highest level of safety features available for safe, effective treatments with minimal discomfort.

**The Newa®**, EndyMed's new breakthrough energy based home use device, provides noninvasive treatment of mild to moderate facial wrinkles. Using the Newa triggers your skin to produce new collagen, lifting your face and neck and leaving your skin smooth, firm and refreshed.



Before treatment



After 6 treatments



9 months follow up

Photographs courtesy of  
I. Vider, M.D., Medical OR, Hertzlia, Israel





Before treatment



After 6 treatments

Photographs courtesy of  
Sharon Kshetry M.D., F.A.C.O.G., MN, USA



Before treatment



After 6 treatments

Photographs courtesy of  
Dr. Mathew Jafarzadeh, Instant Laser Clinic, Melbourne, Australia





Before treatment



After 8 treatments

Photographs courtesy of  
Fiona Wright, M.D., Plano, Texas, USA



Before treatment



After 6 treatments

Photographs courtesy of  
Dr. Mathew Jafarzadeh, Instant Laser Clinic, Melbourne, Australia





Before treatment



After 4 treatments

Photographs courtesy of

Martine Baspeyras M.D., Board Certified Dermatologist, Bordeaux, France



Before treatment




3 Months after 12 treatments

Photographs courtesy of

Isabelle Rousseaux, M.D., Board Certified Dermatologist, Loos, France



3DEEP Facial and Neck



Before treatment



6 months follow up after 6 treatments

Photographs courtesy of  
I. Vider, M.D., Medical OR, Hertzlia, Israel

3DEEP Body



Before treatment



8 months follow up after 6 treatments

Photographs courtesy of  
I. Vider, M.D., Medical OR, Hertzlia, Israel



3DEEP Body



Before treatment




5 months follow up after 6 treatments

Photographs courtesy of  
I. Vider, M.D., Medical OR, Hertzlia, Israel

3DEEP Body



Before treatment



8 months follow up

Photographs courtesy of  
Klaus Fritz, M.D., Landau, Germany





3DEEP Facial



Before treatment




2 years after 8 treatments

Photographs courtesy of  
Isabelle Rousseaux, M.D., Board Certified Dermatologist, Loos, France

3DEEP Facial



Before treatment



After 1 treatment

Photographs courtesy of  
Isabelle Rousseaux, M.D., Board Certified Dermatologist, Loos, France



3DEEP Facial



Before treatment



Immediately after 1 treatment

Photographs courtesy of  
Flor Mayoral, M.D., Miami, Florida

3DEEP Facial and Neck



Before treatment




After 6 treatments

Photographs courtesy of  
Prof. Daniel Cassuto, Milan, Italy



3DEEP Facial and Neck



Before treatment



After 6 treatments

Photographs courtesy of  
Prof. Daniel Cassuto, Milan, Italy

3DEEP Facial and Neck



Before treatment



3 months follow up after 6 treatments



5 months follow up

Photographs courtesy of  
I. Vider, M.D., Medical OR, Hertzlia, Israel



3DEEP Facial and Neck



Before treatment



6 months follow up after 6 treatments

Photographs courtesy of  
I. Vider, M.D., Medical OR, Hertzlia, Israel

3DEEP Facial and Neck



Before treatment



6 months follow up after 6 treatments

Photographs courtesy of  
I. Vider, M.D., Medical OR, Hertzlia, Israel



3DEEP Facial and Neck



Before treatment




After treatment

Photographs courtesy of  
Fiona Wright, M.D., Plano, Texas, USA

3DEEP Neck



Before treatment



7 months after 1 RF treatment

Photographs courtesy of  
Isabelle Rousseaux, M.D., Board Certified Dermatologist, Loos, France



## The iFine Treatment Handpiece by EndyMed

The iFine is designed to target fine lines and wrinkles found in delicate, thin-skinned, hard to reach areas such as around the eyes and mouth. The iFine provides:

- Smoother, tighter skin around the eyes.
- Reduced under eye sagging and improved eye contour.
- Reduced hooding to reveal more of the eyelid.
- More youthful, less tired looking eyes.
- Softening of fine lines enabling easier application of make-up



iFine Handpeice

## ↔ 3DEEP iFine Wrinkle Reduction



Before treatment



After 6 weekly treatment

Photographs courtesy of  
Li-Lin Liu, M.D., Taiwan



↔ 3DEEP iFine Wrinkle Reduction




Before treatment



After 3 treatments

Photographs courtesy of  
I.Vider, M.D., Medical OR, Herzlia, Israel

↔ 3DEEP iFine Wrinkle Reduction



Before treatment




Immediately after treatment


Photographs courtesy of  
Yoram Harth, M.D., Hertzlia, Israel



↔ 3DEEP iFine Wrinkle Reduction



Before treatment



Three weeks after one treatment

Photographs courtesy of  
I. Vider, M.D., Medical OR, Hertzlia, Israel

↔ 3DEEP iFine Wrinkle Reduction



Before treatment



After 3 treatments

Photographs courtesy of  
I. Vider, M.D., Medical OR, Hertzlia, Israel





## Physicians Testimonials



*"The EndyMed PRO 3DEEP technology uses a new sophisticated way to control RF energy in the tissue. Using the attraction and repulsion forces of multiple RF sources working in tandem it forces energy to the depth of the dermis with minimal heat on the surface. My clinical experience with the EndyMed PRO shows consistent skin tightening especially in the lower face and neck areas with 100% satisfied patients."*

**Prof. Daniel Cassuto, Professor of Plastic and Reconstructive Surgery at the University of Catania and a private practice in Milan, Italy**



*"I found the EndyMed PRO a very attractive option for painless skin tightening on the face and body. Its ease of use and lack of consumables makes it a perfect option for long lasting treatment of wrinkles, lax skin and cellulite."*

**Klaus Fritz, M.D., Board Certified Dermatologist, President of the European Society of Laser Dermatology (ESLD), Head Dermatology and Laser Clinic in Landau, Germany**



*"Our patients are delighted with the results of the 3DEEP treatment and I have very little difficulty in recommending this safe and effective treatment for skin tightening and facial scarring."*

**Christopher Inglefield, M.D., BSc FRCS (Plast), London Bridge Plastic Surgery, London, UK**

## Physicians Testimonials



*"EndyMed's 3DEEP RF fractional skin resurfacing applicator is one of the new, exciting technologies used for performing this increasingly popular procedure."*

**Neil Sadick, MD, Board Certified Dermatologist, Clinical Professor of Dermatology at Weill Cornell Medical College, NYC, New York, USA**



*"3DEEP technology enabled us to achieve focused energy in the depth of the tissue with reduced heat and less subsequent damage to the surface of the skin. This results in painless procedures with the need for cooling or topical anesthetics. Now we are able to use, with EndyMed PRO, RF energy more efficiently and effectively."*

**Flor Mayoral, M.D., Board Certified Dermatologist, Coral Gables, Florida**



*"I have used virtually every laser on the market worldwide and can take them apart and put them back together overall and I have to say that EndyMed's 3DEEP is the best total body tightening technology I have seen so far. 3DEEP is both user and patient friendly. Patients, male and female, are beating down my doors worldwide for us to treat them, and that only happens if a technology delivers what it promises with correct value as well. And hallelujah it is a workhorse that I don't have to nonstop troubleshoot! I am a satisfied customer and a hard one to please."*

**Julia T. Hunter, M.D., Skin Fitness Plus, Beverly Hills, California, USA**

# EndyMed Medical

EndyMed Medical is a cutting-edge medical device company that develops energy based aesthetic treatment systems for the professional and home use markets. EndyMed's products are powered by our proprietary, FDA cleared 3DEEP® multisource radiofrequency technology, the most advanced RF treatment technology available today for high efficacy wrinkle reduction and painless skin treatments.



EndyMed™ Medical Ltd., 7 Bareket St., POB 3582, North Industrial Park, Caesarea 30889, Israel  
[ t ] +972 4 630 9100 [ f ] +972 4 630 9101 [ e ] info@endymed.com [ w ] www.endymed.com

**EndyMed USA Toll Free Number - 800-753-8862**

© 2014, EndyMed Medical Ltd. All rights reserved.



Like us on Facebook at [www.facebook.com/EndyMed](http://www.facebook.com/EndyMed)

