

Step-by-Step Guide: Pairing Hearing Aids with iPhone

Step 1: Turn on Bluetooth

- Open the **Settings** app on your iPhone.
- Tap **Bluetooth**.
- Turn the switch on so it turns green.

Step 2: Open Hearing Devices Menu

- Go back to the main **Settings** page.
- Scroll down and tap **Accessibility**.
- Tap **Hearing Devices** (or MFi Hearing Devices).

Step 3: Put Hearing Aids in Pairing Mode

- **For rechargeable aids:** Turn them off and back on by placing them in their charger and taking them out, or by holding the power button.
- **For battery-door aids:** Open the battery doors for three seconds, then close them again.

Step 4: Connect and Pair

- Keep your hearing aids close to your iPhone.
- Wait for your name or the hearing aid model to appear on the screen.
- Tap the name of the hearing aids.
- Tap **Pair** on the pop-up message.
- Tap **Pair** a second time if you are connecting a left and right pair.

This is for informational purposes only. For medical advice or diagnosis, consult a professional. AI responses may include mistakes.