



Product Name: Particulate Respirator KN95

Packing:

2 pcs/polybag (Polybag independent packaging)

40 pcs/carton (Size (LWH): 22*12.2*14cm)

24 boxes/carton (Size (LWH): 50.8*46*42cm)

Classification:

- FFP2 NR D
- CE 2703

Feature:

1. Safe and environmental, no abnormal odor, excellent breathability

2. Three-dimensional tailor and design

a. Folding design

Folding design is easy to carry and store.

b. Nasal splint

Easy to fix, improved airtightness of the mask, thickened round corner design, easy to shape and durable.

c. Ear strap design

Improve comfort level while convenient to wear and take off.

d. Filtration technology

Outer filter layer technology.

3. Four-layer filtration technology

a. Skin-friendly non-woven fabric

b. Melt-blown fabric

c. Filter cotton

d. High quality non-woven fabric

4. One-time filtration efficiency > 95%