



## Infusion Pump Analyzer Chamber

3.5 mL Chamber  
MC3.5-1000

35 mL Chamber  
MC35-1000



8 pin Mini-Din  
connector for  
connection to  
IPA-1000



8 pin Mini-Din  
connector for  
connection to  
IPA-1000

Luer Lock  
Connector for  
Fluid In

### MC3.5-1000 & MC35-1000 CHAMBERS

FLOW		
RANGE	3.5 mL Chamber	0 to 999.9 mL/hr
	35 mL Chamber	15.0 to 9999.9 mL/hr
MEDIA	Distilled Water	
MEDIA CONNECTION	Luer Lock Fitting	

ELECTRICAL	
POWER	Powered by IPA-1000 Unit
CONNECTIONS	8-pin Female Mini-DIN Receptacle Connects to IPA-1000 Base Unit using BC20-40607 Mini-DIN Straight-Thru Interconnect Cable

PHYSICAL & ENVIRONMENTAL		
CONSTRUCTION	ENCLOSURE	ABS Plastic
	FACE PLATE	Back-printed Lexan
SIZE	14.00 x 5.00 x 3.50 Inches (355.6 x 127.0 x 88.9 mm)	
WEIGHT	< 2.2 Lbs (1.00 kg)	
OPERATING RANGE	15 to 30 °C (59 to 86 °F) 10 to 80% RH, Non-Condensing	
STORAGE RANGE	-40 to 60 °C ( -40 to 140 °F)	

Metphil Medical Company  
3F Insular Life Bldg., Abanao Extension  
Baguio City, Philippines  
Phone: +63 74 422 5751  
Email: sales@metphil.com  
Website: www.metphil.com

