The DigitQuant ProstaSense (DPS) Assay

To be comfort, accurate, and trustworthy

Who Should be Tested



Men

Age > 50

Blood PSA Level: 4-10ng/mL

The DPS Procedure




Non-DRE First Catch Urine Sample


Automated Target Specific RNA Isolation

Digital RNA Amplification for Absolute Quantification

What does the DPS Score mean?



The DPS Benefits



The DPS Advantage

	Sensitivity	Specificity	NPV	PPV	False Negative Rate	Avoid Biopsy
DPS	93.0%	63.2%	94.1%	58.8%	2.5%	42.9%
Competitor A	77.5%	57.1%*	90.0%	33.6%	4.9%	49.6%*
Competitor B	88.5%	38.0%	91.8%	29.9%	2.6%	29.4%
Competitor C	92.5%	30.2%	90.4%	36.1%*	2.3%*	23.4%
Competitor D	96.9%*	27.4%	95.9%*	33.7%	3.1%	N/A