

Patient Handling

Novis offers an extensive range of patient handling solutions, designed to provide safe and secure transfer. Focussing on user ergonomics, safety and therapeutic benefits, Novis specialises in transfer aids that provide the best in pressure care treatment and prevention whilst reducing the risk of injury involved in providing care.

- Patient Lifters
- Slings
- Horizontal Positioning
- Seated Positioning
- Seated Transfer
- Horizontal Transfer
- Sit to Stand Transfer
- Patient Support Aids
- Special Solutions



Manual handling represents one of the most important OHS problems facing the nursing profession.

Traditionally, in nursing, manual handling has been thought of as lifting and transferring patients.

A key approach in these standards is to acknowledge that manual handling includes much more than activities involving direct patient care.

It is important to utilise the correct handling equipment when transferring and moving a patient to ensure the safety and comfort for both patient and caregiver.

Manual Transfer Products – an Assessment Guide

Always identify whether the patient is 'able to assist' or is 'dependent', according to particular objective criteria. For these purposes, a patient would be identified as 'able to assist' if they are able to understand instructions, are cooperative and would be physically able to assist the process.

If they do not meet these criteria, they would be determined as 'dependent'.

Any special needs of the person that may affect the move should also be identified. For example, they may be resistant, prefer to move to a particular side or may be deaf.

Determine the patient's ability to assist through:

- Cognitive signs
- Physical signs
- Behavioural signs
- Clinical constraints

Select the safest moving/handling method. Once the initial assessment has been carried out, identify the lowest risk method appropriate to the situation.



Lifters

These Patient Handling devices assist in the care of someone with limited mobility. They provide a safe and effective method of moving, transferring and repositioning a user in a safe and controlled manner.

The LiftAid 180 and 360 are lightweight, easy to operate mobile lifters that allow transfers from the floor, bed or chair with ease and safety.

With low lifting height and soft start/stop function, this lifter has been designed to provide patient and carer with optimum comfort, safety and security.



Slings

Versatile enough for patient transfers, toileting or bathing, the ProSling range includes both pivot and yoke styles with durable embroidered information label. Slings are developed to suit the majority of lifter products on the market, and include three loop support straps and keyhole clips to suit yoke or pivot connections.




Transfer Aids

Designed to provide assistance, facilitate transfer and frequent repositioning, Transfer Aids improve patient care, while reducing strain on carers and improving quality of life for users.

Transfer Aids provide assistance to both the user and carer, improving ergonomics and reducing the risk of injury. The Novis range includes a large variety of exclusive and specialist transfer aids with a focus on ergonomics and facilitating care.

Air Transfer Mat

Novis Air Transfer Mat is inflated with low pressure and a high volume of air. Hover holes on the underside of the mat allow air to escape and the escaping air acts as a gentle air cushion layer, reducing friction, making the lateral movement of the mattress and the patient effortless.



Greater Staff Safety and Patient Comfort

- 80-90% less physical effort needed to perform transfers, greatly reducing injury risks to carers.
- The Novis Air Transfer Mat cradles and supports the patient, maintaining a stable position while transferring.
- Proven design eliminates shearing and friction forces providing more comfortable patient transfers.

Infection Control and Cost Effectiveness

- Single patient use mats deliver greater infection control management.
- Save time and costs on decontamination as no laundering is required.
- Number of staff required is significantly reduced during manual handling transfers representing cost and time efficiency.

Effective for a Wide Range of Patient Transfer Requirements

- Accommodates various patient sizes up to 450kg.
- Vapour permeable nylon helps manage microclimate while helping provide added protection against shear and friction.
- Heavy duty patient safety straps and handles.
- Radiolucent/MRI safe.

	AVAILABLE IN 3 SIZES (W X L)	CODE
Air Transfer Mat	864 x 1980 mm	ATM102D
	990 x 1980 mm	ATM103D
	1270 x 1980 mm	ATM104D
Blower Pump	Air Transfer Mat Blower Pump	ATMP001

*All Novis Air Transfer Mats are compatible with many of the most commonly available pump units



LiftAid 180

- Lightweight, easy to operate mobile patient lifter
- Powder coated aluminium frame with 3 adjustable height settings
- Quick release yoke spreader bar
- Swivel style front and rear castors; foot operated rear brake locks
- Electrically operated 'U' shaped base with wide flanged front legs
- Linak actuator systems with intelligent diagnostics and status indicators
- Handheld controller with simple push button operation
- Rechargeable battery featuring power status indicator
- Emergency stop and manual emergency lowering device
- Optional pivot frame
- Compatible with a range of slings (refer to pages 49-50)

LiftAid180

MHPH0418



LiftAid 320

- Heavy duty, fully featured mobile patient lifter
- Powder coated aluminium frame with 3 height settings
- Quick release yoke double spreader bar
- Swivel style front and rear castors; foot operated rear brake locks
- Electrically operated 'U' shaped base with wide flanged front legs
- Sturdy fold down armrests to assist with gait training
- Linak actuator systems with intelligent diagnostics and status indicators
- Handheld controller with simple push button operation
- Rechargeable battery featuring power status indicator
- Emergency stop and manual emergency lowering device
- Compatible with a range of slings (refer to pages 49-50)

LiftAid 320

MHPH0419




Pivot Frame

- Uses pivot slings with keyhole fittings
- Enhanced patient positioning and comfort
- Improved balance point
- Compatible with a range of pivot slings (refer to page 50)

Pivot Frame

MHPH0408



ProSling Range



Versatile enough for patient transfers, toileting or bathing, the ProSling range includes both pivot and yoke styles with durable embroidered information label. Slings are developed to suit the majority of lifter products on the market, and include three loop support straps and keyhole clips to suit yoke or pivot connections.

All ProSling products are subject to rigorous testing and standards (ISO 10535:2006) to guarantee the best quality product possible.

Slings are available in a range of sizes (with colour coded straps) to accommodate variety in patient height, shape and weight.

ProSling Yoke General Purpose

- General purpose sling for patients with basic upper body control

ProSling Yoke General Purpose with Head Support

- General purpose sling for patients with less control
- Inbuilt head support

ProSling Yoke Hygiene

- Designed for toilet transfers
- Open area from mid back to mid thigh
- Buckle waist strap for added security
- Generous under arm padding

ProSling Yoke Hygiene with Head Support

- Designed for toilet transfers for patients with less control
- Open area from mid back to mid thigh
- Inbuilt head support
- Buckle waist strap for added security
- Generous under arm padding

ProSling Yoke Disposable

- Disposable, general purpose sling for patients with less control
- Single patient use only (cannot be washed and re-used)
- Inbuilt head support
- Non-woven, spun bond polypropylene
- Instruction label with space for patient details

PROSLING YOKE
GENERAL PURPOSE



PROSLING YOKE GENERAL
PURPOSE WITH HEAD SUPPORT



PROSLING YOKE HYGIENE



PROSLING YOKE HYGIENE
WITH HEAD SUPPORT



PROSLING YOKE
DISPOSABLE



ProSling Range



ProSling Pivot General Purpose with Head Support

- General purpose for patients with basic upper body control
- Reinforced nylon resin pivot clips and buckles
- Adjustable straps for position control
- Inbuilt head support

ProSling Pivot Hygiene

- Designed for toilet transfer
- Open area from mid back to mid thigh
- Reinforced nylon resin pivot clips and buckles
- Increased positioning control for correct seating

PROSLING PIVOT HYGIENE



PROSLING PIVOT GENERAL PURPOSE WITH HEAD SUPPORT



Yoke	General Purpose	Small	SL147181
		Medium	SL147281
		Large	SL147381
		Extra Large	SL147481
		Extra Extra Large	SL147581
Yoke	General Purpose with Head Support	Small	SL146181
		Medium	SL146281
		Large	SL146381
		Extra Large	SL146481
		Extra Extra Large	SL146581
Yoke	Hygiene	Small	SL267181
		Medium	SL267281
		Large	SL267381
		Extra Large	SL267481
		Extra Extra Large	SL267581
Yoke	Hygiene with Head Support	Small	SL246181
		Medium	SL246281
		Large	SL246381
		Extra Large	SL246481
		Extra Extra Large	SL246581
Yoke	Disposable General Purpose with Head Support	Medium	SL156281
		Large *	SL156381
Pivot	Pivot General Purpose with Head Support	Small	SL146191
		Medium	SL146291
		Large	SL146391
Pivot	Pivot Hygiene	Small	SL247191
		Medium	SL247291
		Large	SL247391

* Maximum capacity 205 kg



Immedia SatinSheet System



- Complete 'in-bed system' for turning and positioning
- Carefully designed system of base sheets and draw sheets that work together in various combinations, depending on patient needs
- 'In-bed system' means no application or removal of product required
- Reduces friction between user and bed for reduced risk of injury and reduced shear
- Low friction satin material encourages movements
- Can be used to replace regular bed sheets



2Direction Flat Sheet	Mini	1950 x 1050 (Satin 600 x 1050) mm	IM4133S
		1950 x 1050 (Satin 700 x 1050) mm	IM4108S
		2300 x 1050 (Satin 850 x 1050) mm	IM4175S
	Midi	2000 x 1400 (Satin 600 x 1400) mm	IM4119S
		2000 x 1400 (Satin 700 x 1400) mm	IM4113S
		2300 x 1400 (Satin 850 x 1400) mm	IM4176S
	Maxi	2000 x 2000 (Satin 600 x 2000) mm	IM4181S
		2000 x 2000 (Satin 700 x 2000) mm	IM4114S
		2300 x 2000 (Satin 850 x 2000) mm	IM4177S
2Direction Flat Sheet - Double	Mini	2300 x 1050 (Satin 960 x 1050) mm	IM4178S
	Midi	2300 x 1400 (Satin 960 x 1400) mm	IM4179S
	Maxi	2300 x 2000 (Satin 960 x 2000) mm	IM4180S
2Direction Flat Sheet - Offset Double	Maxi	2800 x 2000 (Satin 600 x 2000) mm	IM4122S
2Direction Corner Sheet	Maxi	850 x 2000 (Satin 650 x 2000) mm	IM4107S
		1050 x 2000 (Satin 800 x 2000) mm	IM4147S
		1200 x 2000 (Satin 900 x 2000) mm	IM4137S
2Direction Fitted Sheet	Maxi - Fitted	900 x 2000 (Satin 700 x 2000) mm	IM4112S
	Maxi - Elastic Fit	900 x 2000 (Satin 700 x 2000) mm	IM4142S
2Direction - DrawSheet	Mini	2000 x 1000 mm	IM4116S
	Midi	2000 x 1400 mm	IM4118S
4Direction - DrawSheet	Midi	2000 x 1400 (Satin 1470 x 1400) mm	IM4115S
	Midi with Handles	2000 x 1400 (Satin 1470 x 1400) mm	IM4115SH
	Maxi	2000 x 2000 (Satin 1470 x 2000) mm	IM4121S
	Maxi with Handles	2000 x 2000 (Satin 1470 x 2000) mm	IM4121SH
	Mini with Wings	1850 x 900 mm	IM4110S
In2Sheet with Incontinence Pad - 4Direction	Mini with Wings/Handles	1850 x 900 mm	IM4110SH
	Midi with Wings	2200 x 1400 mm	IM4140S
	Midi with Wings/Handles	2200 x 1400 mm	IM4140SH
SatinSheet Pillow Case		600 x 700 mm	IM4117S
SatinSheet Sets	Mini Set - IM4108S + IM4116S		IM4201S
	Midi Set - IM4113S + IM4118S		IM4200S
	Maxi Set - IM4114S + IM4118S		IM4202S
	Maxi/Midi FIT 4D Set - IM4112S + IM4115S		IM4203S
	Midi 4D Set - IM4113S + IM4115S		IM4204S
	Midi Small Set - IM4119S + IM4118S		IM4205S
	Mini In2Sheet Set - IM4108S + IM4110S		IM4210S
	Midi/Mini In2Sheet Set - IM4119S + IM4110S		IM4209S
	Midi 4D Small Set - IM4119S + IM4115S		IM4207S
	Midi Corner Set - IM4107S + IM4118S		IM4208S
	Midi 4D Corner Set - IM4107S + IM4115S		IM4211S
	Maxi 4D Corner Set - IM4107S + IM4121S		IM4212S



Immedia 4WayGlide Midi and Mini



- Complete 'in-bed system' for turning and repositioning for independent or partially independent users
- Nylon Sheet fixed to mattress and Glide Mattress on top
- Non-slip nylon sheets minimises the risk of sliding out of bed
- 4WayGlide Mini for short in-bed system
- 4WayGlide Midi for short in-bed system for heavy users with partial mobility in legs

4WAYGLIDE MIDL



4WAYGLIDE MINI

4WayGlide Mini	Glide Mattress	1400 x 1000 mm	IM140/100
	Nylon Sheet with Non-Slip Sides	1400 x 1000 mm	IM140/100NS
	Set includes IM140/100 + IM140/100NS		IM4301
4WayGlide Midi	Glide Mattress	1400 x 1500 mm	IM140/150
	Nylon Sheet with Non-Slip Sides	1400 x 1500 mm	IM140/150NS
	Set Includes IM140/150 + IM140/150NS		IM4313
Accessories	Positioning Wedge with Non-Slip Cover	300 x 200 x 530 mm	IM53NS
		300 x 200 x 800 mm	IM53/80NS

Immedia 4WayGlide System



- Complete 'in-bed system' for turning and positioning for users with limited mobility
- Ultra low friction surface for movement in all directions
- Non-slip section (covered or uncovered)
- Elastic corner straps
- Standard and Maxi sizes
- Optional nylon sheet locking mechanism for greater control
- Optional glide mattress




4WayGlide LPL	Glide Mattress, Cotton/Polyester	1400 x 2000 mm	IM140/200LPL
	Glide Mattress with Length Handles, Cotton/Polyester	1400 x 2000 mm	IM140/200LH
	Glide Mattress, Polyurethane	1400 x 2000 mm	IM140/200LPLG
	Nylon Sheet with Non-Slip Lock and Non-Slip Sides	850 x 2000 mm	IM85/200LPLNS
	Nylon Sheet with Non-Slip Lock	850 x 2000 mm	IM85/200LPL
	Nylon Sheet with No Lock and Non-Slip Sides	850 x 2000 mm	IM85/200NS
	Nylon Sheet	850 x 2000 mm	IM85/200
	Nylon Sheet with 3 x Non-Slip Lock and Non-Slip Sides	850 x 2000 mm	IM85/2003LPL
	Nylon Sheet with 2 x Non-Slip Locks	850 x 2000 mm	IM85/2002LPL
4WayGlide LPL Maxi	Glide Mattress, Cotton/Polyester	1800 x 2100 mm	IM180/200LPL
	Glide Mattress, Polyurethane	1800 x 2000 mm	IM180/200GLPL
	Nylon Sheet with Non-Slip Lock and Non-Slip Sides	1100 x 2000 mm	IM110/200LPLNS
	Nylon Sheet with 3 x Non-Slip Lock and Non-Slip Sides	1200 x 2000 mm	IM120/2003LPL
Accessories	Positioning Wedge with Non-Slip Cover	300 x 200 x 530 mm	IM53NS
		300 x 200 x 800 mm	IM53/80NS
	LPL Cover	1400 x 2200 mm	IM140/200PU
4WayGlide Set	Standard Set - Includes IM140/200LPL + IM85/200LPLNS		IM4309
	3-Lock Set - Includes IM140/200LPL + IM85/2003LPLNS		IM4300

Immedia TwinSheet 4Glide system



- Complete 'in-bed system' for turning and positioning incontinent users with severely reduced mobility
- TwinSheet4Glide system is a 2-part glide mattress with incontinence pad and handles
- Also suited for users suffering dementia, brain damage or pressure injuries
- Good access for care of pressure injuries, hygiene, dressing and bedding changes
- For use on top of 4WayGlide or similar in-bed systems



TwinSheet 4Glide	Top Part Underpad	1500 x 750 mm	IM150/75N
	Lower Part Underpad	1500 x 1400 mm	IM150/140N
	Long Underpad	1500 x 2000 mm	IM150/200N
	LPL Nylon Sheet with Non-Slip Lock and Non-Slip Sides	850 x 2000 mm	IM85/200LPLNS
	LPL Nylon Sheet with Non-Slip Lock	850 x 2000 mm	IM85/200LPL
	LPL Nylon Sheet with 3 x Non Slip Locks and Non-Slip Sides	850 x 2000 mm	IM85/2003LPL
Accessories	Positioning Wedge with Non-Slip Cover	300 x 200 x 530 mm	IM53NS
		300 x 200 x 800 mm	IM53/80NS
TwinSheet 4Glide Sets	1 Lock Set - including IM150/75N + IM150/140N + IM85/200LPLNS		IM4312
	3 Lock Set - including IM150/75N + IM150/140N + IM85/2003LPL		IM4314
	Underpad Set - including IM150/75N + IM150/140N		IM4315
	3 Lock Set with Long Underpad - including IM150/200N + IM85/2003LPL		IM4316
	1 Lock Set with Long Underpad - including IM150/200N + IM85/200LPLNS		IM4317

In2Sheet Incontinence Pad



- Incontinence pad with satin low friction underside
- Ideal for turning incontinent user with reduced mobility
- 2 direction low friction surface
- Available with side handles for better grip



In2Sheet Incontinence Pad	900 x 850 mm	IM4109
In2Sheet Incontinence Pad with Handles	1050 x 1250 mm	IM4128

Slide Sheet with Handles



- Tubular high-low friction slide sheet with integrated handles for repositioning in bed
- High friction on outside; ultra low friction inside
- Waterproof barrier on outside
- Designed to facilitate transfers up and down in bed, X-ray table or MT-scanner
- Loop strap for easy storage and colour codes for differentiation



Multiglide Tubular with Handles	Pink Strap	1300 x 800 mm	IM130DR
	Light-Blue Strap	1600 x 800 mm	IM160DR

MultiGlide Slide Sheet



- ❑ Tubular high-low friction slide sheet for repositioning in bed
- ❑ High friction on outside, ultra low friction inside
- ❑ Waterproof barrier on outside
- ❑ Suited for transferring and turning in bed; between bed and wheelchair
- ❑ Loop strap for easy storage and colour codes for differentiation



MultiGlide Tubular	Blue Strap	400 x 400 mm	IM40
	Silver Strap	600 x 250 mm	IM62
	Yellow Strap	600 x 600 mm	IM66
	Brown Strap	720 x 700 mm	IM77
	Red Strap	720 x 1000 mm	IM80
	Light-Brown Strap	700 x 2000 mm	IM82
	Green Strap	720 x 1300 mm	IM83
	Silver Strap	800 x 800 mm	IM88
	White Strap	1000 x 700 mm	IM100
	Silver Strap	1000 x 2000 mm	IM102
	Green Strap	1000 x 1350 mm	IM103/1
	Grey Strap	1080 x 800 mm	IM108
	White Strap	1200 x 700 mm	IM120
	Pink Strap	1300 x 800 mm	IM130
	Light-Blue Strap	1600 x 800 mm	IM168
	Purple Strap	2000 x 700 mm	IM200
MultiGlide Open		1000 x 1400 mm	IM100/140

SG Slide Sheet



- ❑ Tubular low-low friction slide sheet for repositioning in bed
- ❑ Ultra low friction on outside, ultra low friction inside
- ❑ Use to minimise friction under pressure points
- ❑ Waterproof barrier on outside
- ❑ Easy to apply and remove
- ❑ SG tubular available in several sizes or cut to size roll
- ❑ Loop strap for easy storage and colour codes for differentiation



MultiGlide Tubular	Black Strap	400 x 400 mm	IM0040
		600 x 600 mm	IM0066
		720 x 700 mm	IM0077
		800 x 800 mm	IM0088
		720 x 1000 mm	IM0080
		720 x 1300 mm	IM0083
		1000 x 700 mm	IM00100
		1080 x 800 mm	IM00108
		1200 x 900 mm	IM00120/1
		1300 x 800 mm	IM00130
		1300 x 900 mm	IM00139
		1400 x 1200 mm	IM00170
		2000 x 700 mm	IM00200
MultiGlide Nylon Sheet (on a roll)		1.5 x 10 m	IM2245

BetterLiving Single Patient Use Slide Sheet



- ❑ Disposable design for improved infection control practices
- ❑ Small spills can be wiped down with a mild soapy solution
- ❑ Polyurethane coated polyester construction
- ❑ Heal sealed edges
- ❑ Large and practical size of 1.5 x 2 m

Single Patient Use Slide Sheet	BL1000
--------------------------------	--------



BetterLiving Washable Slide Sheet



- ❑ Machine wash and tumble dry (low)
- ❑ Silicone coated polyester construction
- ❑ Heal sealed and hemmed edges
- ❑ Large and practical size of 1.5 x 2 m

Washable Slide Sheet	BL1001
----------------------	--------




Immedia GlideCushion



- ❑ Padded glide cushion
- ❑ Suited for turning or moving in bed, or transfer to wheelchair
- ❑ Supports independent movement or aids assisted movement

GlideCushion	Polyester / Cotton	600 x 500 mm	IM60/50
		600 x 900 mm	IM60/90
		690 x 700 mm	IM69/70
		900 x 900 mm	IM90/90
GlideCushion X-Small Set – Includes IM60/50 + IM6000N			IM4303



Immedia TurningMattress



- ❑ Padded polyester/cotton turning mattress
- ❑ Suited for turning in bed
- ❑ Patient remains on top mattress in bed
- ❑ Ultra low friction helps users with reduced mobility to turn independently in bed

TurningMattress	920 x 1950 mm	IM92/195
-----------------	---------------	----------



Immedia AutoTurn



- Soft swivel cushion to facilitate car entry / exit
- Low friction surface allows user to slide and turn
- Side handles for assisted positioning
- Belt and buckle for attachment

AutoTurn	450 mm	IM99999
----------	--------	---------




1



2



3



4

Immedia EasyTurn



- Soft swivel cushion to facilitate bed or car entry / exit
- Portable and easy to use in any location
- Non-slip underside with soft velour top
- Low friction movement lowers strain on body

	400 mm	IM999405
EasyTurn	450 mm	IM990454
	500 mm	IM999505



Immedia PediTurn




- Floor swivel cushion that facilitates transfers between seated positions
- For users with reduced movement, especially in feet
- Reduces strain on knee and foot joints
- Non-slip top and bottom for stability

	Soft	IM99836GT
PediTurn	Hard	IM99700



PEDITURN HARD

PEDITURN SOFT



Immedia Glide Pad

- ▣ Suited for transfer between bed and wheelchair; in and out of car
- ▣ Used independently or in combination with transfer board

Swan Glide Pad with Handles 450 x 1400 mm IM2698121




Immedia OneWayGlide

- ▣ Promotes controlled and easy positioning
- ▣ One way slide direction; resists sliding in other direction
- ▣ OneWayGlide Non-Slip for tubular cushion with non-slip outer that grips to any surface
- ▣ OneWayGlide Velour for tubular cushion with velour material for added comfort
- ▣ OneWayGlide Pad for one layer pad with velour top for use directly on textile surface.
- ▣ OneWayGlide Long for patients requiring repeated chair repositioning



ONEWAYGLIDE VELOUR TUBULAR CUSHION



ONEWAYGLIDE NON-SLIP TUBULAR CUSHION



ONEWAYGLIDE PAD VELOUR

	430 x 370 mm	IM39/8
OneWayGlide Tube Non-Slip	450 x 450 mm	IM38/8
	600 x 800 mm	IM37/8
OneWayGlide Tube Velour	430 x 370 mm	IM45/8
	500 x 420 mm	IM50/8
OneWayGlide Pad Velour	430 x 370 mm	IM44/8
	400 x 1300 mm	IM34/8
OneWayGlide Long	450 x 1300 mm	IM51/8
	550 x 1300 mm	IM35/8
OneWayGlide Long, Velour, Tie Strings	400 x 1300 mm	IM33/8



Immedia 3B Board



- ▣ Traditional curved, low friction transfer board
- ▣ Longer length ideal for wider gaps
- ▣ Position curve forward or backward for multiple transfer situations
- ▣ Suitable for patients with sufficient upper body strength and balance for self transfer



3B Board	230 x 770 mm	IM403
----------	--------------	-------

Immedia Butterfly Transfer Board



- ▣ Long, wide, stable transfer board with wheel cut outs for wheelchair positioning
- ▣ Smooth top surface with large sitting area
- ▣ Non-slip underside for stability
- ▣ Butterfly GlideWing as optional accessory for reduced transfer friction



Transfer Board		IM401
GlideWing	Large, 350 x 400 mm	IM4001
	X-Large, 350 x 475 mm	IM4002

Immedia E-Board



- ▣ Thin, flexible transfer board for easier application and removal
- ▣ Lightweight and compact for an easy to carry solution
- ▣ E-Board with Wings provides added stability and protection against wheel of wheelchair



E-Board	Small, 200 x 450 mm	IM410
	Medium, 250 x 600 mm	IM408
E-Board with Wings	Wide, 330 x 600 mm	IM404
	Long, 330 x 750 mm	IM405

Immedia Dolphin Toilet Transfer Board



- ▣ Transfer board designed for toilet transfers
- ▣ Smooth fibreglass top
- ▣ Unique cut-out design allows board to remain during toileting
- ▣ Non-slip underside
- ▣ Cut-outs for secure wheelchair positioning



Dolphin Board	690 mm	IM409
---------------	--------	-------

Immedia 2Move Transfer Board



- Transfer board for horizontal transfer
- Used to bridge gap and even height differences between two horizontal surfaces
- Simple to transport and store
- Available in 2 lengths; foldable and wide options available
- Used with a low friction cover



DISPOSABLE COVER IN ROLL



2MOVE, LONG, TRANSFER BOARD, INCLUDING NYLON COVER

2Move Transfer Board	Small	500 x 950 mm	IM8050095
	Large Foldable	500 x 1900 mm	IM8050190
	Large Non-Foldable	500 x 1900 mm	IM8150190
	X-Large Foldable	700 x 1900 mm	IM8070190
2Move Covers	Nylon Cover, Waterproof	for IM8050095	IM4047
		for IM8050190 & IM8150190	IM4048
		for IM8070190	IM4049
	Disposable Cover	for IM8050095, roll (100 pieces)	IM4035
		for IM8050190 & IM8150190, roll (40 pieces)	IM4034
		for IM8070190, roll (40 pieces)	IM4036
2Move Set	Soft Padded Cover	for IM8050190 + IM8150190	IM60/195
	Small Standard Set – Includes IM80500905 + IM4047		IM4319
	Large Standard Set – Includes IM8050190 + IM4048		IM4318
	Large Non-Foldable Set – Includes IM8150190 & IM4048		IM4320
	X-Large Foldable Set – Includes IM8070190 & IM4049		IM4321

Immedia TransferMattress



- Soft padded mattress for horizontal transfer
- Transfer between two supine positions or for turning in bed
- Padded cotton/polyester material
- Combine with waterproof nylon or disposable cover

Narrow	600 x 1950 mm	IM60/195
Wide	700 x 1950 mm	IM70/195
Nylon cover for IM60/195	600 x 2200 mm	IM60220N
Nylon cover for IM70/195	700 x 2000 mm	IM7000N
Narrow Set – Includes IM60/195 + IM60220N		IM4304
Wide Set – Includes IM70/195 + IM7000N		IM4305



Etac Turner PRO



- Turning and transfer aid for weight bearing patients
- Lightweight and compact at only 7.5 kilograms
- Minimal adjustment between users
- Reduces risk of falls and workplace injury
- Ergonomic working position for carers, while maintaining eye contact
- Enhanced compact design with moulded parts for easy cleaning
- Contrasting orange coloured handle for high visibility
- Soft, comfortable individual angled leg supports

Turner PRO

16090106



Immedia Sling with Handles



- Flexible support aid for bed positioning and transfers
- Four handles each side for easy placement and removal
- Low-friction outer padding
- Non-slip, inner padding

Sling with Handles, Long

160 x 1900 mm

IM428

Sling with Handles, Short

160 x 1200 mm

IM428K



Immedia SupportBelt



- Mobilisation and transfer aid to help patients stand or sit from a wheelchair, chair or bed
- General walking and exercise support
- Non-slip inner surface with unique locking system
- Ergonomically angled handles for a better grip

Support Belt	XXS, 4 Handles	450 to 900 mm at waist	IM8009009
	XS, 4 Handles	650 to 1050 mm at waist	IM8010509
	S, 4 Handles	650 to 1050 mm at waist	IM8010513
	M, 6 Handles	900 to 1300 mm at waist	IM8013013
	L, 8 Handles	1250 to 1500 mm at waist	IM8015013
	XL, 10 Handles	1450 to 1700 mm at waist	IM8017015



Immedia SmartMove



- Support pad with handles
- Non-slip material for safe positioning
- Acts as arm extension for support without sacrifice to ergonomic grip
- Used for a range of transfers

SmartMove

210 x 600 mm

IM426



Immedia BedString

- An aid for the active user to independently turn or sit up in bed
- Long strap with handles for self-propelled movement
- Style attachment suited to most bed frames
- Useful for application of sling or for dressing

BedString with 8 Handles	85 x 1700 mm	IM429
BedString with 2 Handles	50 x 1800 mm	IM9910001



Immedia HandyGrip

- Adjustable grip to support independent turning in bed
- Useful for self transfer from wheelchair to bed
- Long strap with adjustable buckle
- Both length and position varied for optimum grip

HandyGrip	50 x 1600 mm	IM461
-----------	--------------	-------



Immedia SlingOn

- Aid to facilitate the application of a sling with minimal movement
- Flexible tool allowing ergonomic grip for caregiver with 4 handles each side
- Used to support various transfers and positioning

SlingOn	Small	800 x 1000 mm	IM437/S
	Medium	800 x 1700 mm	IM438
	Large	1200 x 1500 mm	IM437



Immedia WiseMover

- In-bed positioning for bedridden patients needing frequent positioning and transfer
- Store in bed head and fold out when required
- Suited for turning, positioning and horizontal transfer
- Also for insertion of X-ray cassette when patient immobilised

WiseMover	850 x 2000 mm	IM441
-----------	---------------	-------



Immedia Mini Board



- ❑ Soft transfer stretcher
- ❑ Suited for transfer between bed and operating table, or lifting from the floor
- ❑ Board combined with low friction material for easy insertion

4 Handles, Rectangular	540 x 1100 mm	IM8011054
6 Handles, Rectangular	540 x 1500 mm	IM8015054
7 Handles, with Head Support	540 x 1240 mm	IM8012454
11 Handles, with Head Support	540 x 2000 mm	IM8020054



Immedia OneManSling



- ❑ Unique tool for repositioning patients higher up in bed with minimal effort
- ❑ Padded nylon transfer strap with 6 handles
- ❑ Ideal for use with one carer only
- ❑ Prevents friction and shear

OneManSling with 6 Handles	75 x 3000 mm	IM460PU
----------------------------	--------------	---------



Immedia MultiGlide Gloves



- ❑ Multi-purpose product for different manual handling situations
- ❑ Ideal for small re-positioning adjustments and movements
- ❑ Low friction external surface material
- ❑ One size fits all

MultiGlide Gloves, Waterproof, Pair	IM35
-------------------------------------	------



Immedia Thorax SlideSecure



- ❑ Aids X-ray cassette insertion for bedridden patient
- ❑ Fixed to bed head for easy storage and access
- ❑ Can be used with SlideSecure Bag for easier application


Thorax SlideSecure	900 x 800 mm	IM4008SECURE
Thorax SlideSecure Bag	560 x 470 mm	IM4017



Immedia Mole X-Ray Board

- ❑ X-ray glideboard for easy X-ray cassette positioning in bed
- ❑ Transparent top board to control cassette positioning
- ❑ Lower board to hold cassette in place
- ❑ Can be used with SlideSecure Bag for easier application

Mole X-Ray Board	400 x 800 mm	IM407
------------------	--------------	-------



Immedia RescueSheet

- ❑ Emergency evacuation slide system
- ❑ Remains under mattress for use in emergency
- ❑ Transfer patient and mattress out of room and down stairs
- ❑ Waterproof

RescueSheet	Standard	750 x 2000 mm	IM108088
	Extra Wide	850 x 2000 mm	IM108089



1



2



3

Immedia HipTurner

- ❑ Aid to turn and support patient in static position on side
- ❑ Ideal for bed pan insertion or wound dressing
- ❑ Support good carer positioning without strain
- ❑ Can be used in conjunction with patient lifter

HipTurner	400 x 200 x 45 mm	IM462
-----------	-------------------	-------



Immedia NonSlip



- ▣ High friction pad that can be used when increased friction is desired
- ▣ Use on floor for increased resistance or on mattress for increased manoeuvrability
- ▣ Roll can be cut to desired size

NonSlip Square	380 x 480 mm	IM2248
	1.5 x 5 m	IM2246
NonSlip Roll	1.5 x 30 m	IM2247



Immedia CloverFoot



- ▣ Stool for elevating foot position during transfer
- ▣ Can be used to elevate and support feet during toileting
- ▣ Angled shape to allow placement close to obstacle

CloverFoot	380 x 90 x 320 mm	IM2001
------------	-------------------	--------

