


CERTIFIED  
PRODUCTS



FOREMOST  
**ANESTHESIA**  
VENTILATOR



CERTIFIED  
PRODUCTS

MOUNT HERE



**EASILY ADJUSTABLE  
WITH ANY UNIVERSAL  
BOYLE ANESTHESIA  
MACHINE**

**SAFE DURABLE & COMPACT DESIGN**

# TECHNICAL SPECIFICATIONS




FOREMOST  
**ANESTHESIA**  
VENTILATOR

* AVAILABLE MODES	:	VCV, PCV, VACV, VSIMV, PSIMV, PACV, CPAP
* PEEP	:	0 TO 30 cm H <sub>2</sub> O (ELECTRONIC)
* PRESSURE SUPPORT	:	0 TO 50 cm H <sub>2</sub> O
* RESPIRATORY RATE	:	1 TO 80 BPM
* INSPIRATORY PRESSURE	:	5 TO 100 cm H <sub>2</sub> O High pressure limit Low pressure limit
* TIDAL VOLUME	:	40 TO 2000 ml
* INSPIRATORY TIME	:	0.1 TO 9.9 sec
* I : E RATIO	:	1:9.9 -2.5:1
* TRIGGER	:	PRESSURE (0 TO -20CM H <sub>2</sub> O)
* FIO <sub>2</sub>	:	21% TO 100%
* DIMENSIONS (W*D*H)	:	723 mm* 408 mm* 132 mm
* APNEA BACKUP	:	YES
* DETACHABLE AND AUTO CLAVABLE EXPIRATORY VALVE :		NO
* RECHARGABLE Li-ion BATTERY:		YES WITH UPTO 12 HOURS BACKUP

## A true life saver in the battle of OT




FOREMOST  
**ANASTHESIA**  
VENTILATOR



**EASILY ADJUSTABLE  
WITH ANY UNIVERSAL  
BOYLE ANESTHESIA  
MACHINE  
SAFE, DURABLE & COMPACT DESIGN**

FOREMOST MEDITECH with its delivery reach to almost all serviceable pin codes in India is continuously striving for Product Excellence, Innovation, Precision and Perfection. We are among the top 25 brands of ventilators in the world and the top most leading brand in India.



Mfd. By



Mktd. By

**Medical Engineering And Devices**

# City Off: 4/18 Harimohan Dutta Road, Kolkata 700028 Factory : 20/19, Jessore Road By Lane, Habra, North 24 Parganas , 743263

visit us: [www.med-india.co.in](http://www.med-india.co.in)

email us: [info@med-india.co.in](mailto:info@med-india.co.in)

contact nos: +91- 74777 63529

### Our Certifications:

