

Syringe & Volume Infusion Pump

MED
A COMPANY BY MEDICAL PROFESSIONALS



*Designed to Perform
Uninterrupted*

Sflow-10

SYRINGE PUMP

MED
A COMPANY BY MEDICAL PROFESSIONALS

A company By Medical Professionals

www.med-india.co.in



Follow us on



WWW.MED-INDIA.CO.IN



**Compact, Accurate,
Most Easy to use**

Product Specifications:

Product Name	MED Sflow-10 Syringe Pump
Brand	MED
Model	Sflow-10
Syringe Size	5,10,20,30,50/60ml
Flow Rate Range	5ml Syringe 0.1ml/h-150ml/h 10ml Syringe 0.1ml/h-300ml/h 20ml Syringe 0.1ml/h-600ml/h 30ml Syringe 0.1ml/h-900ml/h 50/60ml Syringe 0.1ml/h-1200ml/h
Syringe Capacities	5 ml, 10 ml, 20 ml (30 ml), 50 ml/60 ml
Accuracy	±2%
Infusion Mode	Rate Mode, Rate-Time, Rate-VTBI, Time-VTBI, Body weight
Volume to be Infused	0-9999.9ml
Bolus or Purge Rate	50 ml syringe: From 0.1 to 1800.0 ml/h 30 ml syringe: From 0.1 to 1200.0 ml/h 20ml syringe: From 0.1 to 800.0 ml/h 10ml syringe: From 0.1 to 400.0 ml/h 5 ml syringe: From 0.1 to 200.0 ml/h Increment: 0.1 ml/h for all syringe types
Bolus Volume	From 1.0 to 100.0 ml, 0.1 ml increment
Drug Library	20 drug list with drug code display
KVO Rate	0.1-5.0ml/h Adjustable
patient weight	from 0.01 to 99.9 kg, 0.1 increment from 100 to 300kg, 1 increment
Solution volume	from 0.01 to 199.9 ml, 0.1 ml increment
pressure Model	Low, Medium, High
Display	LCD Screen
Stackability	Yes
Dimension weight	32x21.5x20cm Approx.3 kg
Battery	Lithium Battery 11,1/2000mAh
Classification	Class I, type Cf, IPX4
Power Spply	AC : 100-240V, 50Hz/60Hz

Medical Engineering And Devices

Regd. Office: 4/18, Harimohan Dutta Road,
Kolkata-700028 Website: www.med-india.co.in M: 74777 63529
Email: info@med-india.co.in



Channel Partner



Branches: Odisha (7439531744), Jharkhand (6296380007), Guwahati (9679762035), Bihar (9679761803), Siliguri (6295027297)

© 2021 Medical Engineering And Devices All rights reserved. Medical Engineering And Devices reserves the right to make changes in product features, specifications, aesthetics and/or to discontinue the same at any time without notice or obligation