

Volume Infusion Pump

MED

A COMPANY BY MEDICAL PROFESSIONALS



*Designed to Perform
Uninterrupted*

Vflow-10

VOLUME INFUSION PUMP

MED

A COMPANY BY MEDICAL PROFESSIONALS

A company By Medical Professionals

www.med-india.co.in



Follow us on



WWW.MED-INDIA.CO.IN



**Compact, Accurate,
Most Easy to use**

Product Specifications:

Product Name	MED Vflow-10 Infusion Pump
Brand	MED
Model	Vflow-10
Dimension & Weight	13x17.5x23cm: 2kg
Occlusion	High: 40 kpa+-20kpa Middle: 40 kpa+-20kpa Low: 100kpa+-20kpa
Total volume infused	0-9999ml
Accuracy	±2%
Infusion Mode	Rate Mode, Rate-Time, Rate-VTBI, Time-VTBI, Body weight, Drop Volume, Drop Time
Volume to be Infused	0-9999.9ml
Bolus or Purge Rate	1-1200ml/h default, 800ml/h, step by 1ml/h
Bolus Volume	1.0-10ml adjustable, default 3ml
Drug Library	20 drug list with drug code display
KVO Rate	1ml/h-5ml/h, default value 1ml/h, can be programmed by user step 0.1ml/h
patient weight	from 0.01 to 99.9 kg, 0.1 increment from 100 to 300kg, 1 increment
Solution volume	from 0.01 to 199.9 ml, 0.1 ml increment
pressure Model	Low, Medium, High
Display	LCD Screen
Stackability	Yes
Dimension weight	32x21.5x20cm Approx.3 kg
Battery	Lithium Battery 11,1/2000mAh
Classification	Class I, type Cf, IPX4
Power Supply	AC : 100-240V, 50Hz/60Hz

Medical Engineering And Devices

Regd. Office: 4/18, Harimohan Dutta Road,
Kolkata-700028 Website: www.med-india.co.in M: 74777 63529
Email: info@med-india.co.in



Channel Partner



Branches: Odisha (7439531744), Jharkhand (6296380007), Guwahati (9679762035), Bihar (9679761803), Siliguri (6295027297)

© 2021 Medical Engineering And Devices All rights reserved. Medical Engineering And Devices reserves the right to make changes in product features, specifications, aesthetics and/or to discontinue the same at any time without notice or obligation