

NEOMED

A Complete Neonatal Solution

**Radiant Warmer
Phototherapy
Open Care System
Incubator
Bubble CPAP
Neonatal Ventilator**

Experience the new generation

MED



Follow us on



www.med-india.co.in

Babylus Warmer

Model: 02/03/PRO




- Made of FRP material
- Long lasting Ceramic heater
- Parabolic reflector
- Shock & Fire Proof *
- Servo & Manual mode
- Skin & Air Temperature
- Dazzle free LED examination lamp
- Display Range 15 to 40 Deg C
- Set temp range(Servo) 32 to 38 C
- Manual mode 0% to 100%
- Attached Stand for patient monitor
- Attached IV stand
- Tiltable bed for easy procedure
- Robust and Intuitive easy to use
- Audio visual alarm for safety
- Long lasting probe easy to clean
- Operating voltage 240V Ac, 50Hz

Model: PT100/200

Lumino Phototherapy

- Made of FRP material
- With intensity controller
- Four examination LED
- Uniform intensity > 40 uW/cm
- Affective treatment area 50cm X 25cm
- Peak wavelength 450-455nm
- Tiltable head 0 to 90 deg
- Adjustable height
- LED life > 20000 hours
- Available with optional trolley
- Available with optional dual surface
- Robust and Intuitive easy to use
- Long lasting No fan cooling
- Operating voltage 240V
- Max power consumption



Medical Engineering And Devices

Regd. Office: 4/18, Harimohan Dutta Road,
Kolkata-700028 Website: www.med-india.co.in
M: 74777 63529
Email: info@med-india.co.in



Branches: Odisha (7439531744), Jharkhand (6296380007), Guwahati (9679762035), Bihar (9679761803), Siliguri (6295027297)

© 2021 Medical Engineering And Devices All rights reserved. Medical Engineering And Devices reserves the right to make changes in product features, specifications, aesthetics and/or to discontinue the same at any time without notice or obligation