


VentTrax 360

Ventilator Alarm Alert and Compliance Monitoring

The VentTrax Display provides alarm type, location, duration, and monitoring status. In addition, patient alarm activity is captured and displayed. The self supervising systems also displays transmitter battery condition, connectivity, and system status. The alarm alert information is easy to read and changes color as the event time increases.

VentTrax Battery Operated Wireless Room Transmitters Support All Ventilators



VentTrax Display Monitors Located at Nursing Station and Hallways

System Status Display

Fault "RED"
Low Battery "YELLOW"

Alarm Display

1-20 Seconds "GREEN"
21-40 Seconds "YELLOW"
41-60 Seconds "RED"
61 > Seconds "ALT RED/WHITE"

Room Status Display

Not In Use "GRAY"
Monitoring "GREEN"
Frequent Alarms "ALT GREEN/WHITE"

System Status	Alarm Type / Location / Duration	Room Status
Power OK		101A 109A
Com OK		101B 109B
108B Fault	VENT 101A	102A 110A
115A Low Batt	0:47	102B 110B
		103A 111A
	VENT 105B	103B 111B
	0:26	104A 112A
		104B 112B
	VENT 102B	105A 113A
	0:14	105B 113B
		106A 114A
	VENT 113A	106B 114B
	0:02	107A 115A
		107B 115B
		108A 116A
		108B 116B
		VENTS 29

VentTrax / Sai Systems Technologies

VentTrax 360

Ventilator Alarm Alert and Compliance Monitoring

Benefits

Administration

Expandable systems, add transmitters as you grow

No hard wiring required

Improved patient safety

Reduce facility noise

Comprehensive assessment tools to measure compliance and success

Cloud-based reporting

As low as \$399/month

Nursing

No Alarm Fatigue from numerous remote boxes

Easily Identify alarm location from up to 200 feet away

Ventilator clears, alarm clears, No wasted steps

Pagers display alarm location even in parking lot

Alarms are easily identified from **ANY** location in the facility

Patient

No annoying alarms

Discrete display boards

Improved care

Improved safety

Please email us at pat.devaney@predictivehealthdevices.com to learn more



VentTrax is proudly made in the USA