

**Medline Hudson RCI Air-Eze Incentive Deep Breathe Exerciser**

- Wide flow rate range from 200 to 1200 cc/sec
- Low resistance design

- Sliding incentive indicator
- Compact design

[Report Incorrect Product Information](#)**Hudson RCI + Medline**The Hudson RCI brand is now part of the Medline family. [Learn More.](#)

Packaging	12 Each / Case
Manufacturer	Teleflex Medical
Manufacturer #	8884715700
Comparable Item(s)	HUD8884719011

SKU HUD715700

**Specifications**

CONNECTOR TYPE	Universal
LATEX FREE	No
PRODUCT TYPE	Flowmeter
SIZE	Universal
UNSPSC	42271604

Unit of Measure	Conversion	Net/Gross Weight (lbs)	Volume (cubic ft)	Shipping Dimensions (inch) L x W x H	GTIN
Each (EA)	1.0 Each	0.001 / 0.371	0.072	5.002 x 5.002 x 5.002	00000000000000
Case (CS)	12.0 Each	0.012 / 4.45	0.869	17.3 x 12.4 x 7.0	00000000000000

**HCPCS Information****HCPCS Code: NO GOVERNMENT OR SUPPLIER HCPCS CODE VERIFICATION**

**HCPCS Disclaimer:** HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.

**Documents/SDS**

No documents/SDS are available for this product