

# Personal Protective Apparel

Head-to-Toe Protection for Healthcare



**Know  
you're  
covered.**



# Protect your people **and** **your bottom line.**

You want to know you're protected. Our protective apparel goes through rigorous testing and meets strict FDA guidelines. We follow the highest quality standards to make sure you get the protection you need.

As a leader in protective apparel, we also have the most comprehensive product line and low acquisition costs. You can find the exact products you need—and maintain cost efficiencies.



AAMI Level  
isolation gowns  
are color-coded  
by neck binding

**Recommended Gowns** 4-5

**AAMI Isolation Gowns** 6-7

**Cover Gowns** 8-9

**Protective Gowns** 10-11

**Footwear** 12-13

**Headwear and Beard Cover** 14-15

**Scrubs** 16

**Coveralls** 17

**Lab Coats** 18-19

**Breast Cancer Month  
Awareness PPE** 20-21

**Face Masks** 22

**Patient Wear** 23



# Understanding gown materials

## Spunbond polypropylene

For basic coverage, this non-woven fabric bonds fibers together to form a single layer.

## SMS material

Spunbond/meltblown/spunbond (SMS) is a multi-layer fabric composed of inner layers of meltblown polypropylene between outer layers of spunbond polypropylene.

## Coated polypropylene (poly-coated)

Soft, spunbond polypropylene is coated with a layer of polyethylene (plastic) film.

## Polyethylene film

Polyethylene (plastic) film is composed of a single layer of film and depending on thickness can be appropriate to low to moderate fluid exposure.



# Choose the right gown

Our recommendations for your workplace environment.

## A Housewide, Med/Surg and ICU

Our QuickComply™ AAMI Level 2 gowns feature an over-the-head style which improves compliance while reducing donning and removing time. Side ties eliminate the difficulties of tying behind the back, while thumb loop wrists allow for a good gown and glove fit.

Item No.	Description
NONLV240	Regular/Large
NONLV240XL	X-Large

*Additional items on page 7*

## B Endoscopy and Emergency

AAMI Level 3 gowns are made from comfortable and breathable premium multi-ply material. Features a classic tie waist, hook and loop neck and knit cuffs.

Item No.	Description
NONLV325	Regular/Large
NONLV325XL	X-Large

*Additional items on page 7*

## C Decontamination

Made from premium polyethylene coated materials that have been tested against ASTM F1670 and ASTM F1671. An over-the-head, open back style reduces donning and removing time, while thumb loop wrists allow for a good gown and glove fit.

Item No.	Description
NON26276	Regular/Large
NON26276XL	X-Large
NON26276XXL	XX-Large

*Additional items on page 11*

## D Chemotherapy

Protective, polyethylene-coated materials & seams tested against top 11 chemotherapy drugs per ASTM F739-12. Features hook and loop neck, waist ties and knit cuffs. Taped seams to meet USP 800 guidelines.

Item No.	Description
NON25457M	Medium
NON25457	Regular/Large
NON25457XL	X-Large
NON25457XXL	XX-Large

*Additional items on page 11*





**A** Housewide, Med/Surg and ICU



**B** Endoscopy and Emergency



**C** Decontamination



**D** Chemotherapy

## Why AAMI?

### AAMI Guidelines: A Reliable Benchmark for Protection

Under AAMI Guidelines, protective apparel must meet minimum performance standards for strength, barrier protection and fluid resistance to ensure the safety of healthcare workers and patients.



## Levels of Barrier Protection – AAMI PB70:2012

AAMI Level	Test Method	Requirements at 4% AQL	Anticipated Risk of Fluid Exposure
<b>Level 1</b> 💧	AATCC 42 - Water Resistance: Impact Penetration	≤ 4 g	Low
<b>Level 2</b> 💧💧	AATCC 42 - Water Resistance: Impact Penetration	≤ 1.0 g	Low to Moderate
	AATCC127 - Water Resistance: Hydrostatic Pressure	≥ 20 cm	
<b>Level 3</b> 💧💧💧	AATCC 42 - Water Resistance: Impact Penetration	≤ 1.0 g	Moderate to High
	AATCC127 - Water Resistance: Hydrostatic Pressure	≥ 50 cm	



## AAMI Level 1 Isolation Gown

Recommended for: Med/Surg, Visitors, Contact Precautions

- Made from multi-ply material
- Side ties eliminate the difficulties of tying behind the back
- Choose between over-the-head or neck tie styles

Item No.	Description	Pkg.
<b>Medium Weight, Neck Ties, Elastic Wrists</b>		
NONLV100	Regular/Large, Yellow	100/cs
NONLV100XL	X-Large, Yellow	100/cs
<b>Medium Weight, Over-the-Head, Thumb Loop Wrists</b>		
NONLV140	Regular/Large, Yellow	100/cs
NONLV140XL	X-Large, Yellow	100/cs

## AAMI Level 2 Isolation Gown

Recommended for: Med/Surg, ICU, Emergency, Airborne and Contact Precautions

- Made from multi-ply material with side ties to eliminate the difficulties of tying behind the back
- Choose between over-the-head or tape tab neck styles
- Available with elastic wrists or thumb loops

Item No.	Description	Pkg.
<b>Medium Weight, Tape Tab Neck, Elastic Wrists</b>		
NONLV200M	Medium, Yellow	100/cs
NONLV200	Regular/Large, Yellow	100/cs
NONLV200XL	X-Large, Yellow	100/cs
NONLV200XXL	XX-Large, Yellow	100/cs
NONLV200XXXL	XXX-Large, Yellow	100/cs
<b>Medium Weight, Over-the-head, Overlapping Back panels, Thumb loops</b>		
NONLV235	Regular/Large, Yellow	100/cs
NONLV235XL	X-Large, Yellow	100/cs
<b>Medium Weight, Over-the-Head, Thumb Loops</b>		
NONLV240	Regular/Large, Yellow	100/cs
NONLV240XL	X-Large, Yellow	100/cs
NONLV240XXL	XX-Large, Yellow	100/cs
NONLV240XXXL	XXX-Large, Yellow	100/cs

## AAMI Level 2 Polyethylene-Coated Isolation Gown

Poly Coated recommended for: Decontam/SPD, Endoscopy, Trauma and ER

- Available over-the-head
- Available with thumb loops or knit cuffs

Item No.	Description	Pkg.
<b>Polyethylene Coated, Over-the-Head, Knit Cuff</b>		
NON27340	Regular/Large, Yellow	100/cs
NON27340XL	X-Large, Yellow	100/cs
<b>Polyethylene Coated, Over-the-Head, Thumb Loops</b>		
NON27360	Regular/Large, Yellow	100/cs
NON27360XL	X-Large, Yellow	100/cs
<b>Polyethylene Coated, Over-the-Head, Thumb Loops with Elastic Wrists</b>		
NON27380	Regular/Large, Blue	100/cs
NON27380XL	X-Large, Blue	100/cs

## AAMI Level 3 Isolation Gown

Recommended for: ICU, Emergency/Trauma, Airborne, Contact and Droplet Precautions

- Comfortable and breathable materials with classic tie waist
- Choose between tape tab, hook and loop and over the head styles
- Available with elastic wrists, knit cuffs, or thumb loops

Item No.	Description	Pkg.
<b>Heavyweight Multi-Ply, Tape Tab Neck, Elastic Wrists</b>		
NONLV315	Regular/Large, Yellow	100/cs
NONLV315XL	X-Large, Yellow	100/cs
<b>Premium Multi-Ply, Hook and Loop Neck, Knit Cuffs</b>		
NONLV325	Regular/Large, Blue	50/cs
NONLV325XL	X-Large, Blue	50/cs
<b>Medium Weight Multi-Ply, Over-The-Head, Thumb Loop Wrists</b>		
NONLV340	Regular/Large, Yellow	50/cs
<b>Microporous, Tape Tab Neck, Knit Cuffs</b>		
NONLV350	Regular/Large, White	50/cs
NONLV350XL	X-Large, White	50/cs

## Gown Features



### Side ties

Available on AAMI Level 1, 2 and 3 gowns



### Neck tie style

Available on AAMI Level 1 gowns



### Thumb loop wrists

Available on AAMI Level 1, 2 and 3 gowns



### Tape tab neck style

Available on AAMI Level 2 and 3 gowns



CRI4000



Neck ties help ensure the gown remains securely in place.

### Classic Cover Gown

Recommended for: Visitors, Food Services

- Made from lightweight polypropylene and multi-layer materials
- Generously sized for full coverage and flexibility
- Side ties eliminate difficulties of tying behind the back
- Neck ties and elastic wrists

Item No.	Description	Pkg.
<b>Lightweight Polypropylene</b>		
CRI4000	Regular/Large, Yellow	50/cs
CRI4001	X-Large, Yellow	50/cs
CRI4003	XXX-Large, Yellow	50/cs
<b>Medium Weight Multi-Layer</b>		
NON27SMS2	Regular/Large, Yellow	100/cs
NON27SMS2XL	X-Large, Yellow	100/cs
NON27SMS2XXXL	XXX-Large, Yellow	100/cs





NON27SMS5



Convenient over-the-head, open back style reduces donning and removing time.

### Over-the-Head, Open Back Cover Gown

Recommended for: Visitors, Food Services

- Available in medium weight multi-layer
- Generously sized for coverage and flexibility
- Over-the-head style reduces donning and removing time
- Thumb loop wrists

Item No.	Description	Pkg.
<b>Medium Weight Multi-Layer</b>		
NON27SMS5	Regular/Large	100/cs
NON27SMS5XL	X-Large	100/cs

## Polyethylene Gowns

Recommended for: Decontamination, Emergency, ICU, Food Services, Airborne and Contact Precautions and Droplet Precautions

- For blood and bodily fluid penetration
- Over-the-head style reduces donning and removing time
- Thumb loop wrists help to prevent sleeves from sliding
- Available in bulk or individually folded in a dispenser box

Item No.	Description	Pkg.
<b>Polyethylene, Blue</b>		
CRI5000	Regular/Large	75/cs
CRI5000B	Regular/Large	100/cs
CRI5010	Universal	75/cs
CRI5001	X-Large	75/cs
CRI5001B	X-Large	100/cs
CRI5003	XXX-Large	75/cs
<b>Premium Polyethylene, Blue</b>		
NONTH150	Regular/Large	75/cs
NONTH180	Universal	75/cs
NONTH200	X-Large	75/cs
NONTH100B	Regular/Large	100/cs
NONTH200B	X-Large	100/cs

## Over-the-Head Microporous Gown

Recommended for: Endoscopy, GI Lab, ICU, Emergency, Maternity

- Made of microporous breathable material tested against ASTM F1670
- Over-the-head style reduces donning and removing time
- Thumb loop wrists prevent sleeve slide
- Individually folded and in a dispenser box

Item No.	Description	Pkg.
NONTH350	Regular/Large	40/cs
NONTH400	X-Large	40/cs

Learn more  
about bouffant  
caps and shoe  
covers.

See pages 12–14



CRI5010

NONTH350



## Polyethylene-Coated Gowns

Recommended for: Decontamination, ICU, Endoscopy, GI Lab

- Neck and side ties
- Available with elastic wrists or knit cuffs

Item No.	Description	Pkg.
<b>Polyethylene Coated, Yellow</b>		
CRI4010	Elastic Wrists, Regular/Large	50/cs
CRI4011	Elastic Wrists, X-Large	50/cs
CRI4010CF	Knit Cuffs, Regular/Large	50/cs
CRI4011CF	Knit Cuffs, X-Large	50/cs
CRI4012CF	Knit Cuffs, XX-Large	50/cs

## Chemotherapy Gowns

Recommended for: Oncology, Radiology, Lab, Pharmacy

- Laminate ChemoAire™ material tested against ASTM F1670 & F1671
- Materials and seams tested against top 11 chemotherapy drugs per ASTM F739-12 and ASTM D6978
- Hook and loop neck, waist ties and knit cuffs
- Taped seams to meet USP 800 guidelines

Item No.	Description	Pkg.
<b>Protective Polyethylene Coated AAMI Level 2</b>		
NON25457M	Medium	30/cs
NON25457	Regular/Large	30/cs
NON25457XL	X-Large	30/cs
NON25457XXL	XX-Large	30/cs
<b>Protective Premium Breathable Material</b>		
NON27457M	Medium	30/cs
NON27457	Regular/Large	30/cs
NON27457XL	X-Large	30/cs
NON27457XXL	XX-Large	30/cs
<b>ChemoAire™ Premium Breathable Material</b>		
NON29457	Regular/Large	30/cs
NON29457XL	X-Large	30/cs

Chemotherapy gowns tested against the following chemotherapy drugs per ASTM F739-12:

- |                              |                       |
|------------------------------|-----------------------|
| • Dacarbazine                | • Cisplatin           |
| • Carmustine (BiCNU)         | • Etoposide (Toposar) |
| • Cyclophosphamide (Cytoxan) | • Paclitaxel (Taxol)  |
| • Doxorubicin Hydrochloride  | • Thio-Tepa           |
| • 5-Fluorouracil             | • Mitomycin C         |
|                              | • Vincristine Sulfate |

## Over-the-Head Polyethylene-Coated Gowns

Recommended for: Decontamination, Emergency, ICU, Maternity

- Over-the-head style reduces donning and removing time
- Available with thumb loop wrists, elastic wrists with thumb loop, or knit cuffs

Item No.	Description	Pkg.
CRI4500	Thumb Loops, Regular/Large, Yellow	100/cs
CRI4600	Elastic wrists with Thumb loops, Universal	100/cs
<b>Premium Polyethylene Coated</b>		
NON26278	Knit Cuffs, Regular/Large, Blue	50/cs
NON26278XL	Knit Cuffs, X-Large, Blue	50/cs
NON26278XXL	Knit Cuffs, XX-Large, Blue	50/cs
NON26276	Thumb Loops, Regular/Large, Blue	50/cs
NON26276XL	Thumb Loops, X-Large, Blue	50/cs
NON26276XXL	Thumb Loops, XX-Large, Blue	50/cs

Meets USP  
800 Require-  
ments

Meets  
USP 800  
Requirements





NON29859

### Multi-Layer Shoe Covers

- Made of multi-layer material
- Heavyweight material features stripe non-skid pattern for increased traction
- Range of sizes to fit all shoe types

Item No.	Description	Pkg.
NON29858	Heavyweight, Regular/Large	100/bx, 300/cs
NON29859	Heavyweight, X-Large	100/bx, 200/cs
NON28858	Medium Weight, Regular/Large	100/bx, 300/cs
NON28859	Medium Weight, X-Large	100/bx, 200/cs
NON28852	Medium Weight, XX-Large	100/bx, 200/cs



CRI2002



CRI2000

### Spunbond Polypropylene Shoe Covers

- Made of spunbond polypropylene
- Choice of non-skid pattern or smooth bottom

Item No.	Description	Pkg.
NON28758	Non-Skid, Regular/Large	100/bx, 300/cs
NON28759	Non-Skid, X-Large	100/bx, 200/cs
NON28752	Non-Skid, XX-Large	100/bx, 200/cs
CRI2002	Non-Skid, Regular/Large	100/bx, 300/cs
CRI2003	Non-Skid, X-Large	100/bx, 200/cs
CRI2000	Smooth Bottom, Regular/Large	100/bx, 300/cs
CRI2001	Smooth Bottom, X-Large	100/bx, 200/cs

Size	Shoe Size Conversion	
	Women's	Men's
Regular/Large	Fits up to size 13.5	Fits up to size 12
X-Large	Fits up to size 16.5	Fits up to size 15



CRI2010

### Polyethylene Shoe Cover

- Smooth bottom for situations where traction is not essential
- Made from polyethylene material designed to withstand high levels of fluids

Item No.	Description	Pkg.
CRI2010	One Size Fits Most	1,000/cs

### Polyethylene-Coated Shoe Covers

- Medium-weight coated spunbond polypropylene
- Available with seam on bottom or seamless
- Non-skid treated for better traction

Item No.	Description	Pkg.
<b>Non-Skid/Bottom Seam</b>		
NON30858	Regular/Large	300/cs
NON30859	X-Large	300/cs
<b>Non-Skid/Seamless</b>		
NON31858	Regular/Large	300/cs
NON31859	X-Large	300/cs
NON31852	XX-Large	300/cs



NON31852



## Knee-High Boot Covers

- Choose between non-skid bottom or foam strips for increased traction
- Available with hook and loop closure or elastic ankle

Item No.	Description	Pkg.
NON27348P	Premium Breathable, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27348PXL	Premium Breathable, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27349PXL	Thigh High Premium Breathable, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27200	Blue Laminate, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27200XL	Blue Laminate, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27144	Poly-coated, Non-skid, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27144XL	Poly-coated, Non-skid, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27150	Poly-coated, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27150XL	Poly-coated, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27143	Poly-coated to the Ankle, SMS Material to the Knee, Non-skid, Elastic, Regular/Large	50/bx, 150/cs
NON27143XL	Poly-coated to the Ankle, SMS Material to the Knee, Non-skid, Elastic, X-Large	50/bx, 150/cs



NON27348P



NON27144



NON27200




NON27143



NON28229




CRI1009




NON30233B



NON28226



CRI1004




CRI1005




NON28233F



NON28230R



CRI1009



CRI1207

### Multi-Layer Bouffant Caps

- Variety of prints and sizes; all come 100/bx, 500/cs
- Choice of heavy or medium weight multi-ply or heavyweight single-ply

Item No.	Description	Pkg.
NON30232B	21", Blue	500/cs
NON30233B	24", Blue	500/cs
NON30223B	28", Blue	500/cs
NON28229	24", Brush Stroke	500/cs
NON28230R	24", Breast Cancer Ribbon	500/cs
NON28533B	24", Blue, Heavyweight	500/cs

### Polypropylene Bouffant Caps

- Variety of colors and sizes, all come 100/bx, 500/cs
- Choice of medium weight polypropylene or lightweight polypropylene

Item No.	Description	Pkg.
CRI1000	21", Green	500/cs
CRI1001	21", Blue	500/cs
CRI1002	21", White	500/cs
CRI1003	24", White	500/cs
CRI1004	24", Blue	500/cs
CRI1005	24", Green	500/cs
CRI1009	24", Red	500/cs
CRI1202	24", Blue, Medium Weight	500/cs
CRI1203	24", White, Heavyweight	500/cs
CRI1207	24", Purple, Medium Weight	500/cs



**Surgical headwear  
that meets  
AORN Guidelines  
for Surgical  
Headwear**

### Surgical Headwear

- Better coverage from the back of the neck to sideburn area
- One size fits most—provides all-day breathable comfort and a secure fit
- No elastic band in front means cap won't leave a mark
- Quality, longer-lasting design creates less waste than bouffant caps

Item No.	Description	Pkg.
NON61950	Premium Blue Material	600/cs
NON61990	Scrim Material	500/cs
NON61980	Multi-layer Material	500/cs



NON61980



NON61950



NON61990

### Surgeon Caps

- Available in scrim, multi-layer or spunlace materials
- One size fits most
- Ties in the back
- Polypropylene crown for breathability

Item No.	Description	Pkg.
NON28625	Scrim Material for Absorbency	500/cs
NON28626	Multi-layer Material for Fluid-resistance	500/cs
NON28628	Spunlace Material for Comfort	500/cs



NON28625



NONSH100C



NONSH200



NONSH600



NONSH400

### Surgeon Hood and Beard Covers

- Made of multi-layer or polypropylene material
- One size fits most

Item No.	Description	Pkg.
NONSH100C	Surgeons Hood with Multi-Layer Material, Polypropylene Crown	300/cs
NONSH200	Multi-layer Hood	300/cs
NONSH600	Polypropylene Hood	300/cs
NONSH400	Polypropylene Beard Cover	1,000/cs

# Scrubwear

## ComfortSure™ Scrubs

- V-neck top with 3 pockets
- Drawstring pant with pockets
- Made from premium ComfortSure™ material that provides soft, comfortable protection

Item No.	Description	Size	Pkg.
<b>ComfortSure™ Disposable Scrubs, V-Neck, Blue</b>			
NON57202XS	Scrub Shirt	XS	50/cs
NON57203XS	Scrub Pants		50/cs
NON57202S	Scrub Shirt	S	50/cs
NON57203S	Scrub Pants		50/cs
NON57202M	Scrub Shirt	M	50/cs
NON57203M	Scrub Pants		50/cs
NON57202L	Scrub Shirt	L	50/cs
NON57203L	Scrub Pants		50/cs
NON57202XL	Scrub Shirt	XL	50/cs
NON57203XL	Scrub Pants		50/cs
NON57202XXL	Scrub Shirt	XXL	50/cs
NON57203XXL	Scrub Pants		50/cs



NON57202M

NON27212L

NON37202M

## Traditional Scrubs

- Choice of v-neck or round neck top with chest pocket
- Choice of elastic waist or drawstring pant with two pockets
- Made from multi-ply material

### Scrub Top, Round Neck, Blue

NON27212S	Small	30/cs
NON27212M	Medium	30/cs
NON27212L	Large	30/cs
NON27212XL	X-Large	30/cs
NON27212XXL	XX-Large	30/cs
NON27212XXXL	XXX-Large	30/cs

### Scrub Top, V-Neck, Blue

NON27202S	Small	30/cs
NON27202M	Medium	30/cs
NON27202L	Large	30/cs
NON27202XL	X-Large	30/cs
NON27202XXL	XX-Large	30/cs

### Scrub Pant, Elastic Waist, Blue

NON27213S	Small	30/cs
-----------	-------	-------

NON27213M	Medium	30/cs
NON27213L	Large	30/cs
NON27213XL	X-Large	30/cs
NON27213XXL	XX-Large	30/cs
NON27213XXXL	XXX-Large	30/cs

### Scrub Pant, Drawstring Waist, Blue

NON27203S	Small	30/cs
NON27203M	Medium	30/cs
NON27203L	Large	30/cs
NON27203XL	X-Large	30/cs
NON27203XXL	XX-Large	30/cs

### Scrub Top, Round Neck, Wine Red

NON37212S	Small	30/cs
NON37212M	Medium	30/cs
NON37212L	Large	30/cs
NON37212XL	X-Large	30/cs
NON37212XXL	XX-Large	30/cs
NON37212XXXL	XXX-Large	30/cs

### Scrub Top, V-Neck, Wine Red

NON37202S	Small	30/cs
-----------	-------	-------

NON37202M	Medium	30/cs
NON37202L	Large	30/cs
NON37202XL	X-Large	30/cs
NON37202XXL	XX-Large	30/cs
NON37202XXXL	XXX-Large	30/cs

### Scrub Pant, Elastic Waist, Wine Red

NON37213S	Small	30/cs
NON37213M	Medium	30/cs
NON37213L	Large	30/cs
NON37213XL	X-Large	30/cs
NON37213XXL	XX-Large	30/cs
NON37213XXXL	XXX-Large	30/cs

### Scrub Pant, Drawstring Waist, Wine Red

NON37203S	Small	30/cs
NON37203M	Medium	30/cs
NON37203L	Large	30/cs
NON37203XL	X-Large	30/cs
NON37203XXL	XX-Large	30/cs



# Coveralls

## Prohibit™ Coveralls

Recommended for: Lab, Facilities Management, Maintenance and Staff

- Comfortable yet protective for high fluid or viral hazards

Item No.	Description	Pkg.
NONCVP2M	Small/Medium	25/cs
NONCVP2L	Large	25/cs
NONCVP2XL	X-Large	25/cs
NONCVP2XXL	XX-Large	25/cs
NONCVP2XXXL	XXX-Large	25/cs
NONCVP2XXXXL	XXXX-Large	25/cs

## Microporous Coveralls

Recommended for: Lab, Maintenance, and Visitors

- Choose between elastic wrists and elastic ankles or straight wrists and straight ankles

Item No.	Description	Pkg.
<b>Elastic Wrists and Ankles</b>		
NONCV700L	Large	25/cs
NONCV700XL	X-Large	25/cs
NONCV700XXL	XX-Large	25/cs
NONCV700XXXL	XXX-Large	25/cs

## Multi-Layer Coveralls

Recommended for: Staff and Visitors

- Made from multi-layer polypropylene
- Elastic wrist and straight ankle

Item No.	Description	Pkg.
NONCV250XL	X-Large	25/cs
NONCV250XXL	XX-Large	25/cs



NONCVP2L

NONCV150L

## Spunbond Coveralls

Recommended for: Staff and Visitors

- Available in heavy weight polypropylene material
- Choose between elastic wrist and straight ankle or straight wrist and straight ankle

Item No.	Description	Pkg.
<b>Heavyweight, Elastic Wrists, Straight Ankles</b>		
NONCV150L	Large	25/cs
NONCV150XL	X-Large	25/cs
NONCV150XXL	XX-Large	25/cs
NONCV150XXXL	XXX-Large	25/cs
NONCV150XXXXL	XXXX-Large	25/cs

- | Item No.    | Description | Pkg.  |
|-------------|-------------|-------|
| NONRP602S   | Small       | 30/cs |
| NONRP602M   | Medium      | 30/cs |
| NONRP602L   | Large       | 30/cs |
| NONRP602XL  | X-Large     | 30/cs |
| NONRP602XXL | XX-Large    | 30/cs |

## Premium Lab Jackets with Knit Collar and Cuffs

- | Item No.    | Description | Pkg.  |
|-------------|-------------|-------|
| NONRP601S   | Small       | 30/cs |
| NONRP601M   | Medium      | 30/cs |
| NONRP601L   | Large       | 30/cs |
| NONRP601XL  | X-Large     | 30/cs |
| NONRP601XXL | XX-Large    | 30/cs |





### Classic Lab Coats with Knit Collar and Cuffs

- Multi-layer material provides comfort and protection
- Three reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
<b>White</b>		
NONSW500S	Small	30/cs
NONSW500M	Medium	30/cs
NONSW500L	Large	30/cs
NONSW500XL	X-Large	30/cs
NONSW500XXL	XX-Large	30/cs
NONSW500XXXL	XXX-Large	30/cs
<b>Blue</b>		
NONSW600S	Small	30/cs
NONSW600M	Medium	30/cs
NONSW600L	Large	30/cs
NONSW600XL	X-Large	30/cs
NONSW600XXL	XX-Large	30/cs

### Classic Lab Jackets with Knit Collar and Cuffs

- Multi-layer material provides comfort and protection
- Three reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
<b>White</b>		
NONRP500S	Small	30/cs
NONRP500M	Medium	30/cs
NONRP500L	Large	30/cs
NONRP500XL	X-Large	30/cs
NONRP500XXL	XX-Large	30/cs
<b>Blue</b>		
NONRP600S	Small	30/cs
NONRP600M	Medium	30/cs
NONRP600L	Large	30/cs
NONRP600XL	X-Large	30/cs
NONRP600XXL	XX-Large	30/cs

### Classic Lab Coats with Traditional Collar and Knit Cuffs

- Multi-layer material provides comfort and protection
- Three reinforced pockets
- Convenient loop for hanging

Item No.	Description	Pkg.
<b>White, Snap Front</b>		
NONSW100S	Small	30/cs
NONSW100M	Medium	30/cs
NONSW100L	Large	30/cs
NONSW100XL	X-Large	30/cs
NONSW100XXL	XX-Large	30/cs
NONSW100XXXL	XXX-Large	30/cs
<b>Blue, Snap Front</b>		
NONSW400S	Small	30/cs
NONSW400M	Medium	30/cs
NONSW400L	Large	30/cs
NONSW400XL	X-Large	30/cs
NONSW400XXL	XX-Large	30/cs

### Protective Sleeve Covers

- Provides maximum amounts of protection in high fluid scenarios, or if extra protection is needed
- Available in premium breathable laminate or microporous material
- Elastic bands at both ends secure the covers at the wrists and upper arms



NON18SLEEVE1

Item No.	Description	Pkg.
NON18P2SLEEVE	20" (50.8 cm) Long Blue Breathable Laminate Material	200/cs
NON18SLEEVE1	20" (50.8 cm) Long White Microporous Breathable Material	200/cs

# Breast Cancer Awareness Month Pink PPE

Protection with a cause.



## Staff Apparel

### Lab Coat

- Pink lab coat made from multi-layer material
- Hits right below the knee
- Snap front closure
- Knit collar and cuffs

Item No.	Description	Pkg.
NONSW700S	Small	30/cs
NONSW700M	Medium	30/cs
NONSW700L	Large	30/cs
NONSW700XL	X-Large	30/cs
NONSW700XXL	XX-Large	30/cs



NONSW700

## Patient Apparel

### Patient Cape

- Patient cape made from multi-layer material
- Sleeveless with tape tab neck
- 31" x 21"

Item No.	Description	Pkg.
NON27038	One Size Fits Most	100/cs

### Patient Gown

- Patient gown made from multi-layer material
- Sleeveless with neck and waist ties

Item No.	Description	Pkg.
NON27048	Regular/Large	50/cs



NON27048



NON27038



## Bouffant Caps

### Pink Ribbon Bouffant Cap

- Bouffant cap made with lightweight polypropylene material

Item No.	Description	Pkg.
NON28230R	Size 24"	500/cs



NON28230R

## Procedure Face Masks

### Face Mask with Ear Loops

- Face mask with ear loops

Item No.	Description	Pkg.
NON27712EL	Pink/White	300/cs



NON27712EL

## Surgical Face Masks

### Face Mask with Ties

- Face mask with ties

Item No.	Description	Pkg.
NONRIBBON	Pink Ribbon	300/cs
NON27712	Pink/White	300/cs



NONRIBBON



NON27712

# Face Masks

## Procedure Face Masks



**NON27122**

### Basic Procedure Mask

Spunbond Polypropylene outer, thermalbond polypropylene Inner.

- Earloop Style
- 50/bx, 6 bx/cs



**NON27378**

### Basic Procedure Mask

Spunbond Polypropylene outer, thermalbond polypropylene Inner.

- Earloop Style
- 50/bx, 6 bx/cs



**NON27375**

### Basic Procedure Mask

Spunbond Polypropylene outer, cellulose inner.

- Earloop Style
- 50/bx, 6 bx/cs

## Surgical Face Masks



**NON27402**

### Basic Surgical Mask

Spunbond Polypropylene outer, cellulose inner.

- Ties
- 50/bx, 6 bx/cs



**NON27385**

### Hypoallergenic Surgical Mask

Cellulose inner and outer.

- Ties
- 50/bx, 6 bx/cs



**NON27386**

### Sensitive Skin Surgical Mask

Spunbond Polypropylene outer, thermalbond polypropylene Inner.

- Ties
- 50/bx, 6 bx/cs



# Patient Wear

Thoughtfully designed in a variety of sizes to fit every shape.

## Multi-Layer Patient Gowns

- Disposable material has a fabric-like feel for comfort
- Opening in front or back for easy access
- Attached ties for a custom fit

Item No.	Description	Pkg.
NON27046	Gown, Sleeveless, Waist Tie, Regular/Large, Blue	50/cs
NON27048	Gown, Sleeveless, Waist Tie, Regular/Large, Pink	50/cs
NON27146	Gown, Sleeveless, Neck and Waist Tie, Regular/Large, Blue	50/cs
NON27146SL	Gown, Short Sleeves, Neck and Waist Tie, Regular/Large, Blue	50/cs
NON27146SLXL	Gown, Short Sleeves, Neck and Waist Tie, XL, Blue	50/cs
NON27346SL	3 Arm Hole Gown, Sleeveless, Regular/Large, Blue	50/cs



NON27048

NON27146SL

## Reinforced Tissue Patient Gowns

- 3-ply tissue-poly-tissue (TPT) material aids in fluid absorption and protection
- Scrim material provides added strength and absorption
- Opening in front or back for easy access

Adjustable ties for a custom fit

Item No.	Description	Material	Pkg.
NON24356	Gown, Waist Tie, 30" x 42", Blue	TPT	50/cs
NON24355	Gown, Waist Tie, 30" x 42", White	TPT	50/cs
NON24252	Gown, Waist Tie, 30" x 42", Mauve	TPT	50/cs
NON24354	X-Ray Gown, Waist Tie, 30" x 42", Blue	TP	50/cs
NON254572	Bariatric Gown, Waist Tie, 72" x 45", Blue	Scrim	25/cs



NON24356



NON24355



NON24252

# Medline PPE Companion™

Protection that's always by your side.

Streamline caregiving. Move PPE closer to the point of care. This unique isolation cart has all you need to protect yourself—and your patients.

## Isolation Cart

### A Precaution Identification

Multilingual instructions show the clinically appropriate level of personal protection.

### B Organized Storage

Multilingual labels on each drawer.

### C Hand Protection

Three glove sizes. More storage behind.

### D Face Protection

At-a-glance supply of face masks with shields and earloop masks, ample storage for N95 respirators and more.

### E Body Protection

Gowns dispense one at a time with a touch and a pull.

### F Gown Storage

Additional storage for a minimum of 50 gowns.

### G Quick Disinfection

Sanitize hands. Disinfect surfaces.



NONISOORG CART

# Disposable Isolation Station

- Keeps hallways free from isolation cart congestion
- Organizer stores flat-packed isolation gowns, three gloves boxes, and one box of face masks
- 3 rows of pockets with eyelets for hanging
- Made of SMS material



NONISOCADDY

**Know you're protected.** Contact your Medline representative for more information.



**Medline Industries, Inc.**  
Three Lakes Drive, Northfield, IL 60093  
Medline United States | 1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

**Medline Canada**  
1-800-268-2848 | medline.ca | canada@medline.com  
**Medline México**  
01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    

#### REFERENCE:

1. Considerations for selecting protective clothing used in healthcare for protection against microorganisms in blood and body fluids. Centers for Disease Control and Prevention website. Table 3: ANSI/AAMI PB 70:12 Classification of barrier performance of surgical gowns, other protective apparel, surgical drapes and drape accessories. Available at: <http://www.cdc.gov/niosh/npptl/topics/protectiveclothing/>. Accessed May 3, 2016.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.

© 2019 Medline Industries, Inc. Companion and QuickComply are trademarks and Medline is a registered trademark of Medline Industries, Inc. MKT19W1334766 / LIT171 / 1.5M / TC / 75