



Dr. Device & Lab
D.D.L. Medical Supplies



Dr. Device & Lab

Company Profile

Founded in 2020, DDL Medical Supplies is a Malaysia based company that specializes in distributing, marketing and supplying scientific & laboratory supplies to various hospitals, diagnostics laboratories, and private sectors with a range of products and laboratory supplies such as laboratory consumables products, laboratory plasticware, laboratory glassware and blood collection tubes. We are certified with GDPMD (Good Distribution Practice For Medical Device), compliant distribute medical devices with MDA Malaysia approval.

Mission

To deliver quality medical laboratory equipment supplies, services and solutions through our commitment to integrity, honesty, responsibility and care for our customers. To add value and to create an ideal long-term customer-supplier relationship.

Vision

To expand our company's footprint and become one of the most leading medical supplies company in Malaysia with variety of medical and laboratory supplies.



Certificate No.:
20431EN44020230325



License No: MDA-4599-WDP123



Gel and Clot Activator Tube >>>

Usage: Applying to the acquisition of high-quality serum, specimen for biochemical, immunological and serum clinical examination. Serum separation gel is a kind of semi-solid inert gel, composed of several compounds and synthesized by medical polymeric material. After blood centrifugal, the specimen can be separated into serum top, separation gel middle and blood corpuscle lower. In the middle of the tube, a thick layer of separation gel can completely isolate serum from blood corpuscle, so as to ensure the purity of

Code	Material	Additive	Size(mm)	Vol	Packing
TYK05	PET	Gel & Clot Activator	13*100	5-7ml	100pcs/rack 1800pcs/ctn
TYK06	PET	Gel & Clot Activator	16*100	8-10ml	100pcs/rack 1200pcs/ctn



Edta Tube >>>

Usage: applying to clinical haematology test, such as blood routine examination, all kinds of blood cells analysis (such as blood cell morphology examination, blood cell count, hematocrit, hemoglobin, sicking test, etc.)

Code	Material	Additive	Size(mm)	Vol	Packing
TYK11	PET	EDTA K2/K3	13*75	1-5ml	100pcs/rack 1800pcs/ctn



Heparin Tube >>>

Usage: applying to clinical biochemistry, blood specimen collection for emergency biochemical tests and anticoagulation(blood rheology measurement, biochemical test, genetic diagnosis, enzymology, enzyme immunoassay, etc.)

Code	Material	Additive	Size(mm)	Vol	Packing
TYK14	PET	Heparin lithium	13*75	3-5ml	100pcs/rack 1800pcs/ctn



PT Tube >>>

Usage: Applying to blood clotting (PT, APTT, coagulation factor) test for clinical examination. The concentration of sodium citrate is 109mol/L and the accurate ratio of blood and anticoagulant is 9:1.

Code	Material	Additive	Size(mm)	Vol	Packing
TYK10	PET	Sodium Citrate1:9	13*75	1.8-4.5ml	100pcs/rack 1800pcs/ctn



Glucose Tube >>>

Usage: Applying to clinical determination of blood glucose. The original characters of blood specimen can be maintained for a long time.

Code	Material	Additive	Size(mm)	Vol	Packing
TYK17	PET	Sodium fluoride /potassium oxalate	13*75	2-5ml	100pcs/rack 1800pcs/ctn

Blood Collection Needle



Code	Description	Packing
TYJ09	Sterilization: gamma radiation Gauge: 18g-24g	100pcs/box, 5000pcs/ctn

Sample Handling Series



Code	Description	Packing
NWM05	Biohazard Bag 15x25cm	100pcs/bag, 2500pcs/ctn

Centrifuge Tube




Code	Description	Packing
TYA03	Centrifuge Tube 15ml, Conical Bottom, Bulk Pack	100pcs/bag, 2000pcs/ctn
TYA01	Centrifuge Tube 50ml, Free Standing, Bulk Pack	50pcs/bag, 500pcs/ctn
TYA02	Centrifuge Tube 50ml, Conical Bottom, Bulk Pack	50pcs/bag, 500pcs/ctn

Microscope Slide & Cover Glass



Code	Description	Packing	
TYI-7105	Polished Edges, Frosted One End On One Side	50pcs/box, 50box/ctn	
TYI-7101	Polished Edges, Plain	50pcs/box, 50box/ctn	
TYI-2424C	Cover Glass 24x24mm	100pcs/box, 500box/ctn	
TYI-2450C	Cover Glass 24x50mm	100pcs/box, 250box/ctn	




Code	Description	Packing
B-JZH01	Inoculation Loop 1ul, EO Sterile	10pcs/bag, 6000pcs/ctn
B-JZH10	Inoculation Loop 10ul, EO Sterile	10pcs/bag, 6000pcs/ctn




Code	Description	Packing	
TYC07	Eppendorf Pipette Tip 200ul	1000pcs/bag, 40000pcs/ctn	
TYC08	Eppendorf Pipette Tip 1000ul	500pcs/bag, 15000pcs/ctn	
TYC02	Gilson Pipette Tip 200ul	1000pcs/bag, 40000pcs/ctn	



Code	Description	Packing	
TYC22	1ml Transfer Pipette	500pcs/bag, 5000pcs/ctn	
TYC24	3ml Transfer Pipette	500pcs/bag, 5000pcs/ctn	

Tourniquet Series



Code	Description	Packing
TYT-1	Disposable Tourniquet, 50pcs pre-cut	50pcs/box, 60boxes/ctn




D.D.L. Medical Supplies

Ground Floor, Lot 348, Lorong 7A, Jalan Ang Cheng Ho,
93100 Kuching, Sarawak, Malaysia
E-mail: drdevicemed@gmail.com



Certificate No.:
20431EN44020230325



License No: MDA-4599-WDP123