

Manual S.M.P.

Sistema Multi Plataforma



eeee[®]

DENTAL COMPONENTES



Exports to Germany since 2011

Aperfeiçoando tecnologias



CONEXÕES EXTERNAS

S.M.P. EFF	Aplicação	Implante	Plataforma Protética	Página
A1.1	NEO	3.3/3.5	3.3	13
A1.1	STERI-OSS®	3.3/3.5	3.3	13
A1.1	LIFECORE®	3.3/3.5	3.3	13
A1.1	RESTORE®	3.3/3.5	3.3	13
A1.1	BIONOVATION®	3.3/3.5	3.3	13
A1.1	IMPLAC®	3.3/3.5	3.3	13
A1.1	SERSON®	3.3/3.5	3.3	13
A1.2	NEO	3.75/4.0	4.1	9
A1.3	NEO	4.3	4.3	14
A1.3	BIOHORIZONS®	4.3	4.3	14
A1.4	NEO	4.5/5.0	5.0	15
A1.4	RESTORE®	4.5/5.0	5.0	15
A1.4	GLOBTECK®	4.5/5.0	5.0	15
A1.4	INP®	4.5/5.0	5.0	15
A1.4	INTRA-LOCK®	4.5/5.0	5.0	15
A1.4	IMPLAC®	4.5/5.0	5.0	15
A1.4	DENTOFLEX®	4.5/5.0	5.0	15
A1.4	DIAMOND®	4.5/5.0	5.0	15
B1.1	SIN®	3.25/4.5	3.4	11
B1.2	SIN®	3.25/3.75/4.0	4.1	9
B1.3	SIN®	5.0	5.0	12
C1.1	CONEXÃO®	3.3/3.5	3.3	13
C1.2	CONEXÃO®	3.75/4.0	4.1	9
C1.3	CONEXÃO®	5.0	5.0	15
K1.1	NOBEL BIOCARE®	3.3	3.5	8
K1.1	IMPLACIL®	3.3	3.5	8
K1.1	IMPLAMED®	3.3	3.5	8
K1.1	SWEED®	3.3	3.5	8
K1.1	GLOBTEK®	3.3	3.5	8
K1.1	ACE®	3.3	3.5	8
K1.2	NOBEL BIOCARE®	3.75/4.0	4.1	9
K1.2	IMPLACIL®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	DENTOFLEX®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	INP®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	STERI-OSS®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	GLOBTEK®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	IMPLAMED®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	TITANIUMFIX®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	RESTORE®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	INTRA-LOCK®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	SERSON®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	IMPLAC®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	DIAMOND®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	BIOHORIZONS®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.2	BIOFIX®	3.25/3.5/3.75/3.8/4.0/4.1	4.1	9
K1.3	NOBEL BIOCARE®	5.0/6.0	5.1	10
k1.3	IMPLACIL®	5.0	5.1	10
K1.3	IMPLAMED®	5.0	5.1	10
K1.3	GLOBTEK®	5.0	5.1	10
L1.1	BIOMET 3I™	3.25	3.4	11
L1.2	BIOMET 3I™	3.25mini/3.75/4.0	4.1	9
L1.3	BIOMET 3I™	5.0	5.0	12
L1.3	ACE®	4.75	5.0	12
L1.3	DENTOFLEX®	5.0	5.0	12

*Os produtos indicados com ® são nomes de marcas registradas dos seus respectivos fabricantes.

CONEXÕES INTERNAS

S.M.P. EFF	Aplicação	Implante	Plataforma Protética	Página
A2.1	NEO	3.75/4.3/5.0	4.3/5.0	17
B2.1	SIN®	3.8	3.8	11
B2.2	SIN®	4.0/4.3	4.5	12
B2.3	SIN®	4.0/5.0	5.5	13
C2.1/C2.2/C2.3	CONEXÃO®	3.75/4.0/4.3/5.0	4.1/4.3/5.0/6.0	18
C2.1/C2.2/C2.3	IMPLAC®	3.75/4.0/4.3/5.0	4.1/4.3/5.0/6.0	18
D2.1	IMPLACIL®	3.3/3.5	3.5	19
D2.1	DIAMOND®	3.3/3.5	3.5	19
D2.2	IMPLACIL®	3.75/4.3	4.0	20
D2.2	DIAMOND®	3.75/4.3	4.0	20
D2.3	IMPLACIL®	4.75/5.0	5.0	21
D2.3	DIAMOND®	4.75/5.0	5.0	21
G2.3	DENTOFLEX®	3.75	3.75	22
K2.1	NOBEL BIO CARE® NP	3.5	3.5	8
K2.1	DERIG®	3.5	3.5	8
K2.2	NOBEL BIO CARE® RP	4.3	4.3	9
K2.2	DERIG®	4.3	4.3	9
K2.3	NOBEL BIO CARE® WP	5.0	5.0	10
L2.1	BIOMET 3i™	3.25	3.4	14
L2.2	BIOMET 3i™	4.0	4.1	15
L2.3	BIOMET 3i™	5.0	5.0	16
O2.2	FRIALIT®	3.8	3.8	11
O2.3	FRIALIT®	4.5	4.5	12
O2.4	FRIALIT®	5.5	5.5	13



CONE MORSE

S.M.P. EFF	Aplicação	Implante	Plataforma Protética	Página
A3.1	NEO CM	3.5/3.75/4.0/4.3/5.0	3.5/3.75/4.0/4.3/5.0	20 e 21
A3.1	SIN UNITITE®	3.5/3.75/4.0/4.3/5.0	4.8	20 e 21
A3.1	IMPLACIL®	3.5/3.75/4.0/4.3/5.0	3.5	20 e 21
A3.1	INTRAOSS®	3.5/3.75/4.0/4.3/5.0	3.5/3.75/4.0/4.3/5.0	20 e 21
A3.1	PROSS®	3.5/3.75/4.0/4.3/5.0	3.3/4.5	20 e 21
A3.1	DENTOFLEX®	3.5/3.75/4.0/4.3/5.0	3.3/4.5	20 e 21
A3.1	KOPP®	3.5/3.75/4.0/4.3/5.0	3.3/4.5	20 e 21
A3.1	INP®	3.5/3.75/4.0/4.3/5.0	3.3/4.5	20 e 21
A3.2	NEO EGM	3.5/3.75/4.0/4.3/5.0	3.0	22 e 23
C3.1	SIN STRONG®	3.5/3.8/4.5/5.5	3.5/3.8/4.5/5.5	29
C3.2	CONEXÃO®	3.5/3.75/4.0/4.3/5.0	3.5	9
J3.1	STR RN	3.3/4.1/4.8	4.8	11
J3.1	ZIMMER® Swiss Plus	3.3/4.1/4.8	3.5	11
J3.1	LIFECOR®	3.3/4.1/4.8	4.8	11
J3.1	INTRA-LOCK®	3.3/4.1/4.8	4.8	11
J3.1	GLOBTECK®	3.3/4.1/4.8	4.8	11
J3.1	TITANIUM FIX®	3.3/4.1/4.8	4.8	11
J3.1	BIOMET 3i™	3.3/4.1/4.8	4.8	11
J3.2	STR WN	4.8	6.5	12
J3.5	STR SC	2.9	2.9	14
J3.3	STR NC	3.3	3.3	15
J3.3	Dérig® CMI	3.3	3.3	15
J3.4	STR RC	4.1/4.8	4.1	16
J3.6	STR X RB	3.5	3.5/3.75/4.0/4.5/5.0/5.5	17
J3.7	STR X WB	4.5	3.5/3.75/4.0/4.5/5.0/5.5	18
K3.1	NOBEL BIOCARE®	3.0	3.0	8
K3.1	Dérig® CMH 3.0P	3.0	3.0	8
K3.2	NOBEL BIOCARE® NP	3.5	3.5	9
K3.3	NOBEL BIOCARE® RP	4.3/5.0	3.9	10
K3.4	ANKYLOS®	3.5/4.0/5.0	3.8/4.5/5.5	31
Q3.1	SIGNOVINCES® Duocon	3.5/4.0	3.5/4.0	30
Q3.1	TITANIUMFIX® BlackFix	3.5/4.0	3.5/4.0	30

PILAR

S.M.P. EFF	Aplicação	Plat Protética	Especificação	Catálogo EFF	Página
M.P.	TODOS	Ø4.8	MINI PILAR	CM/CE/CI	32/16/23
M.P.L. 1	NOBEL BIOCARE®	Ø5.4	MINI PILAR LARGA	CE	16
M.P.L.2	NEO	Ø5.4	MINI PILAR LARGO (NEO) CE		17
SYRN	STR RN	Ø4.1/4.8	PILAR SynOcta®	CM	13
SYWN	STR WN	Ø4.8/6.5	PILAR SynOcta®	CM	13
MB 3.50	STR NC	Ø3.5	MULTI BASE®	CM	19
MB 4.50	STR RC	Ø4.5	MULTI BASE®	CM	19
SRA 3.50	STR NC	Ø3.5	SRA®	CM	19
SRA 4.60	STR RC	Ø4.6	SRA®	CM	19
EST.	BRANEMARK®	Ø5.4	ESTHETICONE	CM/CE/CI	32/16/23
ESTN.	CONEXÃO®	Ø5.4	ESTHETICONE	CM/CE/CI	32/16/23
PGM	NEO	Ø4.8	GRAN MORSE	CM	24
MPGM	NEO	Ø3.5	GRAN MORSE MICRO	CM	24
MPGT	NEO	Ø4.8	GT PILAR	CM	25
F1.1/F1.2	DSP®	Ø3.0/Ø4.0	FLEXCONE®	CE	19
1NP	TODOS	Ø3.5/Ø4.3	ESTHETICONE	CM/CE/CI	32/18/23