



Accesso Biotech

Accessible Cell Analysis Solutions




POC CD4 Test for AIDS

Introducing:
ACCESSO VITA™

CD4 test for point of care, solving unmet needs in low and mid income Countries

(whole blood samples from HIV+ patients)



- Portable, external rechargeable battery
- No cold chain storage required
- 20 uL whole blood sample
- <2 min hands-on time for each test
- <2 min measurement after staining
- Total assay time: <35 min
- Throughput: up to 20 tests/hour/instrument



<https://www.accessobio.com/cd4>



Exclusively Distributed by:

www.lumiquick.com



LumiQuick
DIAGNOSTICS, INC.

Contact us for further technical information and product availability :
sales@lumiquick.com - Tel.: +1(408)855-0061