

VA National Contract #36E79720D0029  
Available via PPV McKesson



## STERILE GAUZE BANDAGE ROLLS

These Gauze Bandage Rolls are conforming dressings that can serve as primary or secondary dressings. They have six-ply construction with finished edges reducing loose ends and lint. An open weave design provides fast wicking and superior aeration. These Gauze Rolls are absorbent to reduce the risk of maceration, as well as offering cushioning and wound conformability to provide better wound site protection.



### Specifications

- 4.5 inches wide by 4.1 yards in length
- 100% Cotton, 6-ply
- Sterilized by irradiation
- Individually wrapped peel packs
- Not made with natural rubber latex
- Single Use only
- Made in Taiwan



**Part #B140142500PH**  
**1 EA/PK, 150 PK/CS**

Customer Service: 1-888-581-6922  
Email: [jjaquez@goldcoastmedicalsupply.net](mailto:jjaquez@goldcoastmedicalsupply.net)