



# Scissors Bandage Utility

IAS-1803-50

Supplier: [Gold Coast Medical Supply](#)

DAPA# SP0200-15-H-0030

POC: Jaime Jaquez

Email: [jjaquez@goldcoastmedicalsupply.net](mailto:jjaquez@goldcoastmedicalsupply.net)

Phone: (949) 293-6882



- ◆ This product line includes bandage utility scissors that are reusable, not made with natural rubber latex, made with two stainless steel blades, one serrated blade and one smooth blade, each blade with an angle and moving on a pivot. They are non-sterile multipurpose instrument used for cutting materials such as dressings, casts, and surgical drapes. Rings for fingers and thumb are made from a polypropylene or plastic variant. One plastic finger ring may be larger than the other for maximized control while cutting.
- ◆ Displayed line item and associated image exemplify this product line; contact above POC for additional information.

| Mfr Part #              | Item Description                             | Cardinal PVON | O&M PVON    |
|-------------------------|----------------------------------------------|---------------|-------------|
| <a href="#">GC111SC</a> | STAINLESS STEEL BLACK BANDAGE SCISSORS 7.5IN | GC111SC       | 2850GC111SC |

**Medical Materiel Enterprise Standardization Offices**