

CONTOURWAVE




AXICON
TECHNOLOGIES

PRINCIPLE

- Immerse yourself in the pinnacle of rejuvenation with cutting-edge micro-vibration compression, achieved through a sophisticated 360° rotation of the inner ball within the roller.
- Experience the refined sensation of "pulsation compression," as the rotating ball precisely applies pressure to the skin, creating a continuous push-pull dynamic for unparalleled results.
- Revitalize cellular activity, enhance blood flow, and elevate oxygen levels with expertly targeted pressure and lifting techniques applied to the skin's tissues.
- Effortlessly diminish the appearance of cellulite by expertly pressurizing and breaking down fat deposits, leading to their refined decomposition.
- Attain a smoother, more luxurious skin texture through this advanced wellness technology.





FUNCTION

- Transforms skin texture, significantly reducing the visibility of orange peel-like irregularities.
- Harmonizes the lymphatic system, promoting superior detoxification.
- Enhances skin smoothness, revealing a more luminous and refined complexion.
- Optimizes microcirculation, supporting overall skin health and vibrancy.
- Provides deep muscle relaxation, effectively relieving tension and stress.
- Alleviates discomfort, ensuring a truly luxurious and soothing experience.

Each treatment head is meticulously crafted to address different stages of orange peel cellulite, delivering a customized sensory experience:

- For dense or stubborn cellulite, a dual approach combining mechanical rotation followed by mechanical vibration is recommended for maximum effectiveness.
- For milder cellulite conditions, direct vibration offers a highly effective treatment.



TREATMENT RANGE

BODY

- Legs, Arms, Buttocks, Back, Abdomen, Feet
- Boosts local blood circulation
- Enhances lymphatic drainage
- Eases minor muscle aches and pains
- Accelerates muscle recovery post-exercise



FACE

- Mouth, Eyes, Forehead, Cheeks, Cheekbone, Glabella, Neck, and Lower Neck Area
- Stimulates collagen production
- Enhances vascularization and oxygenation
- Relaxes contracted muscles, smoothing expression lines
- Reduces the appearance of wrinkles




DESIRED OUTCOMES

- Enhanced oxygenation and pain relief.
- Increased blood circulation, resulting in better tissue nourishment.
- Toxin removal, promoting efficient lymphatic drainage.
- Smoother, more defined contours with diminished wrinkles.
- Toned muscles and tightened tissues for enhanced firmness.



HANDPIECES



Mechanical Vibration Handpiece



Mechanical Rotation Handpiece

SPECIFICATIONS

SPEED LEVEL	10 - 100%
ROTATING SPEED	30 - 390mp
RED LIGHT WAVELENGTH	650 - 660nm
MOTOR	Motor with constant speed under load 10Kg
VOLTAGE	220V±10% 50Hz or 110V±10% 60Hz
POWER	135W
MECHANICAL VIBRATION HANDLE SIZE	480.5*108.5mm
MECHANICAL ROTATION HANDLE SIZE	350*57mm



Tomorrow's Technologies. Today.

www.axicontechnologies.com

© 2024 Axicon Technologies LLC - USA
All rights reserved.

