



SIGVARIS  
GROUP

# Medical Compression Therapy Solutions



2024

sigvaris

**mobilis**  
by SIGVARIS GROUP

# Our mission

Established in 1864, SIGVARIS GROUP is, and has always been, a family-owned company. Founded and headquartered in Switzerland, our mission for helping people feel their best is why we've been a global market leader for more than 155 years.

Our compression portfolio includes Medical, Consumer Health Care, products for Medically Complex Edema, and Soft Orthopedic Goods that are specifically designed to restore, assist and enhance people's lives. We are passionate about making exceptional products that combine advanced technology and quality craftsmanship with diverse styles that fit our modern customers' lives.

When heritage, competence, and dedication come together, there's no limit to where we'll go.

Live your life. Be yourself.

For more information, please visit:  
**[sigvaris.com](https://sigvaris.com)**

WE'VE GOT YOU COVERED  
FOR **ALL YOUR COMPRESSION**  
**THERAPY NEEDS!**

**We produce the best compression garments to enhance our customers' health and quality of life.**

We are dedicated to helping legs look and feel their best.



**RECIPROCATED  
TOE PRODUCTS**

Crafted with a smooth, flat toe seam and contoured toe box for added comfort.



**THE MOST COMFORTABLE  
FABRIC OPTIONS AVAILABLE**

Including the finest cottons, wools, microfibers and double covered yarns.



**WE ARE DEDICATED  
TO INNOVATION**

We continue to invest in medical research to help people feel their best. Every day.



**WE OFFER  
MORE VARIETY**

More materials, colours, patterns, and styles than any other medical compression brand.



**THE LARGEST SELECTION OF  
PREMIUM NATURAL FIBERS**

Including Sea Island Cotton, Supima Cotton, Linen and Merino Wool.



**THE HIGHEST  
QUALITY STANDARDS**

More than 20 people handle our products from start to finish, ensuring every garment is of the highest quality.



**DESIGNED FOR  
YOUR BEST FIT**

Using studies on human data measurements, we fit 99.49% of the population.

**WHY  
SIGVARIS  
GROUP?**



# Customer Care

## DEALER PORTAL

The Sigvaris Group Dealer Portal is your online location for securely ordering products and managing your account.

**24/7 access – order anytime!**

View order history  
and tracking

Live inventory status

Sizing Wizard

SKU level comments  
for customer names

Obtain invoice copies

UPC code for a product

Shop the catalogue

Quick order entry

Want to activate your account? Email: [canadaorders@sigvaris.com](mailto:canadaorders@sigvaris.com)



## ONLINE CHAT

Found on dealer portal

**Real time convenience with a live agent**

**Instant responses**

**Fast customer inquiry resolution**

## To Order:

Online Dealer Portal: [Sigvaris-online.com](http://Sigvaris-online.com)

Phone: **1.800.363.4999**

Email: [canadaorders@sigvaris.com](mailto:canadaorders@sigvaris.com)



# Become a Certified Fitter

Elevate your knowledge with our **Sigvaris Online Fitter Certification Course**! Dive into the world of compression wear and become a certified expert in the art of comfort and wellness.



Upon completion you will receive:

- ✓ Certificate
- ✓ Fitter pin
- ✓ Complimentary pair of compression socks

**Questions?** Email [canada.training@sigvaris.com](mailto:canada.training@sigvaris.com)



		PRODUCT	STYLES	COMPRESSION LEVEL	PAGE
MEDICAL COMPRESSION THERAPY SOLUTIONS	ESSENTIAL	<b>COTTON</b> 230	<ul style="list-style-type: none"> <li>● calf, calf w/grip-top, thigh w/grip-top</li> <li>○ calf, thigh w/grip-top</li> </ul>	20–30 mmHg 30–40 mmHg	<b>8</b>
		<b>OPAQUE</b> 860	<ul style="list-style-type: none"> <li>● calf, petite calf, calf plus, calf w/grip-top, thigh w/grip-top, petite thigh w/grip-top, extra-large thigh w/grip-top, pantyhose, petite pantyhose, plus-sized pantyhose</li> <li>○ calf, calf plus, calf w/grip-top, thigh w/grip-top, extra-large thigh w/grip-top, thigh w/waist attachment</li> </ul>	20–30 mmHg 30–40 mmHg	<b>10, 12</b>
	STYLE	<b>SHEER</b> 780 <i>(for women only)</i>	<ul style="list-style-type: none"> <li>● calf, thigh w/grip-top, pantyhose,</li> <li>○ calf, thigh w/grip-top, pantyhose</li> </ul>	15–20 mmHg 20–30 mmHg 30–40 mmHg	<b>16</b>
		<b>MEDIUM SHEER</b> 750 <i>(for women only)</i>	● calf, thigh w/grip-top, pantyhose	20–30 mmHg	<b>18</b>
		<b>SOFT OPAQUE</b> 840 <i>(for women only)</i>	<ul style="list-style-type: none"> <li>● calf, calf w/grip-top, thigh w/grip-top, pantyhose, maternity pantyhose</li> <li>○ calf, thigh w/grip-top, pantyhose</li> </ul>	15–20 mmHg 20–30 mmHg 30–40 mmHg	<b>20</b>
		<b>SEA ISLAND COTTON</b> 220	● calf	20–30 mmHg	<b>22</b>
		<b>MERINO WOOL</b> 240	● calf	20–30 mmHg	<b>24</b>
		<b>LINEN</b> 250	● calf	20–30 mmHg	<b>26</b>
		<b>MICROFIBER</b> 820 <i>(for men only)</i>	● calf, calf w/grip-top, thigh w/grip-top (extended foot sizes available)	15–20 mmHg 20–30 mmHg 30–40 mmHg	<b>28</b>
		<b>MICROFIBER PATTERNS</b> 830	● calf	20–30 mmHg	<b>30</b>
	MOTION	<b>CUSHIONED COTTON</b> 360	● calf	20–30 mmHg	<b>34</b>
		<b>HIGH TECH</b> 412C	● calf	20–30 mmHg	<b>36</b>
		<b>THERMOREGULATING WOOL</b> 422	● calf	20–30 mmHg	<b>38</b>
	SPECIALTY	<b>NATURAL RUBBER</b> 500	○ calf, thigh, thigh w/grip-top, thigh w/waist attachment, pantyhose	30–40 mmHg 40–50 mmHg 50–60 mmHg	<b>42</b>
		<b>ADVANCE ARMSLEEVE</b> 910	○ armsleeve w/grip-top w/gauntlet and without gauntlet	15–20 mmHg 20–30 mmHg	<b>44</b>
		<b>INELASTIC WRAPS</b>			<b>46</b>

Closed Toe = ●      Open Toe = ○

	PRODUCT	STYLES	COMPRESSION LEVEL	PAGE
CONSUMER HEALTH	<b>CASUAL COTTON</b> 146/186	● calf	15–20 mmHg	<b>66</b>
	<b>TRAVENO TRAVEL SOCKS</b> 146/T - 186/T	● calf	15–20 mmHg	<b>67</b>
	<b>SEA ISLAND COTTON</b> 151/191	● calf	15–20 mmHg	<b>68</b>
	<b>ALL-SEASON MERINO WOOL</b> 152/192	● calf	15–20 mmHg	<b>70</b>
	<b>LINEN</b> 153/193	● calf	15–20 mmHg	<b>71</b>
	<b>MICROFIBER SHADES</b> 143/183	● calf	15–20 mmHg	<b>72</b>
	<b>SHEER FASHION</b> 120	● calf, thigh w/grip-top, pantyhose ○ calf	15–20 mmHg	<b>74</b>
	<b>SOFT SILHOUETTE LEGGINGS</b> 170	○ open foot	15–20 mmHg	<b>76</b>
	<b>ATHLETIC RECOVERY SOCKS</b> 401	● calf	15–20 mmHg	<b>78</b>
	<b>PERFORMANCE SLEEVES</b> 412V	○ open foot	20–30 mmHg	<b>79</b>
	<b>MERINO OUTDOOR SOCKS</b> 421	● calf	15–20 mmHg	<b>80</b>
	<b>EVERSOFT DIABETIC SOCKS</b> 160	● calf	8–15 mmHg	<b>81</b>
	<b>ACCESSORIES</b>			<b>82</b>
RESOURCES	<b>SYMPTOMS &amp; INDICATIONS</b>			<b>90</b>
	<b>GLOSSARIES OF ICONS &amp; TERMS</b>			<b>91-92</b>
	<b>HOW TO MEASURE &amp; ORDER</b>			<b>93-95</b>
MOBILIS	<b>mobilis</b> ORTHOPEDIC BRACING			<b>96-117</b>
	<b>ELBOW</b>			<b>100</b>
	<b>WRIST &amp; THUMB</b>			<b>103</b>
	<b>BACK</b>			<b>107</b>
	<b>KNEE</b>			<b>110</b>
	<b>ANKLE</b>			<b>114</b>
	<b>HOW TO ORDER MOBILIS</b>			<b>117</b>



**sigvaris** ESSENTIAL

**Medical Compression Therapy Solutions**





## Support Yourself

### Why choose **Essential**?

- Reduce swelling and discomfort
- Long-distance travel
- Everyday comfort and reliability

#### ESSENTIAL PRODUCTS

Series	Product	Page
<b>230</b>	Cotton	8
<b>860</b>	Opaque (Women)	10
<b>860</b>	Opaque (Men)	12

#### BEST SUITED FOR

- Containment of various types of edema, including venous edema, phlebolymphe<sup>d</sup>ema, and early stages of lymph<sup>e</sup>dema
- Ideal for dress or casual shoes
- Suited for post-vein treatment

Featured  
colour is  
Light Beige

## COTTON 230

for women & men

### Compression levels

20–30 mmHg

30–40 mmHg

### Features and Benefits



Sensitive Skin



Cotton



Breathable

- 100% Supima Cotton against the skin
- Ideal for sensitive skin and allergies
- Breathable and durable
- Excellent containment of various types of edema
- Classic ribbed design for everyday wear

### Fiber Content

- 50% Nylon, 25% Spandex, 25% Supima Cotton (content can vary slightly by style)
- Latex free



**20–30 mmHg**  
Closed Toe Calf w/ Grip-top  
Medium Long in Black

**232CMLM99/S**



Silicone dotted grip-top available  
for calf and thigh



**Click here to see product images online**  
[sigvaris.com](https://sigvaris.com)



## SIZING CHART






	Small		Medium		Large		X Large	
<b>Ankle Circumference</b>	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
<b>Calf Circumference</b>	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
<b>Calf Length</b>	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
<b>Thigh Circumference</b>	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		NA	
<b>Leg Length: Women &amp; Open Toe</b>	Short	Long	Short	Long	Short	Long		
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")		
<b>Leg Length: Men</b>	Short	Long	Short	Long	Short	Long		
	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")		
<b>Size</b>	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

## STYLES AND COLOURS

Closed Toe = ●  
Open Toe = ○

				
WHITE 00	NAVY 10	BLACK MIST 14	LIGHT BEIGE 66	BLACK 99

<b>Calf (C)*</b>	●	●	●	●○	●
<b>Calf w/Grip-Top (C/S)</b> Silicone dotted grip-top band				●**	●
<b>Thigh w/Grip-Top (N)</b> Silicone dotted grip-top band				●○	●

\* 30–40 mmHg extra large calf only available in black with grip-top

\*\* 20–30 mmHg only

## Product Code Example: 232NMLW66

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
23	2	N	ML	W	66

## Product Code Example: 232CLSM66/S

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour	Grip-top
23	2	C	LS	M	66	/S

Featured  
colour is  
Honey

## OPAQUE 860

for women

### Compression Levels

20–30 mmHg

30–40 mmHg

### Features and Benefits



Durable



Largest Size Range



Opaque

- Largest size range, from petite to plus
- Durable construction combined with a smooth opaque look
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema
- Plus size pantyhose are an alternative to maternity pantyhose

### Fiber Content

- 64% Nylon, 36% Spandex  
(content can vary slightly by style)
- Latex free



**20–30 mmHg**  
Closed Toe Calf Medium Long  
in Black

**862CMLW99**



Silicone dotted grip-top available  
for calf and thigh



**Click here to see product images online**  
[sigvaris.com](http://sigvaris.com)

## SIZING CHART

	Small			Medium			Large		X Large	
<b>Ankle Circumference</b>	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
<b>Calf Circumference</b>	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")		41–61 cm (16"–24")	
<b>Calf Circumference: Calf Plus w/Grip-Top</b>	NA			NA	39.5–53 cm (15.5"–21")		45.5–61 cm (18"–24")		51–71 cm (20"–28")	
<b>Calf Length</b>	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤35.5 cm (14")	≤41 cm (16")	>41 cm (16")	≤35.5 cm (14")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
<b>Thigh Circumference</b>	45.5–66 cm (18"–26")			48–73 cm (19"–28.5")			52–84 cm (20.5"–33")		64–95 cm (25"–37.5")	
<b>Leg Length</b>	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤65 cm (25.5")	≤72 cm (28.5")	>72 cm (28.5")	≤65 cm (25.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
<b>Hip Circumference</b>	83–112.5 cm (32.5"–44.5")			87–127 cm (34.5"–50")			98–150.5 cm (38.5"–59.5")		NA	
<b>Plus or Maternity Size Hip Circumference</b>	NA	≤142 cm (56")		NA	≤157.5 cm (62")		≤178 cm (70")		NA	
<b>Size</b>	SP	SS	SL	MP	MS	ML	LS	LL	XS	XL
<b>Calf Plus w/Grip-Top</b>					MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Women's Shoe Size : Petite Leg 4–8, Short Leg 5.5–9.5, Long Leg 8–12

## STYLES AND COLOURS

Closed Toe = ●  
Open Toe = ○

	HONEY 33	GOLDEN 36	LIGHT BEIGE 66	BLACK 99
<b>Calf (C)</b>	●	●	○	●○
<b>Calf Petite Sizes (C)*</b> Not available with dotted grip-top band	●	●		●
<b>Calf w/Grip-Top (C/S)</b> Silicone dotted grip-top band			●○	●
<b>Calf Plus w/Grip-Top (C/P)</b> Only available with silicone dotted grip-top band			●○	●○*
<b>Thigh w/Grip-Top (N)</b> Silicone dotted grip-top band	●	●	○	●○
<b>Thigh w/Grip-Top (N) XS &amp; XL</b> Silicone dotted grip-top band			●○	○
<b>Thigh Petite Sizes (N)*</b>	●	●		●
<b>Pantyhose (P)**</b>	●	●		●
<b>Pantyhose Petite (P)*</b>	●	●		●
<b>Pantyhose Plus and Maternity (P/P)**</b>			●	●
<b>Thigh w/Waist Attachment (W)**</b> Indicate leg LEFT (/L) or RIGHT (/R)			○	

\* 20-30 mmHg only

\*\* Not available in XS or XL sizing

## Product Code Example: 863NMSW33

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
86	3	N	MS	W	33

## Product Code Example: 862CLLO66/P

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Plus size with Grip-top
86	2	C	LL	O	66	/P



Featured  
colour is  
Ribbed Black



## OPAQUE 860

for men

### Compression Levels

20–30 mmHg

30–40 mmHg

### Features and Benefits



Durable



Largest Size Range



Opaque

- Durable construction with a ribbed design on closed toe calf products
- Smooth opaque design on open-toe styles
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema

### Fiber Content

- 64% Nylon, 36% Spandex  
(content can vary slightly by style)
- Latex free



**20–30 mmHg**  
Open Toe Thigh Medium Long  
in Light Beige

**862NMLO66**



Silicone dotted grip-top available  
for calf and thigh



**Click here to see product images online**  
[sigvaris.com](https://sigvaris.com)





## SIZING CHART

	Small		Medium		Large		X Large	
<b>Ankle Circumference</b>	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
<b>Calf Circumference</b>	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
<b>Calf Circumference: Calf Plus w/Grip-Top</b>	NA		39.5–53 cm (15.5"–21")		45.5–61 cm (18"–24")		51–71 cm (20"–28")	
<b>Calf Length</b>	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
<b>Thigh Circumference</b>	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		64–95 cm (25"–37.5")	
<b>Leg Length</b>	Short	Long	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
<b>Size</b>	SS	SL	MS	ML	LS	LL	XS	XL
<b>Calf Plus w/Grip-Top</b>			MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

## STYLES AND COLOURS

Closed Toe = ●  
Open Toe = ○



<b>Calf (C)</b>	●*	●*	○	○
<b>Calf w/Grip-Top (C/S)</b> Silicone dotted grip-top band			○	
<b>Calf Plus w/Grip-Top (C/P)</b> Only available with silicone dotted grip-top band			○	○*
<b>Open Toe Thigh w/Grip-Top (N)</b> Silicone dotted grip-top band			○	○
<b>Thigh w/Waist Attachment (W)**</b> Indicate leg LEFT (/L) or RIGHT (/R)			○	

\* 20–30 mmHg only

\*\* Not available in XS or XL with waist attachment

## Product Code Example: 862CMLM99

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
86	2	C	ML	M	99

## Product Code Example: 862NMLO66

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
86	2	N	ML	O	66



## Medical Compression Therapy Solutions





## Express Yourself

Why choose **Style**?

- Fashionable, modern patterns and designs
- Special occasion and formal wear
- Our most luxurious feel

### STYLE PRODUCTS

Series	Product	Page
<b>780</b>	Sheer	16
<b>750</b>	Medium Sheer	18
<b>840</b>	Soft Opaque	20
<b>220</b>	Sea Island Cotton	22
<b>240</b>	Merino Wool	24
<b>250</b>	Linen	26
<b>820</b>	Microfiber	28
<b>830</b>	Microfiber Patterns	30

### BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment with minimal edema
- Pregnancy
- People who prefer fashionable options
- Ideal for business or dress attire



Featured colour is  
Toasted  
Almond

## SHEER 780

for women

### Compression Levels

15–20 mmHg

20–30 mmHg

30–40 mmHg

### Features and Benefits



Sheer



Elegant



Breathable

- Fashionable for professional and dress wear
- Consistent sheerness throughout the leg
- Excellent breathability and moisture management

### Fiber Content

- 15–20 and 20–30 mmHg: 69% Nylon, 31% Spandex (content can vary slightly by style)
- 30–40 mmHg: 60% Nylon, 40% Spandex (content can vary slightly by style)
- Latex free



**15–20 mmHg**  
Closed Toe Thigh Medium Long  
in Dark Navy

**781NMLW08**



Sensinno® grip-top on thigh



**Click here to see product images online**  
[sigvaris.com](https://sigvaris.com)

## SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

## STYLES AND COLOURS

Closed Toe = ● Open Toe = ○	DARK NAVY 08	HONEY 33	WARM SAND 34*	GOLDEN 36	TOASTED ALMOND 37*	CAFÉ 73*	MOCHA 85*	NIGHT-SHADE 94	BLACK 99
	●	●○	●○	●○	●○	●○	●○	●	●○
Calf (C)	●	●○	●○	●○	●○	●○	●○	●	●○
Thigh w/Grip-Top (N) Sensinnov® grip-top band	●	●○	●○	●○	●○	●○	●○	●	●○
Pantyhose (P)	●	●○	●○	●○	●○	●○	●○	●	●○

\* 20–30 mmHg only


## Product Code Example: 781NMLW36

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
78	1	N	ML	W	36

## Product Code Example: 783PSSO33

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
78	3	P	SS	O	33

Featured  
colour is  
Black



## MEDIUM SHEER 750

for women

### Compression Level

20–30 mmHg

### Features and Benefits



Medium Sheer



Durable



Breathable

- Fashionable, durable, medium sheer hosiery for everyday wear
- Excellent breathability and moisture management
- Comfortable 2" top band on calf style

### Fiber Content

- 68% Nylon, 32% Spandex  
(content can vary slightly by style)
- Latex free



**20–30 mmHg**  
Closed Toe Thigh Medium Long  
in Honey

**752NMLW33**



Silicone dotted grip-top on thigh



**Click here to see product images online**  
[sigvaris.com](https://sigvaris.com)



## SIZING CHART

	Small		Medium		Large	
<b>Ankle Circumference</b>	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
<b>Calf Circumference</b>	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
<b>Calf Length</b>	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
<b>Thigh Circumference</b>	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
<b>Leg Length</b>	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
<b>Hip Circumference</b>	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
<b>Size</b>	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

## STYLES AND COLOURS

Closed Toe = ●			
<b>Calf (C)</b>	●	●	●
<b>Thigh w/Grip-Top (N)</b> Silicone dotted grip-top band	●	●	●
<b>Pantyhose (P)</b>	●	●	●

### Product Code Example: 752NMLW33

Series	mmHg Level	Style	Size	Gender	Colour
75	2	N	ML	W	33



Featured colour  
is Mulberry

# SOFT OPAQUE 840

for women

## Compression Levels

15–20 mmHg

20–30 mmHg

30–40 mmHg

## Features and Benefits



**Opaque**



**Microfiber**



**Soft**

- Fashionable opaque hosiery for professional and everyday wear
- Long-lasting microfiber softness provides all day comfort and allows product to easily glide on and off
- Durable and breathable
- Maternity pantyhose option available

## Fiber Content

- 79% Nylon, 21% Spandex  
(content can vary slightly by style)
- Latex free



**20–30 mmHg**  
Closed Toe Pantyhose Medium  
Long in Black

**842PMLW99**



Silicone dotted grip-top for calf  
and SensinnoV® grip-top for thigh



**Click here to see product images online**  
[sigvaris.com](http://sigvaris.com)

## SIZING CHART

	Small		Medium		Large	
<b>Ankle Circumference</b>	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
<b>Calf Circumference</b>	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
<b>Calf Length</b>	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
<b>Thigh Circumference</b>	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
<b>Leg Length</b>	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
<b>Hip Circumference</b>	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
<b>Maternity hip circumference</b>	≤142 cm (56")		≤157.5 cm (62")		≤178 cm (70")	
<b>Size</b>	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

## STYLES AND COLOURS

Closed Toe = ●  
Open Toe = ○

	MIDNIGHT BLUE 09	CHAI 35	PECAN 40	MULBERRY 60*	ESPRESSO 89*	GRAPHITE 91	BLACK 99
<b>Calf (C)</b>	●	●○	●	●	●	●	●○
<b>Calf w/Grip-Top (C/S)*</b> Silicone dotted grip-top							●
<b>Thigh w/Grip-Top (N)</b> Sensinnov® grip-top band	●	●○	●	●	●	●	●○
<b>Pantyhose (P)</b>	●	●○	●	●	●	●	●○
<b>Maternity Pantyhose (M)*</b>							●

\* 20–30 mmHg only

## Product Code Example: 842CMSW60

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
84	2	C	MS	W	60





Featured colour  
is Black

## SEA ISLAND COTTON 220

for women & men

### Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 68

### Features and Benefits



Soft



Cotton



Breathable

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

### Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free



**20–30 mmHg**  
Closed Toe Calf Small Long  
in Brown

**222CSLW11**



**20–30 mmHg**  
Closed Toe Calf Medium Long  
in Navy

**222CMLM10**



**Click here to see product images online**  
[sigvaris.com](https://sigvaris.com)

## SIZING CHART


	Small		Medium		Large		X Large	
<b>Ankle Circumference</b>	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
<b>Calf Circumference</b>	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
<b>Calf Length</b>	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
<b>Size</b>	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

## STYLES AND COLOURS

Closed Toe = ●



<b>Calf (C)</b>	●	●	●
-----------------	---	---	---

## Product Code Example: 222CLLM10

Series	mmHg Level	Style	Size	Gender	Colour
22	2	C	LL	M	10



Featured colour  
is Black

## MERINO WOOL 240

for women & men

### Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 70

### Features and Benefits



Wool



Thermoregulating



Stylish

- High quality extra-fine Merino wool makes this sock soft and breathable
- Moisture wicking capabilities keep feet dry
- Naturally odour resistant
- Thermoregulating for year-round wear
- Smooth, flat toe seam

### Fiber Content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free



20–30 mmHg  
Closed Toe Calf Medium Long  
in Brown

**242CMLM11**



20–30 mmHg  
Closed Toe Calf Medium Long  
in Navy

**242CMLW10**



*Click here to see product images online*  
[sigvaris.com](https://sigvaris.com)



## SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

## STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---

## Product Code Example: 242CMSW99

Series	mmHg Level	Style	Size	Gender	Colour
24	2	C	MS	W	99

Featured  
colour is  
Denim

## LINEN 250

for women & men

### Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 71

### Features and Benefits



**Hypoallergenic**  
Natural fibers



**Anti-bacterial**



**Moisture  
Management**

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

### Fiber Content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free



**20–30 mmHg**  
Closed Toe Calf Medium Long  
in Light Grey

**252CMLM205**



**20–30 mmHg**  
Closed Toe Calf Medium Long  
in Lavender

**252CMLM207**



*Click here to see product images online*  
[sigvaris.com](https://sigvaris.com)

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12  
Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●	●	●
----------	---	---	---	---

Product Code Example: 252CMLM207

Series	mmHg Level	Style	Size	Gender	Colour
25	2	C	ML	M	207



# MICROFIBER 820

for men

## Compression Levels

15-20 mmHg

20-30 mmHg

30-40 mmHg

## Features and Benefits



Microfiber



Soft



Stylish

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking
- Glides on and off with ease
- Comfortable and durable sock for professional and everyday wear
- Available in extended foot sizes
- Smooth, flat toe seam

## Fiber Content


- 87% Nylon , 13% Spandex  
(content can vary slightly by style)
- Latex free

Featured  
colour is  
Tan Khaki



**20-30 mmHg**  
Closed Toe Thigh Medium Long  
in Black

**822NMLM99**



Silicone dotted grip-top option  
available for calf. Sensinno®  
grip-top for thigh.



**Click here to see product images online**  
[sigvaris.com](https://sigvaris.com)

## SIZING CHART

	Small		Medium		Large		X Large	
<b>Ankle Circumference</b>	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
<b>Calf Circumference</b>	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
<b>Calf Length</b>	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
<b>Thigh Circumference</b>	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		NA	
<b>Leg Length</b>	Short	Long	Short	Long	Short	Long		
	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")		
<b>Size Shoe Size</b>	SS 7–11	SL 9.5–14	MS 7–11	ML 9.5–14	LS 7–11	LL 9.5–14	XS 7–11	XL 9.5–14
				MX* 12.5–17.5		LX* 12.5–17.5		XX* 12.5–17.5

Men's shoe sizes: Short Leg 7–11, Long Leg 9.5–14 (For extended foot sizes: long leg 12.5–17.5)

\*MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15–20 & 20–30 mmHg only. Extended foot sizes not available with grip top

## STYLES AND COLOURS

Closed Toe = ●



<b>Calf (C)</b>	●	●***	●
<b>Calf w/Grip-Top (C/S)**</b> Silicone dotted grip-top band	●		●
<b>Thigh w/Grip-Top (N)**</b> Sensinnov® grip-top band	●		●

\*\* Not available in 15–20 mmHg

\*\*\* Steel Grey only available in 15–20 mmHg and 20–30 mmHg

## Product Code Example: 823CMLM32

Series	mmHg Level	Style	Size	Gender	Colour
82	3	C	ML	M	32

## Product Code Example: 822CMLM32/S

Series	mmHg Level	Style	Size	Gender	Colour	Grip-top
82	2	C	ML	M	32	/S

# MICROFIBER PATTERNS 830


for women & men

## Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 72

## Features and Benefits



Microfiber



Stylish



Patterns

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking.
- Fashionable, durable sock for everyday wear
- Glides on and off with ease
- Smooth, flat toe seam

## Fiber Content

- 87% Nylon, 13% Spandex  
(content can vary slightly by style)
- Latex free

Featured colour is  
Wisteria Windows



20–30 mmHg  
Closed Toe Calf Small Long  
in Royal Blue Argyle

**832CSLM48**



20–30 mmHg  
Closed Toe Calf Small Short  
in Mariner Stripe

**832CSSW02**



*Click here to see product images online*  
[sigvaris.com](https://sigvaris.com)

## SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

## STYLES AND COLOURS

### Women

Closed Toe = ●




Calf (C)	●	●	●	●	●
----------	---	---	---	---	---



Calf (C)	●	●	●	●	●
----------	---	---	---	---	---

### Men

Closed Toe = ●



Calf (C)	●	●	●	●	●
----------	---	---	---	---	---



Calf (C)	●	●
----------	---	---

### Product Code Example: 832CMLW02

Series	mmHg Level	Style	Size	Gender	Colour
83	2	C	ML	W	02



sigvaris MOTION

Medical Compression Therapy Solutions



## Push Yourself

### Why choose **Motion**?

- Cushioned sole for added comfort
- Durable thermoregulating wool (422)
- Moisture-wicking yarns

#### MOTION PRODUCTS

Series	Product	Page
<b>360</b>	Cushioned Cotton	34
<b>412C</b>	High Tech	36
<b>422</b>	Thermoregulating Wool	38

#### BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- People who stand for prolonged periods of time
- Casual and occupational footwear
- Walking and everyday activities

# CUSHIONED COTTON 360

for women & men

## Compression Level

20–30 mmHg

## Features and Benefits



Cushioned Sole



Cotton



Sensitive Skin

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Ideal for sensitive skin and allergies
- Smooth, flat toe seam
- Diabetic friendly
- Ideal for people who stand for prolonged periods of time

## Fiber Content

- 64% Supima Cotton, 28% Nylon, 8% Spandex
- Latex free



**Click here to see product images online**  
sigvaris.com

Featured  
colour is  
White



**20–30 mmHg**  
Closed Toe Calf Medium Long  
in Black

**362CMLW99**



Flat, low friction toe seam

## SIZING CHART


	Small		Medium		Large		X Large	
<b>Ankle Circumference</b>	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
<b>Calf Circumference</b>	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
<b>Calf Length</b>	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
<b>Size</b>	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

## STYLES AND COLOURS

Closed Toe = ●



<b>Calf (C)</b>	●	●
-----------------	---	---

## Product Code Example: 362CMLM00

Series	mmHg Level	Style	Size	Gender	Colour
36	2	C	ML	M	00





Featured colour is Lime

## HIGH TECH 412C

for women & men

### Compression Level

20–30 mmHg

### Features and Benefits



**Moisture Wicking**



**Anatomically Contoured**



**Colours**

- High-Tech fibers for moisture wicking and improved thermoregulation
- Cushioned footbed, toe, and Achilles tendon zones for optimal comfort
- Anatomically contoured to support the natural shape of the foot and calf
- Smooth, flat toe seam
- Perfect for everyday activities and active individuals

### Fiber Content

- 64% Nylon, 17% Polyester, 19% Spandex
- Latex free



**20–30 mmHg**  
Closed Toe Calf Small Medium in Steel Blue

**412CSM62**



**20–30 mmHg**  
Closed Toe Calf Medium Long in Pink

**412CML56**






**Click here to see product images online**  
[sigvaris.com](https://sigvaris.com)

## SIZING CHART

	Small			Medium			Large			X Large		
<b>Ankle Circumference</b>	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")			29–33 cm (11.5"–13")		
<b>Calf Circumference</b>	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")			40.5–58 cm (16"–23")		
<b>Women's Shoe Size</b>	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14	8–10.5	11–13.5	>14
<b>Men's Shoe Size</b>	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5	7–9.5	10–12.5	13–15.5
<b>Size</b>	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

## STYLES AND COLOURS

Closed Toe = ●



<b>Calf (C)</b>	●	●	●	●	●	●	●
-----------------	---	---	---	---	---	---	---

## Product Code Example: 412CML54

Series	mmHg Level	Style	Size	Colour
41	2	C	ML	54



Featured colour is Charcoal

# THERMO-REGULATING WOOL 422

for women & men

## Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 80

## Features and Benefits



**Wool**



**Cushioned Sole**



**Thermoregulating**

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned, padded soles help prevent blisters
- Soft and comfortable
- Smooth, flat toe seam
- Perfect for outdoor work environments and activities

## Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free



**Click here to see product images online**  
sigvaris.com



**20–30 mmHg**  
Closed Toe Calf Large Long  
in Charcoal

**422CLL12**



**20–30 mmHg**  
Closed Toe Calf Medium Short  
in Olive


**422CMS31**

## SIZING CHART

	Small			Medium			Large		
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")		
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")		
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL

## STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●
----------	---	---

## Product Code Example: 422CML12

Series	mmHg Level	Style	Size	Colour
42	2	C	ML	12



**sigvaris** SPECIALTY

**Medical Compression Therapy Solutions**



## Take care of Yourself

### Why choose **Specialty**?

- Specific medical needs
- Need low fatigue fibers
- Natural rubber (500)

#### SPECIALTY PRODUCTS

Series	Product	Page
<b>500</b>	Natural Rubber	42
<b>910</b>	Advance Armsleeve	44

#### BEST SUITED FOR

- Severe varicosities
- Chronic edema post DVT
- Severe venous stasis
- Early stages of lymphedema
- Active venous ulcers and prevention of ulcer recurrence
- People with early stages of arm lymphedema
- Post-mastectomy care





# NATURAL RUBBER 500

for women & men

## Compression Levels

30–40 mmHg

40–50 mmHg

50–60 mmHg


## Features and Benefits



Natural Rubber



Durable



Largest Size Range

- Durable, long-term containment solution with minimal fiber fatigue throughout the day
- Extended sizing for a custom-like fit
- Ideal for patients with more serious indications, including:
  - Severe varicosities
  - Chronic edema post DVT
  - Prevention and treatment of Active Venous Ulcers
  - Severe venous stasis
  - Lymphedema

## Fiber Content

- 70% Nylon, 30% Natural Latex Rubber (Calf & Thigh-high styles)
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibers (Pantyhose style)
- **Caution:** contains natural rubber latex that may cause allergic reaction



**30–40 mmHg**  
Open Toe Thigh Medium Average  
Long in Beige

**503NM2077**



Silicone dotted grip-top available  
on thigh



**Click here to see product images online**  
sigvaris.com

## SIZING CHART

	Small				Medium				Large				X Large			
Ankle Circumference	19–23 cm (7.5"–9")				23–30 cm (9"–12")				30–40 cm (12"–15.5")				40–46 cm (15.5"–18")			
Calf Circumference	Average		Full		Average		Full		Average		Full		Average		Full	
	28–38 cm (11"–15")		≤43 cm (17")		32–43 cm (12.5"–17")		≤48 cm (19")		36–48 cm (14.5"–19")		≤57 cm (22.5")		45–57 cm (18"–22.5")		≤68 cm (27")	
Calf Circumference for 505C (Full only)	505C: full calf only		32–43 cm (12.5"–17")		505C: full calf only		36–48 cm (14.5"–19")		505C: full calf only		45–57 cm (18"–22.5")		505C: not available in XL			
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")
Thigh Circumference	Average		Full		Average		Full		Average		Full		N/A			
	41–54 cm (16"–21")		≤61 cm (24")		48–61 cm (19"–24")		≤71 cm (28")		58–71 cm (23"–28")		≤81 cm (32")					
Thigh Circumference for 505T (Full only)	505T: full thigh only		≤61 cm (24")		505T: full thigh only		≤71 cm (28")		505T: full thigh only		≤81 cm (32")					
Leg Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long				
	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")				
Hip Circumference	<122 cm (48")				<137 cm (54")				<155 cm (61")							
Size	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4

## STYLES AND COLOURS

Open Toe = ○



Caution: this product contains natural rubber latex which may cause allergic reactions. For all indications, cautions, and contraindications, please see page 90 for details. Sigvaris accepts no liability for non-observance of contraindications or cautions.

Calf (C)	○	Calf only available in S, M, and L full circumference for 50–60 mmHg
Thigh w/out Grip-Top (T)	○	
Thigh w/Grip-Top (N) Silicone dotted grip-top band	○	Not available in 50–60 mmHg or XL
Pantyhose (P)	○	Pantyhose only available in 30–40 mmHg and only in full thigh circumference
Thigh w/Waist Attachment (W) Indicate Leg LEFT (/L) or RIGHT (/R)	○	Thigh w/ waist attachment is available in average and full thigh circumferences in 30–40 mmHg and in full thigh circumferences only in 40–50 mmHg

Product Code Example: 504CS2077

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
50	4	C	S2	O	77

Product Code Example: 503WL1077/L

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Leg
50	3	W	L1	O	77	/L





# ADVANCE ARMSLEEVE 910

for women & men

## Compression Levels

15–20 mmHg

20–30 mmHg

## Features and Benefits



Moisture  
Wicking



SensinnoV  
Grip-top Band




Breathable

- Ideal for patients with arm lymphedema or post-mastectomy care
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable and soft
- SensinnoV® grip-top band ensures armsleeve stays in place all day long
- Available with or without gauntlet

## Fiber Content

- 50% Nylon, 28% Rayon, 22% Spandex
- Latex free



With  
Gauntlet

**20–30 mmHg**  
Armsleeve With Gauntlet  
Small Long in Beige

**912ASLG77/S**



SensinnoV® grip-top band

## SIZING CHART

	Small		Medium		Large		X Large	
<b>Palm (A)</b>	17–21 cm (6.5"–8")		18–22 cm (7"–8.5")		19–23 cm (7.5"–9")		20–24 cm (8"–9.5")	
<b>Wrist (C)</b>	14–16 cm (5.5"–6.5")		16–18 cm (6.5"–7")		18–20 cm (7"–8")		20–23 cm (8"–9")	
<b>Upper Mid-Arm (F)</b>	Regular		Regular		Plus		Regular	
	26–34 cm (10"–13.5")		29–37 cm (11.5"–14.5")		37–43 cm (14.5"–17")		32–40 cm (12.5"–15.5")	
<b>Arm Length (C-G)*</b>	Regular	Long	Regular	Long	Regular	Long	Regular	Long
	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")
<b>Size</b>	SR	SL	MR	ML	MR/P	ML/P	LR	LL

\* Should have **slight bend** in arm when measuring length


## STYLES AND COLOURS

Open Sleeve = ●



<b>Armsleeve</b> w/Gauntlet (G) Sensinnov® grip-top band	●
<b>Armsleeve</b> w/out Gauntlet (O) Sensinnov® grip-top band	●

**Note:** all armsleeves come with Grip-top (/S)  
15–20 mmHg not available with gauntlet or plus sizing



## ARMSLEEVE WITHOUT GAUNTLET

**Product Code Example: 912ASLO77/S**

Series	mmHg Level	Style	Size	W/out Gauntlet	Colour	Grip-Top
91	2	A	SL	O	77	/S

## ARMSLEEVE WITH GAUNTLET

**Product Code Example: 912ALRG77/PS**

Series	mmHg Level	Style	Size	W/Gauntlet	Colour	Plus Size	Grip-Top
91	2	A	LR	G	77	/P	S

# sigvaris INELASTIC WRAPS

Ready-to-wear wraps for the active treatment and management of conditions impacting the lymphatic and venous systems



## Why Choose Inelastic compression wraps?

- Adjustable pressure
- Large variety of solutions

### MEDICALLY COMPLEX EDEMA PRODUCTS

Product	Page
Coolflex™	48
Compreflex™	52
Compreflex™ Standard Calf & Foot	54
Compreflex™ Liner	56
Compreboot®	57
Compreflex™ Strap Extender	58
Accutab®	58
Compreflex™ Thigh	59
Compreflex™ Knee	60
Compreknee™	61
Compreshort®	62

*The SIGVARIS GROUP line of products for medically complex edema (MCE) are for the treatment and management of conditions impacting the lymphatic and venous circulatory systems.*

*The product portfolio contains ready-to-wear garments for the active treatment and maintenance of lymphedema and for patients managing chronic venous disorders in the lower extremities. Sigvaris is looking towards a brighter future with new therapies for complex venous and lymphatic diseases. Our inelastic wraps deliver maximum therapeutic results with increased patient compliance.*



## COOLFLEX™

### Wrap Details

#### LIGHTWEIGHT, BREATHABLE COMPRESSION WEAR

Coolflex™ is a lightweight compression device, optimal for individuals with fragile or broken skin. Made of a non-stretch mesh material, the Coolflex™ transfers heat and wicks away moisture, leaving the skin cool and dry.

The features of the closure system, along with the built-in guidelines for application, promote independence and encourage improved patient outcomes.

#### BREATHABLE GARMENT

The Spacer Fabric® TR is designed to be highly breathable for cool-wearing comfort

#### ONE-HANDED DONNING

EZ-Close Pressure System of eight tabs makes one-handed, lateral donning simple, even for patients with dexterity issues

#### INNOVATIVE DESIGN

Durable, metal hook closures for maximum strength and durability



20–50 mmHg



### Ideal for

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Well-managed Lymphedema
- Non-specific or post-surgical edema
- Easy to put on for patients with limited dexterity


### Features

- Lightweight and breathable for cool-wearing comfort
- Durable metal clasps, instead of Velcro closure
- Low profile design allows boot to fit in most footwear
- Spacer Fabric TR is highly breathable for cool-wearing comfort



### ACCURATE COMPRESSION IS EASY WITH THE COOLFLEX™

Each Coolflex™ calf garment includes the unique, EZ-Close™ pressure system. The EZ-Close pressure system is positioned laterally, allowing the patient to apply the product without bending forward. The bands of the closure system are made to an optimal length, making them easy to grip and close using minimal hand strength and little dexterity.



The integrated EZ-Close Pressure System has visual indicators for the accurate delivery of 20–30 mmHg, 30–40 mmHg or 40–50 mmHg on each tab



When pairing with a basic liner, a foot piece is required.



When pairing with a transition liner, a foot piece is not required\*.

\* If a foot piece is not required, please order a transition liner. See page 56.

## Standard Calf

20–50 mmHg



Ideal for both venous and lymphatic limb disorders, the garments are available in standard and max sizing and can accommodate limbs that feature disproportionately small ankles with a larger calf.

### Features

- Lightweight and breathable for cool-wearing comfort
- Easy to put on for patients with limited dexterity
- Durable metal clasps, instead of Velcro closure

### Includes

- Coolflex™ Standard Calf
- 1 Pair of non-compressive basic Liners per leg

### Colour




Black

### COOLFLEX™ STANDARD CALF SIZE CHART (A–C1)

(Available in Black Only)

		Small	Small Max	Medium / Large	Medium / Large Max	X Large / XX Large
	C1	28–39 cm	35–50 cm	34–50 cm	45–64 cm	45–64 cm
	C	28–39 cm	35–50 cm	34–50 cm	45–64 cm	45–64 cm
	B	23–33 cm	27–38 cm	29–43 cm	35–49 cm	37–54 cm
	A	16–26 cm	16–26 cm	23–33 cm	23–33 cm	30–42 cm
Left	Regular	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
	Tall	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L
Right	Regular	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
	Tall	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R



## Standard Foot


Variable compression



Foot compression garments are pictured with the Basic Liner. Basic liner is included with the wrap but not with the footpiece

### Features

- Low profile design allows boot to fit in most footwear
- Spacer Fabric TR is highly breathable for cool-wearing comfort
- Durable metal clasps, instead of Velcro closure



### Foot size measuring Instructions:

- I:** Measure circumference of ankle bend and heel  
**J:** Measure circumference across metatarsal heads  
**K:** Measure **medial** length from heel to 1<sup>st</sup> metatarsal head

## COOLFLEX™ STANDARD FOOT SIZE CHART (I, J, K)

(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–39 cm	N/A	30–43 cm	N/A	34–49 cm	N/A
J	22–26 cm	N/A	25–30 cm	N/A	28–33 cm	N/A
K	14–18 cm	19–21 cm	18–20 cm	20–22 cm	20–22 cm	23–25 cm
Right	2101-BTR-R	2101-BTL-R	2103-BTR-R	2103-BTL-R	2105-BTR-R	2105-BTL-R
Left	2101-BTR-L	2101-BTL-L	2103-BTR-L	2103-BTL-L	2105-BTR-L	2105-BTL-L



## COMPREFLEX™

### Wrap Details

#### ACCURATE, MEASURED THERAPEUTIC COMPRESSION FOR SELF-DIRECTED TREATMENT AND SELF-MANAGED CARE

The Compreflex™ range of products are low profile compression devices designed to deliver moderate to firm compression to the lower extremities. Each Compreflex™ device is made of a soft conforming inelastic material, offering low to minimal stretch for proper fit and the gradient application of compression. The Compreflex™ garments are designed with key features that support the patient, including our patented Accutab® Devices which ensure the accurate application of the prescribed dosage of therapeutic compression, an interior stretch panel for easy donning, a contoured calf for limb shaping, and opposing straps that are rounded for improved fit and patient comfort.

#### INTELLIGENT DESIGN

Garment features a contoured calf, providing a superior fit

#### PREMIUM MATERIALS

Soft, conforming material provides a comfortable fit

#### EASY APPLICATION

With a front stretch panel to secure the garment in place and straps that roll back, the garment is easy to don and adjust, improving compliance

#### SUPERIOR QUALITY

Each garment is double stitched to ensure durability



20–50 mmHg



## Features & Benefits


- High quality inelastic wrap composed of durable breathable inelastic material with moisture wicking capacities
- Contoured design helps the wrap conform to the natural shape of the leg
- Front stretch panel secures the garment in place, while the Velcro straps assist for easy application
- Ability to customize with different liner and boot options

## Conditions

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post surgical edema

## How it works


### PROVEN ACCURACY FOR RESTING PRESSURE



Accutab® has been rigorously tested according to established guidelines in the International Compression Club Consensus documents. Testing was conducted using PicoPress, which is recognized for its high accuracy and reproducibility.

## Effective working pressure

- Inelastic compression reduces leg edema and improves venous blood flow.



## Safe, precise resting pressure

- Under the supervision of a healthcare provider, inelastic compression can be used safely for patients with mild to moderate peripheral arterial disease.
- Accutab® delivers precise dosages for controlled initial resting pressure when donning a Sigvaris Compreflex™ garment.



## Accutab® Devices

(see page 58 for details)



## Compreboot® Standard

Variable compression  
(see page 57 for details)



## Basic Liner

Non-Compressive  
One pair of Liners  
included with each leg  
(see page 56 for details)

## Standard CALF & FOOT

20–50 mmHg



For those who require a calf and foot garment  
Includes:

- Compreflex™ Wrap
- Non-compressive Basic Liner
- Non attached foot piece

## PRODUCT CODE

Size Regular	Black	Beige	Size Tall	Black	Beige
Small	1401-BKR	1411-BKR	Small	1401-BKT	1411-BKT
Medium	1402-BKR	1412-BKR	Medium	1402-BKT	1412-BKT
Large	1403-BKR	1413-BKR	Large	1403-BKT	1413-BKT
X Large	1404-BKR	1414-BKR	X Large	1404-BKT	1414-BKT
XX Large	1405-BKR	1415-BKR	XX Large	1405-BKT	1415-BKT

**Note:** • 1401-BKR, 1411-BKR, 1401-BKT, and 1411-BKT includes a Regular Length Boot

• 1402-BKR, 1412-BKR, 1402-BKT, 1412-BKT, 1403-BKR, 1413-BKR, 1403-BKT, and 1413-BKT includes a Medium/Large Regular Length Boot


• 1404-BKR, 1414-BKR, 1404-BKT, 1414-BKT, 1405-BKR, 1415-BKR, 1405-BKT, and 1415-BKT includes a X-Large/XX-Large Regular Length Boot

## Sizing Chart

### COMPREFLEX™

### CALF SIZING CHART (A-C)

	Small	Medium	Large	X Large	XX Large
<b>C</b>	29–39 cm	34–44 cm	40–50 cm	46–56 cm	54–64 cm
<b>A</b>	19–26 cm	21–30 cm	26–36 cm	31–41 cm	36–46 cm



### G Length

Regular: 30–33.9 cm  
Tall: 34–42 cm

### Foot Size:

**I:** Measure circumference of ankle bend and heel

**J:** Measure circumference across metatarsal heads

**K:** Measure **medial** length from heel to 1<sup>st</sup> metatarsal head

# Compreflex™ Calf options



Standard  
**CALF**

20–50 mmHg

For those who require an alternative foot garment

Includes:

- Compreflex™ Wrap
- Non-compressive Basic Liner

Can be matched with alternative foot pieces\* (see page 57).

\* When paired with a basic liner, a foot piece is required.



Transition  
**CALF**

20–50 mmHg

For those who require a compressive liner for the foot and ankle

Includes:

- Compreflex™ Wrap
- Transition Liner, which provides 10–15 mmHg in the foot and ankle

Ability to wear regular shoe over liner



Complete  
**CALF**

20–50 mmHg

For those who require a compressive liner with compression in the foot and ankle, and graduated compression in the calf

Includes:

- Compreflex™ Wrap
- Complete Liner, which provides 10–15 mmHg in the foot and ankle, and graduated compression in the calf

Ability to wear regular shoe over liner

PRODUCT CODE		
Size Regular	Black	Beige
Small	1401-NF-BKR	1411-NF-BKR
Medium	1402-NF-BKR	1412-NF-BKR
Large	1403-NF-BKR	1413-NF-BKR
X Large	1404-NF-BKR	1414-NF-BKR
XX Large	1405-NF-BKR	1415-NF-BKR

Size Tall		
Small	1401-NF-BKT	1411-NF-BKT
Medium	1402-NF-BKT	1412-NF-BKT
Large	1403-NF-BKT	1413-NF-BKT
X Large	1404-NF-BKT	1414-NF-BKT
XX Large	1405-NF-BKT	1415-NF-BKT

PRODUCT CODE		
Size Regular	Black	Beige
Small	1401-UC-BKR	1411-UC-BKR
Medium	1402-UC-BKR	1412-UC-BKR
Large	1403-UC-BKR	1413-UC-BKR
X Large	1404-UC-BKR	1414-UC-BKR
XX Large	1405-UC-BKR	1415-UC-BKR

Size Tall		
Small	1401-UC-BKT	1411-UC-BKT
Medium	1402-UC-BKT	1412-UC-BKT
Large	1403-UC-BKT	1413-UC-BKT
X Large	1404-UC-BKT	1414-UC-BKT
XX Large	1405-UC-BKT	1415-UC-BKT

PRODUCT CODE		
Size Regular	Black	Beige
Small	1401-FC-BKR	1411-FC-BKR
Medium	1402-FC-BKR	1412-FC-BKR
Large	1403-FC-BKR	1413-FC-BKR
X Large	1404-FC-BKR	1414-FC-BKR
XX Large	1405-FC-BKR	1415-FC-BKR

Size Tall		
Small	1401-FC-BKT	1411-FC-BKT
Medium	1402-FC-BKT	1412-FC-BKT
Large	1403-FC-BKT	1413-FC-BKT
X Large	1404-FC-BKT	1414-FC-BKT
XX Large	1405-FC-BKT	1415-FC-BKT

See page 56 for Liner Sizing chart



## Liner Options



Basic Liner

No compression

- Non-compressive liner for the protection of the skin
- Easy to apply knee-high liners
- drirelease moisture wicking and odour control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Heathered Silver colour only

Size: One size fits all

Pair	10 Pair Pack
1700 - NK1	1700 - NKC
Included With Compreflex™ Standard Calf & Foot, Compreflex™ Standard Calf, Coolflex™ Standard Calf	



Transition Liner

10–15 mmHg

- Mild compression in foot and ankle only
  - drirelease moisture wicking and odour control
  - Extra thick sole from toe to heel for added comfort
  - Silver ionized treatment for bacteriostatic properties
  - Available in Black and Beige
- \* When ordering a Coolflex standard without a foot piece, please order a transition liner.

PRODUCT CODE (Included with Compreflex™ Transition)		
Size	Pair (Black)	Pair (Beige)
Small	1601-CK	1611-CK
Medium	1602-CK	1612-CK
Large	1603-CK	1613-CK
X Large	1604-CK	1614-CK
XX Large	1605-CK	1615-CK



Complete Liner

10–15 mmHg

- Mild compression in foot and ankle with graduated compression in calf
- Comfortable knee-high liners
- drirelease moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige

PRODUCT CODE (Included with Compreflex™ Complete)		
Size	Pair (Black)	Pair (Beige)
Small	1601-CK-F1	1611-CK-F1
Medium	1602-CK-F1	1612-CK-F1
Large	1603-CK-F1	1613-CK-F1
X Large	1604-CK-F1	1614-CK-F1
XX Large	1605-CK-F1	1615-CK-F1

## LINERS

### TRANSITION AND COMPLETE LINERS SIZING CHART

Size	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe
Small	19–26 cm	29–39 cm	<7½	<6½
Medium	21–30 cm	34–44 cm	7.5–9.5	6.5–8.5
Large	26–36 cm	40–50 cm	10–12	9–11
X Large	31–41 cm	46–56 cm	12.5–15	11.5–14
XX Large	36–46 cm	54–64 cm	12.5–15	11.5–14

# Compreboot® Options



Compreboot®  
Lite **FOOT**

Variable  
compression

- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design may be trimmed to fit
- Provides mild, therapeutic compression



Compreboot®  
Standard **FOOT**

Variable  
compression

- Four-way stretch for firm compression around ankle
- Ideal for mild to moderate chronic edema
- Low stretch material conforms to foot shape



Compreboot®  
Plus **FOOT**

Variable  
compression

- Chip foam lining creates high / low pressure ridges
- Removable gripper patches provide traction
- Provides moderate, therapeutic compression

## PRODUCT CODE (Available in Black Only)


Size	Black
<b>Small / Medium / Large</b>	
Regular	1803-BLR
Long	1803-BLL
<b>X Large / XX Large</b>	
Regular	1805-BLR
Long	1805-BLL

## PRODUCT CODE

Size	Black	Beige
<b>Small</b>		
Regular	1801-BTR	1811-BTR
Long	1801-BTL	1811-BTL
<b>Medium / Large</b>		
Regular	1803-BTR	1813-BTR
Long	1803-BTL	1813-BTL
<b>X Large / XX Large</b>		
Regular	1805-BTR	1815-BTR
Long	1805-BTL	1815-BTL

## PRODUCT CODE (Available in Black Only)

Size	Black
<b>Small</b>	
Regular	1801-BPR
Long	1801-BPL
<b>Medium / Large</b>	
Regular	1803-BPR
Long	1803-BPL
<b>X Large / XX Large</b>	
Regular	1805-BPR
Long	1805-BPL



## Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure **medial** length from heel to 1<sup>st</sup> metatarsal head

## COMPREFLEX® LITE

### FOOT SIZING CHART (I, J, K)

	<b>Small / Medium / Large</b>		<b>X Large / XX Large</b>	
	Regular	Long	Regular	Long
<b>I</b>	28–36 cm	28–36 cm	≤44 cm	≤44 cm
<b>J</b>	22–26 cm	22–26 cm	≤33 cm	≤33 cm
<b>K</b>	14–18 cm	19–23 cm	20–22 cm	25–27 cm

## COMPREFLEX® STANDARD AND PLUS

### FOOT SIZING CHART (I, J, K)

	<b>Small</b>		<b>Medium / Large</b>		<b>X Large / XX Large</b>	
	Regular	Long	Regular	Long	Regular	Long
<b>I</b>	28–36 cm	28–36 cm	<39 cm	<39 cm	<44 cm	<44 cm
<b>J</b>	22–26 cm	22–26 cm	<30 cm	<30 cm	<33 cm	<33 cm
<b>K</b>	14–18 cm	19–23 cm	18–20 cm	23–25 cm	20–22 cm	25–27 cm



## Compreflex™ Strap Extender

- Provides up to 10 cm of additional circumferential size range to Compreflex™ calf wraps
- Detaches easily
- Potential solution for patients outside of standard size ranges of Compreflex™ calf wraps or needing custom size

	1 Unit	10 Unit Pack
Black	1401-SE1	1401-SEB
Beige	1411-SE1	1411-SEB



## Accutab® Devices

The Accutab® accurately delivers the level of compression prescribed by the health care provider. An easy to use device, the Accutab® ensures the safe delivery of compression in doses of 20–30 mmHg, 30–40 mmHg, or 40–50 mmHg.





## Standard Thigh

20–50 mmHg



### Features

- The Compreflex™ Standard Thigh is an inelastic compression device for the thigh. The product includes a Belt and Hip Attachment to support the product throughout daily use.
- Pair the Compreflex™ Standard Thigh with the Compreflex™ Standard Knee (sold separately) to create a complete compression system for the upper leg.


### Ideal for

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post surgical edema

### Includes

- Compreflex™ Standard Thigh
- 1 Pair of non-compressive Cotton Liners per leg
- Belt and Hip Attachment
- Accutab® Devices

### Colours



## COMPREFLEX™ STANDARD THIGH SIZE CHART (D-F) & PRODUCT CODE

	Small	Medium	Large	X Large
<b>F</b>	48–58 cm	56–66 cm	64–74 cm	74–84 cm
<b>E</b>	43–53 cm	51–61 cm	58–68 cm	68–78 cm
<b>D</b>	38–48 cm	46–56 cm	53–63 cm	63–73 cm

### BLACK

Left	Regular	1401-TCR-L	1402-TCR-L	1403-TCR-L	1404-TCR-L
	Tall	1401-TCT-L	1402-TCT-L	1403-TCT-L	1404-TCT-L
Right	Regular	1401-TCR-R	1402-TCR-R	1403-TCR-R	1404-TCR-R
	Tall	1401-TCT-R	1402-TCT-R	1403-TCT-R	1404-TCT-R

### BEIGE

Left	Regular	1411-TCR-L	1412-TCR-L	1413-TCR-L	1414-TCR-L
	Tall	1411-TCT-L	1412-TCT-L	1413-TCT-L	1414-TCT-L
Right	Regular	1411-TCR-R	1412-TCR-R	1413-TCR-R	1414-TCR-R
	Tall	1411-TCT-R	1412-TCT-R	1413-TCT-R	1414-TCT-R






## Standard Knee

20–50 mmHg



### Features

- Provides accurate, inelastic compression levels of 20–50 mmHg
- Includes inner donning sleeve for easy application
- Product includes a Butterfly Strap for hard to compress areas



### COMPREFLEX™ STANDARD KNEE SIZE CHART (C1-D)

	Small	Medium	Large	X Large
D	38–48 cm	46–56 cm	53–63 cm	58–73 cm
C1	29–39 cm	34–44 cm	39–49 cm	44–65 cm
<b>Black</b>	1401-KP	1402-KP	1403-KP	1404-KP
<b>Beige</b>	1411-KP	1412-KP	1413-KP	1414-KP

\*Arrows indicate reference points of D and C1. Length measurement is not required.

# COMPREKNEE™

Variable compression




## Standard Knee CLOSED PATELLA



### Features

- Garment covers the knee, providing full coverage
- Straps provide a full-range of adjustment to create a comfortable fit
- Conforming material for a firm secure fit

### Colours



### COMPREKNEE™ SIZE CHART (C1-D)

	Small	Medium	Large / X Large
D	38–48 cm	46–56 cm	53–63 cm
C1	29–39 cm	34–44 cm	39–49 cm
Black	1101-KP	1102-KP	1103-KP
Beige	1111-KP	1112-KP	1113-KP

\*Arrows indicate reference points of D and C1. Length measurement is not required.

## COMPRESHORT®

10-15 mmHg



IDEAL FOR PATIENTS WITH ABDOMINAL OR TRUNCAL EDEMA WHO NEED MILD, SUPPORTIVE COMPRESSION

Product available in short and capri.



### Features

The Compreshort® offers mild, supportive compression for the treatment of abdominal or truncal and upper thigh edema. Available in short and capri lengths, patients with lymphedema and lipedema benefit from the therapeutic compression and support the garment provides.

### Ideal for

- Lymphedema
- Lipedema


### Colour



Black

COMPRESHORT® SIZING CHART (Waist, Hip, Midthigh)				PRODUCT CODE	
	Waist	Hip	Midthigh	Product SKU (Short)	Product SKU (Capri)
Small	69–76 cm	84–91 cm	41–46 cm	1301-CS	1301-CP
Small Max	76–86 cm	91–102 cm	41–46 cm	1301-CM	1301-CPM
Medium	76–86 cm	91–102 cm	46–53 cm	1302-CS	1302-CP
Medium Max	86–97 cm	102–112 cm	46–53 cm	1302-CM	1302-CPM
Large	86–97 cm	102–112 cm	53–61 cm	1303-CS	1303-CP
Large Max	97–112 cm	112–127 cm	53–61 cm	1303-CM	1303-CPM
X Large	97–112 cm	112–127 cm	61–71 cm	1304-CS	1304-CP
X Large Max	112–132 cm	127–147 cm	61–71 cm	1304-CM	1304-CPM
2X Large	112–132 cm	127–147 cm	71–81 cm	1305-CS	1305-CP
2X Large Max	132–152 cm	147–168 cm	71–81 cm	1305-CM	1305-CPM
3X Large	132–152 cm	147–168 cm	81–91 cm	1306-CS	1306-CP
3X Large Max	152–172 cm	168–188 cm	81–91 cm	1306-CM	1306-CPM
4X Large	152–172 cm	168–188 cm	91–101 cm	1307-CS	1307-CP
4X Large Max	172–193 cm	188–208 cm	91–101 cm	1307-CM	1307-CPM
5X Large	172–193 cm	188–208 cm	101–112 cm	1308-CS	1308-CP

**Note:** Available in Black Only





**Compression socks and hosiery for everyday wellbeing**



## Care for Yourself

### Why choose **Consumer Health**?

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease, prevention of swelling and varicose veins
- Simplified sizing

#### CONSUMER HEALTH PRODUCTS

Series	Product	Page
<b>146 / 186</b>	Casual Cotton	66
<b>146/T / 186/T</b>	Traveno Travel Socks	67
<b>151 / 191</b>	Sea Island Cotton	68
<b>152 / 192</b>	All-Season Merino Wool	70
<b>153 / 193</b>	Linen	71
<b>143 / 183</b>	Microfiber Shades	72
<b>120</b>	Sheer Fashion	74
<b>170</b>	Soft Silhouette Leggings	76
<b>401</b>	Athletic Recovery Socks	78
<b>412V</b>	Performance Sleeves	79
<b>421</b>	Merino Outdoor Socks	80
<b>160</b>	Eversoft Diabetic Socks	81



Featured  
colour is  
Black

# CASUAL COTTON 146/186

for women & men

## Compression Level

15–20 mmHg

## Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation

## Fiber Content


- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free

## SIZING CHART

146	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
186	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

## STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●	●	●	●
----------	---	---	---	---	---

Product Code Example:  
146CB99

Series	Style	Size	Colour
146	C	B	99

Featured  
colour is  
Black

# TRAVENO TRAVEL SOCKS 146/T - 186/T

for women & men

## Compression Level

15–20 mmHg

## Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation
- Box cover promotes Traveno as a Travel friendly sock

## Fiber Content


- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free

## SIZING CHART

146/T	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
186/T	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

## STYLES AND COLOURS

Closed Toe = ●



Calf (C)



Product Code Example: 146CB99/T

Series	Style	Size	Colour	Travel
146	C	B	99	/T



# SEA ISLAND COTTON 151/191

for women & men

## Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 22

## Features and Benefits

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

## Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free


Featured colour  
is Black



**151**  
Women's  
Calf Style

**15–20 mmHg**  
Closed Toe Calf Size A in Navy

**151CA10**



**191**  
Men's  
Calf Style

**15–20 mmHg**  
Closed Toe Calf Size B in Brown


**191CB11**

## SIZING CHART

151	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23 cm (7"-9")	19-25 cm (7.5"-10")	20-27 cm (8"-10.5")
191	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25 cm (7.5"-10")	20-27 cm (8"-10.5")	21-29 cm (8.5"-11.5")
	Size	A	B	C

## STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---

## Product Code Example: 151CA99

Series	Style	Size	Colour
151	C	A	99



Featured colour  
is Black

# ALL-SEASON MERINO WOOL 152/192

for women & men

## Compression Level

15–20 mmHg

Also available in 20–30 mmHg,  
see page 24

## Features and Benefits

- High quality extra-fine Merino Wool makes this sock soft and breathable
- Moisture wicking capabilities and naturally odour resistant
- Thermoregulating for year-round wear and keeps feet dry
- Smooth, flat toe seam

## Fiber Content


- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free

## SIZING CHART

152	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
192	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

## STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---

Product Code Example: 152CB99

Series	Style	Size	Colour
152	C	B	99



Featured  
colour is  
Lavender

# LINEN 153/193

for women & men

## Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 26

## Features and Benefits

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

## Fiber Content


- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free

## SIZING CHART

153	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
193	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

## STYLES AND COLOURS

Closed Toe = ●



Product Code Example: 153CB207

Calf (C)	●	●	●	●
----------	---	---	---	---

Series	Style	Size	Colour
153	C	B	207



# MICROFIBER SHADES 143/183

for women & men

## Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 30

## Features and Benefits

- 100% of microfiber sits next to skin for a soft feel
- Fashionable, durable socks for everyday wear
- Excellent breathability and moisture management
- Glides on and off with ease
- Smooth, flat toe seam

## Fiber Content

- 87% Nylon, 13% Spandex
- Latex free

Featured colour is  
Royal Blue Argyle

**143**  
Women's  
Calf Style



**15–20 mmHg**  
Closed Toe Calf Size B  
in Purple Argyle

**143CB46**

**183**  
Men's  
Calf Style



**15–20 mmHg**  
Closed Toe Calf Size C  
in Graphite Stripe

**183CC17**

SIZING CHART

143	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23 cm (7"-9")	19-25 cm (7.5"-10")	20-27 cm (8"-10.5")
183	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25 cm (7.5"-10")	20-27 cm (8"-10.5")	21-29 cm (8.5"-11.5")
	Size	A	B	C

STYLES AND COLOURS

Women Closed Toe = ●							
	Calf (C)	●	●	●	●	●	●

Men Closed Toe = ●					
	Calf (C)	●	●	●	●

Product Code Example: 183CB48

Series	Style	Size	Colour
183	C	B	48



Featured  
colour is  
Honey

# SHEER FASHION 120

for women

## Compression Level

15–20 mmHg

## Features and Benefits

- Incredible sheer design in a variety of fashionable colours
- Sheer toe and heel area
- Comfortable stretch and soft fabric, excellent for everyday wear

## Fiber Content

- 79% Nylon, 21% Spandex  
(Content can vary slightly by style)
- Latex free



**15–20 mmHg**  
Closed Toe Calf Size C  
in Charcoal

**120CC12**



**15–20 mmHg**  
Closed Toe Pantyhose Size A  
in Black

**120PA99**

## PANTYHOSE SIZING CHART

cm/ft	Kg/Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230
145 cm / 4'9"																												
147 cm / 4'10"																												
150 cm / 4'11"																												
152 cm / 5'0"																												
155 cm / 5'1"																												
157 cm / 5'2"																												
160 cm / 5'3"																												
162 cm / 5'4"																												
165 cm / 5'5"																												
167 cm / 5'6"																												
170 cm / 5'7"																												
173 cm / 5'8"																												
175 cm / 5'9"																												
178 cm / 5'10"																												
180 cm / 5'11"																												
183 cm / 6'0"																												

## CALF / THIGH SIZING CHART

120	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23 cm (7"-9")	19-25 cm (7.5"-10")	20-27 cm (8"-10.5")
	Size	A	B	C

Weight Guidelines For Thigh Style: Size A: <135Lbs ; Size B: <180Lbs ; Size C: <230Lbs.

## STYLES AND COLOURS

Closed Toe = ●  
Open Toe = ○

	CHARCOAL 12	TAUPE 29	HONEY 33	GOLDEN 36	CAFÉ 73	BLACK 99
--	----------------	-------------	-------------	--------------	------------	-------------

Calf (C)	●	●○	●○	●○	●	●
Thigh (N)	●	●	●	●	●	●
Pantyhose (P)	●	●	●	●	●	●

## Product Code Example: 120NC29

Series	Style	Size	Colour
120	N	C	29

## Product Code Example: 120PC99

Series	Style	Size	Colour
120	P	C	99





# SOFT SILHOUETTE LEGGINGS 170

for women

## Compression Level

15-20 mmHg

## Features and Benefits

- Advanced soft stretch microfiber fabric shapes legs and smooths skin
- Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort
- Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts

## Fiber Content

- 83% Nylon, 17% Spandex
- Latex free



**15-20 mmHg**  
Legging Size B  
in Black

**170LB99**



**Non-binding ankle cuff**  
stays comfortably in place  
for all day wear

## SIZING CHART

cm / ft	Kg / Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230
147 cm / 4'10"																												
150 cm / 4'11"																												
152 cm / 5'0"																												
155 cm / 5'1"																												
157 cm / 5'2"																												
160 cm / 5'3"																												
162 cm / 5'4"																												
165 cm / 5'5"																												
167 cm / 5'6"																												
170 cm / 5'7"																												
173 cm / 5'8"																												
175 cm / 5'9"																												
178 cm / 5'10"																												
180 cm / 5'11"																												
183 cm / 6'0"																												

Guidelines for Ankle Circumference

**Size A:** 18–23 cm (7"–9"), **Size B & C:** 20–26 cm (8"–10"), **Size D:** 22–29 cm (8.5"–11.5")

## STYLES AND COLOURS

Open Foot = ○

BLACK  
99

Leggings (L)

○

## Product Code Example: 170LB99

Series	Style	Size	Colour
170	L	B	99



Featured  
colour is  
White



15–20 mmHg  
Closed Toe Calf Size M  
in Black

401CM99

# ATHLETIC RECOVERY SOCKS 401

for women & men

## Compression Level

15–20 mmHg

## Features and Benefits

- drirelease® yarns provide moisture wicking and odour control
- Comfortable and easy to put on and take off
- Helps reduce lactic acid buildup
- Mild cushioning in sole, heel, and toe area for added comfort
- Smooth, flat toe seam

## Fiber Content

- 67% drirelease® (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free

## SIZING CHART

401	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15	≥15.5
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14	14.5–17
	Ankle Circumference	18–24 cm (7"–9.5")	19–26 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")	21–29 cm (8.5"–11.5")
	Size	S	M	L	X	2

## STYLES AND COLOURS

Closed Toe = ●

WHITE  
00

BLACK  
99

Calf (C)	●	●
----------	---	---

Product Code Example: 401CM00

Series	Style	Size	Colour
401	C	M	00

# PERFORMANCE SLEEVES

## 412V

for women & men

### Compression Level

20–30 mmHg

### Features and Benefits

- Stabilizes and supports active muscles
- Ideal for running, triathlons, endurance sports, and gym workouts
- Fabric resists deterioration in salt, chlorine, or fresh water
- Quick drying
- High-visibility reflective transfer

### Fiber Content

- 75% Nylon, 25% Spandex
- Latex free

Featured colour is Purple

### SIZING CHART

412V	Ankle Circumference	18–21.5 cm (7"–8.5")	21.5–25.5 cm (8.5"–10")	25.5–29 cm (10"–11.5")	29–33 cm (11.5"–13")
	Calf Circumference	28–39 cm (11"–15.5")	31.5–45 cm (12.5"–17.5")	35.5–51 cm (14"–20")	40.5–58 cm (16"–23")
	Size	S	M	L	X

### STYLES AND COLOURS

Calf Sleeve = ○	WHITE 00	BLUE 50	PURPLE 58*	BLACK 99
Calf Sleeve (V)	○	○	○	○

\* Purple is not available in X-Large

Product Code Example: 412VL99

Series	Style	Size	Colour
412	V	L	99



Featured colour  
is Charcoal

# MERINO OUTDOOR SOCKS 421

for women & men

## Compression Level

15–20 mmHg

Also available in 20–30 mmHg,  
see page 38

## Features and Benefits

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned, padded soles help prevent blisters
- Mesh flex zones for comfort and ventilation
- Smooth, flat toe seam
- Soft and comfortable
- Perfect for outdoor work environments and activities

## Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free

## SIZING CHART

421	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14
	Ankle Circumference	18–24 cm (7"–9.5")	19–26 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	S	M	L	X

## STYLES AND COLOURS

Closed Toe = ●

CHARCOAL  
12

OLIVE  
31

Calf (C)	●	●
----------	---	---

## Product Code Example: 421CS12

Series	Style	Size	Colour
421	C	S	12



# EVERSOFT DIABETIC SOCKS 160

for women & men

## Compression Level

8–15 mmHg

## Features and Benefits

- Padded sole, heel, and toe for comfort and protection
- Smooth, flat toe seam
- Compression-free foot
- drirelease yarns provide moisture wicking and odour control
- Suitable for individuals with limited dexterity or lower tolerance for compression
- Ideal for individuals with sensitive skin

## Fiber Content

- 67% drirelease (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free

## SIZING CHART

160	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14
	Ankle Circumference	18–24 cm (7"–9.5")	19–26 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	S	M	L	X

## STYLES AND COLOURS

Closed Toe = ●

WHITE  
00

Calf (C)



## Product Code Example: 160CS00

Series	Style	Size	Colour
160	C	S	00

**Donning Aids and Product Care Accessories**



## Make Life Easier

### Why choose **Accessories**?

- Aid in putting on and removing stockings
- Stocking cleaning to maximize usage
- Extra liners for added comfort

#### ACCESSORIES

Product	Page
Gloves	84
Cotton Underliner	85
Washing Solution	85
It Stays Roll-On Adhesive	85
Sim-Slide	86
Magnide On/Off	86
Sigvaris Doff N' Donner	87
Sigvaris Doff N' Donner Cone	87
Melany	87

### GLOVES



#### Rubber Gloves\*

Ridged - The practical and protective way to apply and remove your compression stockings.



#### Latex-Free Gloves

Gripping donning gloves recommended for use in a latex-free environment.



#### Textile Gloves\*

These high-quality gloves have a latex surface to provide a good grip when putting on and taking off compression stockings.

#### PACKAGE: 12 PER PACK

SIZE	PRODUCT CODE
Small	591R400S
Medium	591R400M
Large	591R400L
X-Large	591R400X
ASST'D 4X S/M/L	591R400A

#### PACKAGE: 12 PER PACK

SIZE	PRODUCT CODE
Small	592R400S
Medium	592R400M
Large	592R400L
X-Large	592R400X

#### PACKAGE: 12 PER PACK

SIZE	PRODUCT CODE
X-Small	593R400XS
Small	593R400S
Medium	593R400M
Large	593R400L
X-Large	593R400X
ASST'D 4x XS/S/M	593R400A

#### RUBBER GLOVE AND LATEX-FREE GLOVE SIZING CHART

Measure Circumference of the Hand at the Palm

Hand Circumference	16.5–18 cm (6.5"–7")	19–20.5 cm (7.5"–8")	21.5–23 cm (8.5"–9")	24–25.5 cm (9.5"–10")
	Small	Medium	Large	X-Large

#### TEXTILE GLOVE SIZING CHART

Measure Length and Width of the Hand

Hand Length**	20.5–21 cm	21.5–22 cm	22.5–23 cm	23.5–24 cm	24.5–25 cm
Hand Width	7–7.5 cm	7.5–8 cm	8–8.5 cm	9–9.5 cm	10–10.5 cm
	X-Small	Small	Medium	Large	X-Large

\* **Caution:** Rubber Gloves and Textile Gloves contain natural rubber latex that may cause allergic reaction.

\*\* **Hand Length:** measured from the tip of the longest finger to the crease under the palm.

## LINERS



### Cotton Underliner

10–15 mmHg

Can be used as an underlayer to protect fragile skin, secure wound dressings, add 10–15 mmHg of compression, or create an easy glide surface to aid in putting on and taking off compression stockings.

*For sizing, please refer to the Natural Rubber Calf sizing on page 43.*

#### PACKAGE: 2 PER PACK

SIZE	PRODUCT CODE	NATURAL RUBBER SIZING
Small / Short	588USS	S1/S3
Small / Long	588USL	S2/S4
Medium / Short	588UMS	M1/M3
Medium / Long	588UML	M2/M4
Large-XL / Short	588ULS	L1/L3 - X1/X3
Large-XL / Long	588ULL	L2/L4 - X2/X4

## WASHING SOLUTION

### Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for washing delicate garments!

#### PACKAGE OPTIONS

CONTAINS	PRODUCT CODE
Six 16 oz. Bottles (473ml)	586W900
One 16 oz. Bottle (473ml)	586W950
Twelve 2 oz. Bottles (60ml)	586W975



## IT STAYS!

### Roll-on body adhesive

A gentle roll-on body adhesive that washes off with water and may be used to hold compression garments in place against the skin.

IT STAYS! is pliable and moves with the body.



SIZE	PRODUCT CODE
One 2 oz Bottle (60ml)	581I350
Twelve 2 oz. Bottles (60ml)	581I300



## DONNING AIDS



### Sim-Slide

(For open toe Styles)

SIZE	PRODUCT CODE
Small	587SLIDESM
Medium	587SLIDEMD
Large	587SLIDELG
X-Large	587SLIDEXL

#### SIZING CHART For Sim-Slide

Women's Shoe Size	≤4-6	6.5-8	8.5-9.5	≥10
Men's Shoe Size	≤3-4.5	5-6.5	7-8	≥8.5
Size	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>



*Click here to see how to use  
the Sim-Slide®*  
[bit.ly/SIGssyt](https://bit.ly/SIGssyt)



### Magnide On/Off

(For closed toe Styles)

SIZE	PRODUCT CODE
Medium	587MAGMD
Large	587MAGLG
X-Large	587MAGXL

#### SIZING CHART For Magnide On/Off

Women's Shoe Size	≤8.5	9-12	≥ 12.5
Men's Shoe Size	≤7.5	8-11	≥ 11.5
Size	<b>M</b>	<b>L</b>	<b>XL</b>

\* XL suitable for extra-wide feet



*Click here to see how to use  
the Magnide On/Off®*  
[bit.ly/49sxcYX](https://bit.ly/49sxcYX)

## DONNING AIDS



### SIGVARIS Doff N' Donner

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use. Latex free.

SIZE	PRODUCT CODE
One Size	583D100



*Click here to see how to use  
Sigvaris DOFF N' DONNER*  
[bit.ly/3QSBcKU](https://bit.ly/3QSBcKU)



### SIGVARIS Doff N' Donner Cone

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner. The bottom of the base firmly attaches to most clean, flat surfaces via suction.

SIZE	PRODUCT CODE
One Size	583D200



### MELANY

MELANY from Sigvaris makes both donning and doffing compression garments simple. Its innovative one-piece design with adjustable handles was made especially for users with limited range of motion and finger dexterity.\*

SIZE	CALF CIRCUMFERENCE	PRODUCT CODE
X-Small	25–30 cm	579SXS
Small	30–38 cm	579SSM
Medium	38–46 cm	579SMD
Large	46–56 cm	579SLG

\* Individual components not sold separately



*Click here to see how to use  
the MELANY*  
[bit.ly/3SCQYLu](https://bit.ly/3SCQYLu)

**Learn more**



## Help yourself

How can **resources** help?

- Learn about the levels of compression for various indications to better help your clients
- Understand the terms and icons that we use
- Follow specific instructions on how to measure for the perfect fit

### RESOURCES

Section	Page
Symptoms & Indications	90
Glossary of Icons	91
Glossary of Terms	92
How to Measure	93
How to Order	94

SYMPTOMS / INDICATIONS	COMPRESSION LEVEL AT THE ANKLE
<ul style="list-style-type: none"> <li>□ Tired, achy legs</li> <li>□ Occupational or evening edema</li> <li>□ Leg discomfort from long hours of standing or sitting</li> <li>□ Prevent varicosities &amp; swelling during pregnancy</li> <li>□ Preventative for legs predisposed to risk</li> <li>□ Reduce swelling during travel</li> <li>□ Prevent DVT during travel</li> </ul>	<b>15–20 mmHg</b>
<ul style="list-style-type: none"> <li>□ Heavy, fatigued, aching legs</li> <li>□ Mild edema in lower extremities</li> <li>□ Mild varicosities with minimal edema</li> <li>□ Mild varicosities during pregnancy</li> <li>□ Lymphedema upper extremities</li> <li>□ Post-procedure of small veins</li> <li>□ Prevention or management of DVT</li> </ul> <p><i>Note: A prescription is recommended for stockings with the indications of 20 mmHg or higher</i></p>	<b>20–30 mmHg</b>
<ul style="list-style-type: none"> <li>□ Moderate varicosities with mild to moderate edema (also during pregnancy)</li> <li>□ Post-procedure of larger veins to maintain treatment success</li> <li>□ CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations)</li> <li>□ Management and treatment of DVT or Post-Thrombotic Syndrome</li> <li>□ Superficial phlebitis</li> <li>□ Lymphedema after decongestant therapy to maintain reduction</li> <li>□ Orthostatic / postural hypotension</li> <li>□ Post-fracture, post-traumatic edema</li> </ul>	<b>30–40 mmHg</b>
<ul style="list-style-type: none"> <li>□ Lymphedema after decongestant therapy to maintain reduction</li> <li>□ Severe Post-Thrombotic Syndrome</li> <li>□ Severe varicosities and/or edema</li> </ul>	<b>40 mmHg or higher</b>
<ul style="list-style-type: none"> <li>□ Severe arterial insufficiency, intermittent claudication, ischemia</li> <li>□ Uncontrolled congestive heart failure</li> </ul>	<b>Contraindications</b>
<ul style="list-style-type: none"> <li>□ Skin sensitivities or allergies</li> <li>□ Advanced neuropathy with sensory loss in leg or foot</li> <li>□ Diabetes with arterial compromise</li> <li>□ Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician</li> <li>□ Eczema, hypodermatitis, dermatosclerosis, dermatitis</li> </ul>	<b>Cautions</b>
<p>The therapeutic effect of graduated compression stockings is based on the compression level. With the proper CEAP classification, the correct compression level at the ankle can be chosen. For more severe symptoms, consult with a vascular specialist.</p>	
<p><i>No liability accepted for non-observance of contraindications and cautions.</i>            Fronek HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.</p>	





Sensitive Skin



Thermoregulating



Breathable



Moisture Wicking



Soft




Durable



Cotton




Wool



Microfiber



Natural Rubber



Largest Size Range



Sheer



Medium Sheer



Opaque



Stylish



Stylish



Elegant



Colours



Patterns



Sensinno® Grip-top  
Band



Anatomically  
Contoured



Cushioned Sole



Anti-bacterial



Hypoallergenic  
Natural fibers



Moisture  
Management

TERM	DEFINITION	SIGVARIS QUALITY STANDARD
CHRONIC VENOUS DISEASE (CVD)	A disease state comprised of a wide range of disorders including tired, achy legs, varicose veins, Chronic Venous Disease (CVD), skin issues, and venous leg ulcers (open and healed).	SIGVARIS GROUP manufactures MEDICAL and CONSUMER HEALTH CARE products to help reduce the risk of, improve symptoms of, and provide treatment for CVD.
DOUBLE-COVERED YARN	Spandex fibers are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort.	Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort.
DEEP VEIN THROMBOSIS (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.	SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening.
GRADUATED COMPRESSION VS. ANTI-EMBOLISM	Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for non-ambulatory patients.	SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation.
HYPOALLERGENIC	Designed to minimize the chance of an allergic reaction.	SIGVARIS GROUP offers Latex free products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin.
ISO CERTIFIED	SIGVARIS GROUP is ISO (International Organization for Standardization) certified, indicating we have the required quality systems in place.	SIGVARIS GROUP follows ISO Standard 13485:2016 for our products manufactured in North and South America.
MICROFIBER	An extremely fine synthetic fiber knitted into textiles with the texture and drape of natural fiber fabric. This fabric has excellent breathability.	Used to create products with a soft and comfortable feel. Find this fiber in Soft Opaque 840, Microfiber Patterns 830, Microfiber Shades 143/183, and Microfiber 820.
mmHg (MILLIMETERS OF MERCURY)	Defined as the extra pressure generated by a column of mercury one millimeter high; the same unit of measurement used for blood pressure.	SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement.
MERINO WOOL	Merino is the softest and finest wool available in the world and is inherently itch-free.	SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in Thermoregulating Wool 422, Merino Outdoor Socks 421, Merino Wool 240, and All-Season Merino Wool 152/192.
SEA ISLAND COTTON	Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness.	SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our Sea Island 220 and Sea Island Cotton 151/191.
SENSINNOV®	A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day.	This top band is featured in Sheer 780, Soft Opaque 840, and Microfiber 820 to help keep thigh-high stockings in place all day long.
SUPIMA COTTON	Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fiber.	SIGVARIS GROUP uses this premium cotton yarn in the Cotton 230, Sea Island 220, Sea Island Cotton 151/191, Casual Cotton 146/186, and Cushioned Cotton 360, making them comfortable and ideal for those with sensitive skin.

**Please follow these instructions carefully to determine the correct size for different product styles.**

Measurements should always be taken in the following order, and as early in the morning as possible when swelling is minimal.

## CALF STYLE

Open/Closed Toe

<b>1</b>	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. <b>Note: This should be your first measurement.</b>
<b>2</b>	Measure the calf circumference at the fullest part of the calf.
<b>3</b>	Determine <b>the calf length</b> by measuring from the floor up to the bend of the knee. This measurement is taken on the outside of the leg.



**1** Ankle circumference directly above ankle bone



**2** Calf circumference at fullest part of the calf


## PANTYHOSE AND THIGH-HIGH STYLE

Open/Closed Toe

<b>1</b>	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM or LARGE. <b>Note: This should be your first measurement.</b>
<b>2</b>	Measure the calf circumference at the fullest part of the calf.
<b>3</b>	Skip Step 3 (but refer to General Guidelines #6 below).
<b>4</b>	Measure the thigh at its widest circumference.
<b>5</b>	Determine <b>the leg length</b> by measuring from the floor up to the gluteal fold.
<b>6</b>	Determine <b>the hip circumference</b> by measuring the widest part of the hips.



**3** Calf length from the floor to the right-angle bend of the knee



**4** Widest circumference of the thigh

## GENERAL GUIDELINES

<b>1</b>	Measurements should always be taken as early in the morning as possible, when swelling is minimal.
<b>2</b>	Measure against bare skin, not over clothing (Steps 1–5).
<b>3</b>	The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to right for examples).
<b>4</b>	Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
<b>5</b>	Length measurements should be taken from the ground up.
<b>6</b>	Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles.
<b>7</b>	Record shoe size; it may be helpful in determining the best possible fit.



**5** Leg length from the floor up to the gluteal fold



**6** Circumference at the hip

## HOW TO ORDER MEDICAL AND SPECIALTY PRODUCTS



### Product Code Example: 232CMLM99 (page 8)

Series Number	Compression Level	Product Style	Size	Gender	Colour
23 = Cotton <i>(See individual pages, use first 2 digits of series number)</i>	1 = 15–20 mmHg 2 = 20–30 mmHg 3 = 30–40 mmHg 4 = 40–50 mmHg* 5 = 50–60 mmHg*	C = Knee-high or Calf <i>(See individual pages for available styles)</i>	ML = Medium Long <i>(See individual pages for sizing)</i>	M = Men W = Women O = Open toe	99 = Black <i>(See individual pages for available colours)</i>
<b>23</b>	<b>2</b>	<b>C</b>	<b>ML</b>	<b>M</b>	<b>99</b>
* 40–50 mmHg and 50–60 mmHg available in Natural Rubber only					
<b>To order Plus Calf with Grip-top or Plus Size Pantyhose (860P for women only):</b> add “/P” after the colour in the ordering code					
<b>To order Calf Style with Grip-top:</b> add “/S” after the colour in the ordering code <i>(see product pages for availability)</i>					
<b>To order Thigh with Waist Attachment:</b> add “/L” to indicate the Left leg or “/R” to indicate the Right leg after the colour in the ordering code <i>(860W and 500W only)</i>					
<b>To order Motion High Tech and Motion Thermoregulating Wool:</b> The letter for Gender is not used in the ordering code. <b>Example:</b> 412CSM59 = Motion High Tech, 20–30 mmHg, Closed toe Calf, Size Small Medium, Colour Red					

### HOW TO ORDER CONSUMER HEALTH PRODUCTS



#### Product Code Examples:

**151CA99** (page 69)

**191CB99** (page 69)

Series Number	Product Style	Size	Colour
151/191 = SEA ISLAND COTTON  First 3 digits indicate gender 151 = woman 191 = man  <i>(See individual pages for series number)</i>	C = Calf   <i>(See individual pages for available styles)</i>	According to shoe size  A= women's shoe sizes 5-7 B = men's shoe sizes 9-11  <i>(See individual pages for sizing)</i>	99 = Black   <i>(See individual pages for available colours)</i>
<b>151</b>	<b>C</b>	<b>A</b>	<b>99</b>
<b>191</b>	<b>C</b>	<b>B</b>	<b>99</b>
<b>ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, MERINO OUTDOOR SOCKS, and EVERSOFIT DIABETIC SOCKS are UNISEX: the first 3 digits are the same for both genders.</b> <b>Example:</b> 421CL31 = MERINO OUTDOOR SOCK, 15–20 mmHg, Closed toe Calf, Size Large, Colour Olive			



# **mobilis**

by SIGVARIS GROUP

## **SOFT ORTHOPEDIC SUPPORTS AND BRACES**



# INDEX

---

## Intro

98

Mode of action 99  
How to measure 99

## Elbow

100

<b>EPI</b> CARE - Elbow Support	101
<b>EPI</b> BAND - Elbow Strap	101
<b>EPI</b> ACTIVE - Elbow Brace	102

## Wrist and thumb

103

<b>MANU</b> CARE - Wrist Support	104
<b>MANU</b> WRAP - Wrist and Thumb Brace	104
<b>MANU</b> ACTIVE - Wrist Brace	105
<b>MANU</b> SUPPORT - Wrist Brace	105
<b>MANU</b> PLUS SUPPORT - Wrist and Thumb Stabilizer	106
<b>RHIZO</b> SUPPORT - Thumb Brace	106

## Back

107

<b>LUMBO</b> CARE - Back Support	108
<b>LUMBO</b> SUPPORT - Lumbar Brace	108
<b>LUMBO</b> SUPPORT EVO - Lumbar Brace	109
<b>LUMBO</b> INTENSE - Lumbar Brace	109

## Knee

110

<b>GENU</b> CARE - Knee Support	111
<b>GENU</b> WRAP - Knee Wrap	111
<b>GENU</b> BAND - Knee Strap	112
<b>GENU</b> ACTIVE STAB - Stabilizing Knee Brace	112
<b>GENU</b> ACTIVE PAD - Patellar Knee Brace	113

## Ankle

114

<b>MALLEO</b> CARE - Ankle Support	115
<b>MALLEO</b> ACTIVE PAD - Ankle Brace	115
<b>MALLEO</b> ACTIVE BAND - Ankle Brace	116
<b>MALLEO</b> SUPPORT - Ankle Brace	116
<b>MALLEO</b> INTENSE - Ankle Stabilizer	117

## How to order Mobilis

117

# MOBILIS: SOFT ORTHOPEDIC SUPPORTS AND BRACES

A collection of high-quality supports and braces from SIGVARIS GROUP designed to provide support, stability, and pain relief while allowing individuals to stay active.




Mobilis product names are made of two parts:

## 1 Body part where the product is used


	LINE
Elbow	EPI
Hand	MANU
Thumb	RHIZO
Back	LUMBO
Knee	GENU
Ankle	MALLEO

## 2 Support level:




Example 1 + 2 = **MALLEOSUPPORT**

# MODE OF ACTION




## Compression effect

- Decreases swelling and pressure on joints
- Provides pain relief
- Improves mobility



## Massage effect with padding

- Targeted swelling reduction
- Improves blood circulation
- Aids healing



## Joint stabilization

- Rigid elements stabilize the joint and improve muscle activity



## Proprioception

- Aids in perception of movements and body position



## Warming effect


- Increases comfort and may provide pain relief

# HOW TO MEASURE




Measure the arm circumference at 3 cm below the elbow (two fingers)

## Elbow




Measure the elbow circumference with 30 degrees of inclination

## Wrist and thumb




Measure wrist circumference

## Back




Measure waist circumference at belly button

## Knee



Measure knee circumference with 30 degrees of inclination



Measure thigh circumference at 15 cm above the knee

## Ankle



Measure ankle circumference above ankle bone



# Elbow



	EPI CARE	EPI BAND	EPI ACTIVE
INDICATIONS			
Arthritis	○		○
Joint fluid and swelling	○		○
Chronic irritation/ Tendomyopathy	○		○
Post trauma, Post surgery	○		○
Tennis/golfer's elbow (Epicondylopathy)	○	○	○
SUPPORT LEVEL			
	 Support level mild	 Support level medium	 Support level medium



# EPI CARE

## Elbow Support




Contains natural rubber latex that may cause allergic reaction

### INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, tennis/golfer's elbow (epicondylopathy)

### ACTION MODE



### FEATURES

- Mild support
- Maximum comfort
- Discreet look
- Can be worn on either elbow
- Optimal adhesion due to silicone band

### COMPOSITION

64% Polyamide, 27% Elastodiene, 5% Elastane, 4% Silicone



Size	Circumference	SKU
Size XS	20-23 cm	MECE-XS-66
Size S	23-26 cm	MECE-SM-66
Size M	26-29 cm	MECE-MD-66
Size L	29-32 cm	MECE-LG-66
Size XL	32-35 cm	MECE-XL-66

# EPI BAND

## Elbow Strap



### INDICATIONS

Tennis/golfer's elbow (epicondylopathy)

### ACTION MODE




### FEATURES

- Adjustable strap to modify support level
- Silicone pad provides relief for tendon
- Can be worn on either elbow

### COMPOSITION

89% Polyamide, 11% Elastane



Size	Circumference	SKU
One size	20 - 35 cm	MEBE-91

# EPIACTIVE

## Elbow Brace



Support level  
medium



Contains natural  
rubber latex that  
may cause allergic  
reaction


### INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, tennis/golfer's elbow (epicondylopathy)

### ACTION MODE



Compression



Massaging effect



Proprioception

### FEATURES

- Elbow brace with silicone inserts
- Removable elastic strap to adjust support
- Anatomical shape and 3D knitting for maximum comfort
- Can be worn on either elbow

### COMPOSITION

56% Thermoplastic elastomer, 32% Polyamide, 7% Elastodiene, 3% Elastane, 2% Polyester



Measure the elbow  
circumference with  
30 degrees of inclination



Size	Circumference	SKU
Size XS	20 - 23 cm	MEAE-XS-91
Size S	23 - 26 cm	MEAE-SM-91
Size M	26 - 29 cm	MEAE-MD-91
Size L	29 - 32 cm	MEAE-LG-91
Size XL	32 - 35 cm	MEAE-XL-91



Click here to  
watch a video





# Wrist & Thumb



	MANU CARE	MANU WRAP	MANU ACTIVE	MANU SUPPORT	MANU PLUS SUPPORT	RHIZO SUPPORT
INDICATIONS						
Arthritis	○	○	○	○	○	○
Chronic irritation/ Tendomyopathy	○	○	○	○	○	○
Carpal tunnel syndrome			○	○	○	
Instability			○	○	○	○
Post trauma, Post surgery		○	○	○	○	
Rhizarthrosis					○	○
Skier's thumb/Torn thumb ligament (Lateral metacarpo- phalangeal ligament lesion)					○	○
SUPPORT LEVEL	 Support level mild	 Support level moderate	 Support level medium	 Support level strong	 Support level strong	 Support level strong

# MANUCARE

## Wrist Support

### INDICATIONS

Arthritis, chronic irritation/tendomyopathy

### ACTION MODE



Compression



Warming effect



Proprioception

### FEATURES

- Mild support
- Soft compression
- Ultimate comfort
- Discreet look
- Can be worn on either wrist

### COMPOSITION


66% Polyamide, 29% Elastodiene, 5% Elastane



Support level  
**mild**



Contains natural rubber latex that may cause allergic reaction



Measure wrist circumference

Size	Circumference	SKU
Size <b>XS</b>	13 - 15 cm	MMCW-XS-66
Size <b>S</b>	15 - 17 cm	MMCW-SM-66
Size <b>M</b>	17 - 19 cm	MMCW-MD-66
Size <b>L</b>	19 - 21 cm	MMCW-LG-66
Size <b>XL</b>	21 - 23 cm	MMCW-XL-66

# MANUWRAP

## Wrist and Thumb Brace

### INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery

### ACTION MODE



Movement control



Stabilization

### FEATURES


- Moderate adjustable support level
- Breathable fabric
- Can be worn on either hand

### COMPOSITION

96% Polyamide, 4% Elastane



Support level  
**moderate**



Measure wrist circumference

Size	Circumference	SKU
Sizes <b>S/M</b>	13 - 17 cm	MMWW-SM/MD-91
Sizes <b>L/XL</b>	17 - 23 cm	MMWW-LG/XL-91


# MANUACTIVE

## Wrist Brace

### INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery

### ACTION MODE




### FEATURES

- Brace with silicone inserts around cubitus bone
- Metal rod to limit wrist movement
- Elastic strap allows support level adjustment
- Anatomically shaped 3D knit for comfort and fit
- Specific for left or right hand

### COMPOSITION

56% Polyamide, 25% Thermoplastic Elastomer, 9% Elastodiene, 5% Polyester, 5% Elastane



Measure wrist circumference



Contains natural rubber latex that may cause allergic reaction



[Click here to watch a video](#)



Size	Circumference	SKU Left Hand	SKU Right Hand
Size <b>XS</b>	13 - 15 cm	MMAW-XS-L-91	MMAW-XS-R-91
Size <b>S</b>	15 - 17 cm	MMAW-SM-L-91	MMAW-SM-R-91
Size <b>M</b>	17 - 19 cm	MMAW-MD-L-91	MMAW-MD-R-91
Size <b>L</b>	19 - 21 cm	MMAW-LG-L-91	MMAW-LG-R-91
Size <b>XL</b>	21 - 23 cm	MMAW-XL-L-91	MMAW-XL-R-91


# MANUSUPPORT

## Wrist Brace

### INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery

### ACTION MODE



### FEATURES


- Immobilization/stabilization brace: rigid palm structure, external and internal rods to restrain wrist movements
- Three hook-and-loop straps to adjust the fit
- Can be worn on either wrist

### COMPOSITION

92% Polyamide, 5% Elastane, 3% Cotton



[Click here to watch a video](#)



Measure wrist circumference

Size	Circumference	SKU
Sizes <b>S/M</b>	13 - 17 cm	MMSWS-SM/MD-91
Sizes <b>L/XL</b>	17 - 23 cm	MMSWS-LG/XL-91



# MANUPLUS SUPPORT


## Wrist & Thumb Stabilizer



### INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery, rhizarthrosis, skier's thumb/torn thumb ligament (lateral metacarpophalangeal ligament lesion)

### ACTION MODE




### FEATURES

- Immobilization/stabilization brace for wrist and thumb
- The two hook-and-loop straps can be clamped by laces to adjust the fit
- Specific for right or left hand

### COMPOSITION

50% Polyester, 46% Polyamide, 4% Polyurethane



Measure wrist circumference



Size	Circumference	SKU Left Hand	SKU Right Hand
Size <b>XS</b>	13 - 15 cm	MMPSW-XS-L-91	MMPSW-XS-R-91
Size <b>S</b>	15 - 17 cm	MMPSW-SM-L-91	MMPSW-SM-R-91
Size <b>M</b>	17 - 19 cm	MMPSW-MD-L-91	MMPSW-MD-R-91
Size <b>L</b>	19 - 21 cm	MMPSW-LG-L-91	MMPSW-LG-R-91
Size <b>XL</b>	21 - 23 cm	MMPSW-XL-L-91	MMPSW-XL-R-91

# RHIZOSUPPORT


## Thumb Brace



### INDICATIONS

Arthritis, chronic irritation/tendomyopathy, instability, rhizarthrosis, skier's thumb/torn thumb ligament (lateral metacarpophalangeal ligament lesion)

### ACTION MODE



### FEATURES


- Immobilization/stabilization brace for thumb due to external and internal rods
- Two hook-and-loop straps to adjust the fit
- Can be worn on either hand

### COMPOSITION

92% Polyamide, 8% Elastane



Click here to watch a video



Measure wrist circumference




Size	Circumference	SKU
Sizes <b>S/M</b>	13 - 17 cm	MRST-SM/MD-91
Sizes <b>L/XL</b>	17 - 23 cm	MRST-LG/XL-91

**mobilis**  
by SIGVARIS GROUP



# Back



	LUMBO CARE	LUMBO SUPPORT	LUMBO SUPPORT EVO	LUMBO INTENSE
INDICATIONS				
Muscular tension	○	○	○	○
Lower back pain (Lumbar syndrome)	○	○	○	○
Muscular weakness and loss (degeneration)		○	○	○
Post surgery		○	○	○
Sciatica		○	○	○
SUPPORT LEVEL				
	 Support level mild	 Support level strong	 Support level strong	 Support level extra strong

# LUMBOCARE


## Back Support



### INDICATIONS

Muscular tension, lower back pain (lumbar syndrome)

### ACTION MODE




### FEATURES

- Adjustable belt provides mild support
- Comfortable material that provides a warming effect
- Can be worn against the skin

### COMPOSITION

65% Cotton, 31% Polyester, 4% Elastodiene



Size	Circumference	SKU
Size <b>S</b>	70 - 95 cm	MLCB-SM-66
Size <b>M</b>	95 - 115 cm	MLCB-MD-66
Size <b>L</b>	115 - 135 cm	MLCB-LG-66

# LUMBOSUPPORT


## Lumbar Brace



### INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

### ACTION MODE



### FEATURES


- Four stabilizing rods to support and secure
- Adjustable straps for a precise fit
- To be worn over clothing or underwear
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

### COMPOSITION

54% Polyamide, 21% Elastane, 20% Polypropylene, 5% Polyester



[Click here to watch a video](#)



Size	Circumference	SKU Short (S)	SKU Tall (T)
Size <b>S</b>	70 - 95 cm	MLSB-SM-S-91	MLSB-SM-T-91
Size <b>M</b>	95 - 115 cm	MLSB-MD-S-91	MLSB-MD-T-91
Size <b>L</b>	115 - 135 cm	MLSB-LG-S-91	MLSB-LG-T-91



# LUMBO SUPPORT EVO


## Lumbar Brace



### INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

### ACTION MODE



[Click here to watch a video](#)

### FEATURES


- Lumbar spine stabilizing belt with four posterior rods
- The double elastic hook-and-loop harness is removable and permits individual adjustment
- To be worn over the clothes
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

### COMPOSITION

44% Polyamide, 23% Polyester, 8% Viscose, 7% Elastodiene, 6% Ethylene vinyl acetate, 5% Polyvinyl chloride, 4% Polyurethane, 3% Elastane



Contains natural rubber latex that may cause allergic reaction



Size	Circumference	SKU Short (S)	SKU Tall (T)
Size <b>S</b>	56 - 71 cm	MLSEB-SM-S-99	MLSEB-SM-T-99
Size <b>M</b>	68 - 84 cm	MLSEB-MD-S-99	MLSEB-MD-T-99
Size <b>L</b>	81 - 96 cm	MLSEB-LG-S-99	MLSEB-LG-T-99
Size <b>XL</b>	94 - 109 cm	MLSEB-XL-S-99	MLSEB-XL-T-99
Size <b>XXL</b>	107 - 122 cm	MLSEB-2X-S-99	MLSEB-2X-T-99

# LUMBO INTENSE


## Lumbar Brace



### INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

### ACTION MODE



[Click here to watch a video](#)

### FEATURES


- Stabilizing belt with four posterior rods provides stronger pressure on the lumbar spine
- The double elastic hook-and-loop harness is removable and permits individual adjustment
- To be worn over the clothes
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

### COMPOSITION

47% Polyamide, 27% Polyester, 2% Viscose, 4% Elastodiene, 6% Ethylene-acetate, 5% Polyvinyl chloride, 4% Polyurethane, 2% Polypropylene, 3% Elastane








Contains natural rubber latex that may cause allergic reaction



Size	Circumference	SKU Short (S)	SKU Tall (T)
Size <b>S</b>	56 - 71 cm	MLIB-SM-S-99	MLIB-SM-T-99
Size <b>M</b>	68 - 84 cm	MLIB-MD-S-99	MLIB-MD-T-99
Size <b>L</b>	81 - 96 cm	MLIB-LG-S-99	MLIB-LG-T-99
Size <b>XL</b>	94 - 109 cm	MLIB-XL-S-99	MLIB-XL-T-99
Size <b>XXL</b>	107 - 122 cm	MLIB-2X-S-99	MLIB-2X-T-99

# Knee



	GENU CARE	GENU WRAP	GENU BAND	GENU ACTIVE STAB	GENU ACTIVE PAD
<b>INDICATIONS</b>					
Arthritis	○	○		○	○
Joint fluid and swelling	○			○	○
Chronic irritation/ Tendomyopathy	○	○	○	○	○
Post trauma, Post surgery	○	○		○	○
Kneecap (patellar) pain		○	○	○	○
Deep knee pain (Femoro-patellar syndrome)		○	○	○	○
Instability, Hyperlaxity, Meniscus lesion, Gonarthrosis				○	
<b>SUPPORT LEVEL</b>	 Support level <b>mild</b>	 Support level <b>moderate</b>	 Support level <b>medium</b>	 Support level <b>medium</b>	 Support level <b>medium</b>



# GENUCARE

## Knee Support



Contains natural rubber latex that may cause allergic reaction

### INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery

### ACTION MODE



Compression



Proprioception



Warming effect

### FEATURES

- Mild support
- Optimal adhesion due to silicone band
- Discreet look
- Can be worn on either knee

### COMPOSITION

64% Polyamide, 28% Elastodiene  
5% Elastane, 3% Silicone



Size	Circumference	SKU
Size <b>XS</b>	30 - 34 cm	MGCK-XS-66
Size <b>S</b>	34 - 38 cm	MGCK-SM-66
Size <b>M</b>	38 - 42 cm	MGCK-MD-66
Size <b>L</b>	42 - 46 cm	MGCK-LG-66
Size <b>XL</b>	46 - 50 cm	MGCK-XL-66



# GENUWRAP

## Knee Wrap



### INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

### ACTION MODE



Stabilization



Massaging effect




Proprioception

### FEATURES

- Moderate support knee wrap with three straps for optimum fit
- Adjustable size
- Breathe-O-Prene material allows skin to breathe
- Can be worn on either knee

### COMPOSITION

94% Polyamide, 6% Elastane



Size	Circumference	SKU
Size <b>S/M</b>	30 - 42 cm	MGWK-SM/MD-91
Size <b>L/XL</b>	42 - 55 cm	MGWK-LG/XL-91



# GENUBAND

## Knee Strap



### INDICATIONS

chronic irritation/tendomyopathy, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

### ACTION MODE




### FEATURES

- Adjustable knee strap to modify support level
- Silicone insert is located on kneecap tendon
- Can be worn on either knee
- Anatomical shape for extra comfort

### COMPOSITION

89% Polyamide, 11% Elastane



Size	Circumference	SKU
One Size	30 - 55 cm	MGBK-91

# GENUACTIVE STAB


## Stabilizing Knee Brace



### INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome), instability, hyperlaxity, meniscus lesion, gonarthrosis

### ACTION MODE




### FEATURES

- Stabilizing knee brace with lateral supports
- Silicone insert around the kneecap with three Velcro straps
- The anatomically shaped fabric is compressive
- Wearable on both knees

### COMPOSITION

8% Polyamide, 23% Elastomer, 15% Elastodiene, 6% Polyester, 5% Elastane, 2% Polyvinyl chloride, 1% Silicone



Contains natural rubber latex that may cause allergic reaction



Click here to watch a video



Size	Knee Circumference	Thigh Circumference	SKU
Size XS	30 - 34 cm	40 - 44 cm	MGSK-XS-91
Size S	34 - 38 cm	44 - 48 cm	MGSK-SM-91
Size M	38 - 42 cm	48 - 52 cm	MGSK-MD-91
Size L	42 - 46 cm	52 - 56 cm	MGSK-LG-91
Size XL	46 - 50 cm	56 - 60 cm	MGSK-XL-91
Size 2X	50 - 54 cm	60 - 64 cm	MGSK-2X-91

# GENUACTIVE PAD

## Patellar Knee Brace



Support level  
medium



Contains natural  
rubber latex that  
may cause allergic  
reaction

### INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/  
tendomyopathy, post trauma, post surgery, kneecap  
(patellar) pain, deep knee pain (femoro-patellar syndrome)

### ACTION MODE



Compression



Stabilization



Massaging effect




Proprioception

### FEATURES

- Stabilizing knee brace with lateral supports
- Silicone insert around the kneecap
- Anatomically shaped 3D knit for comfort and fit

### COMPOSITION

52% Polyamide, 26% Thermoplastic  
Urethane, 14% Elastodiene,  
5% Polyester, 2% Elastane,  
1% Silicone



Measure knee circumference below the kneecap  
with 30 degrees of inclination, and thigh  
circumference 15 cm above the knee



[Click here to  
watch a video](#)




Size	Knee Circumference	Thigh Circumference	SKU
Size XS	30 - 34 cm	40 - 44 cm	MGAK-XS-91
Size S	34 - 38 cm	44 - 48 cm	MGAK-SM-91
Size M	38 - 42 cm	48 - 52 cm	MGAK-MD-91
Size L	42 - 46 cm	52 - 56 cm	MGAK-LG-91
Size XL	46 - 50 cm	56 - 60 cm	MGAK-XL-91
Size 2X	50 - 54 cm	60 - 64 cm	MGAK-2X-91

**3D Knit**  
Provides breathability and comfort for both  
everyday and physical activities.

**Lateral Supports**  
Provides high stability

**Silicone Insert**  
The anatomical shape helps  
the kneecap to remain stable

**Bend Zone**  
Knitted fabric is less  
compressive for more  
comfort and flexibility







# Ankle

	MALLEO CARE	MALLEO ACTIVE PAD	MALLEO ACTIVE BAND	MALLEO SUPPORT	MALLEO INTENSE
<b>INDICATIONS</b>					
Arthritis	○	○	○	○	○
Joint fluid and swelling	○	○	○		
Chronic irritation/ Tendomyopathy	○	○	○	○	○
Post trauma, Post surgery, Soft tissue lesion	○			○	○
Hyperlaxity				○	○
Chronic instability				○	○
Post trauma, Post surgery, Ligament injury				○	○
Moderate sprains				○	○
<b>SUPPORT LEVEL</b>	 Support level mild	 Support level medium	 Support level medium	 Support level strong	 Support level extra strong



# MALLEOCARE

## Ankle Support

### INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, soft tissue lesion

### ACTION MODE



Compression



Stabilization



Warming effect



Proprioception

### FEATURES

- Mild support
- Adjustable level of stability
- Can be worn on either ankle

### COMPOSITION

70% Polyamide, 12% Elastane,  
9% Elastodiene, 8% Polyester,  
1% Polypropylene



Support level  
**mild**



Contains natural  
rubber latex that  
may cause allergic  
reaction



Measure ankle circumference  
above ankle bone

Size	Circumference	SKU
Size <b>S</b>	18 - 22 cm	MMCA-SM-66
Size <b>M</b>	22 - 25 cm	MMCA-MD-66
Size <b>L</b>	25 - 28 cm	MMCA-LG-66

# MALLEOACTIVE PAD

## Ankle Brace

### INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy

### ACTION MODE



Compression



Stabilization



Massaging effect



Proprioception

### FEATURES

- 2 lateral pads stabilize and support the ankle
- Gentle massaging effect
- Breathable 3D knit for optimal fit and comfort
- Left and right foot option
- Includes a foot slip for easier donning

### COMPOSITION

45% Polyamide, 40% Thermoplastic  
Urethane, 10% Elastodiene, 3% Elastane,  
2% Polyester



Support level  
**medium**



Contains natural  
rubber latex that  
may cause allergic  
reaction



**Click here to  
watch a video**



Measure ankle circumference  
above ankle bone

Size	Circumference	SKU Left Ankle	SKU Right Ankle
Size <b>S</b>	16 - 21 cm	MMAA-SM-L-91	MMAA-SM-R-91
Size <b>M</b>	21 - 26 cm	MMAA-MD-L-91	MMAA-MD-R-91
Size <b>L</b>	26 - 31 cm	MMAA-LG-L-91	MMAA-LG-R-91
Size <b>XL</b>	31 - 37 cm	MMAA-XL-L-91	MMAA-XL-R-91

# MALLEOACTIVE BAND


## Ankle Brace



### INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy

### ACTION MODE



### FEATURES

- The elastic band allows for adjustable support
- The anatomically shaped breathable knit fabric is compressive
- Can be worn in a shoe
- Can be worn on either ankle

### COMPOSITION

62% Polyamide, 18% Elastodiene, 12% Viscose, 5% Polyester, 3% Elastane



Contains natural rubber latex that may cause allergic reaction



Click here to watch a video



Size	Circumference	SKU
Size XS	18 - 20 cm	MMBA-XS-91
Size S	20 - 23 cm	MMBA-SM-91
Size M	23 - 25 cm	MMBA-MD-91
Size L	25 - 27 cm	MMBA-LG-91
Size XL	27 - 30 cm	MMBA-XL-91

# MALLEOSUPPORT


## Ankle Brace



### INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, soft tissue lesion, hyperlaxity, chronic instability, ligament injury, moderate sprains

### ACTION MODE



### FEATURES

- Stabilizing ankle brace with lateral rigid rods
- Two hook-and-loop straps for added support
- The lace closure system allows for an adjustable fit
- Can be worn on either ankle

### COMPOSITION

91% Polyamide, 5% Elastane, 2% Polyester, 2% Cotton



Size	Circumference	SKU
Size S	19 - 22 cm	MMSA-SM-91
Size M	22 - 26 cm	MMSA-MD-91
Size L	26 - 30 cm	MMSA-LG-91



# MALLEOINTENSE

## Ankle Stabilizer



Contains natural rubber latex that may cause allergic reaction




Click here to watch a video



### INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, soft tissue lesion, hyperlaxity, chronic instability, ligament injury, moderate to severe sprains

### ACTION MODE




### FEATURES

- Immobilizing ankle stabilizer based on two rigid shells
- Adaptable fit with two gel cushions to provide comfort and cooling effect
- Two Velcro straps as well as heel support strap provide ability to adjust tightness and support
- Can be worn on either ankle

### COMPOSITION

50% Gel, 22% Polyvinyl chloride, 20% Polypropylene, 7% Polyamide, 1% Neoprene



Size	Height	SKU
One Size	≥ 150 cm	MMIA-91

## HOW TO ORDER MOBILIS PRODUCTS

M	E	A	E	-	MD	-		-	91
Mobilis	Letter matches first part of the product name:  E = Epi M = Manu MP = Manu Plus R = Rhizo L = Lumbo G = Genu M = Malleo	Letter matches line of product:  C = Care W = Wrap A = Active S = Support or Stab B = Band I = Intense SE = Support Evo	Letter matches body part where the product is used:  E = Elbow W = Wrist T = Thumb B = Back K = Knee A = Ankle		Size		When your product has a left or right side, <b>L</b> or <b>R</b> is indicated.  When your product has a short or tall height, <b>S</b> or <b>T</b> is indicated.		Available colours:  <b>66</b> = Beige <b>91</b> = Graphite <b>99</b> = Black

Add **L** to indicate the Left side or **R** to indicate the Right side after the Size in the ordering code. **L** or **R** is not specified in products that fit both left and right. **Example:** MMAW-MD-L-91.

For back braces, add **S** to indicate the Short height or **T** to indicate the Tall height after the Size in the ordering code.

You can find products without a size. **Example:** MEBE-91

You can find products with combined sizes. **Example:** MRST-SM/MD-91

We are dedicated  
to **helping people feel  
their best.** Every day.



**Get more done faster:**

Chat with us live on the  
Online Dealer Portal: [sigvaris-online.com](https://sigvaris-online.com)



Click to see  
digital catalog

**To order:**

[sigvaris-online.com](https://sigvaris-online.com)  
[canadaorders@sigvaris.com](mailto:canadaorders@sigvaris.com)  
1.800.363.4999