

Name		Compounds tested in Healstone Products			
		1	2	3	4
JWH-018	JWH-018 original				
JWH-018	JWH-018 N-propanoic acid; AM2201 N-propanoic acid; JWH 073 N-propanoic acid metabolite	25	25	25	5
JWH-018	5-pentanoic acid metabolite; N-Pentanoic acid metabolite	100	100	50	10
JWH-018	N-4-hydroxypentyl	2000	2000	<2000	500
JWH-019	6-hydroxyhexyl	5000	<2000	<2000	500
JWH-019	5-hydroxyhexyl		<2000	<2000	500
JWH-073	4-butanoic acid metabolite; N-butanoic acid	100	50	50	10
JWH-073	N-2-hydroxybutyl	2000	<2000	<2000	500
JWH-073	N-4-hydroxybutyl	500	500	500	100
JWH-122	N-4-hydroxypentyl	5000	5000	2000	250
JWH-122	N-5-hydroxypentyl	5000	5000	2000	1000
JWH-200	6-hydroxyindole	2000	<2000	2000	250
JWH-210	JWH 210 N-pentanoic acid; JWH 210 N-(5-carboxypentyl) metabolite	<10000	<10000	200	50
JWH-398	N-pentanoic acid metabolite	200		200	50
AM2201	N-(4-hydroxypentyl) metabolite	200	100	>100	>100
MAM2201	N-pentanoic acid metabolite	100	100	100	20
RCS-4	N-(5-carboxypentyl) metabolite	5000		<10000	<10000
JWH-081	JWH 081 N-(4-hydroxypentyl) metabolite	2000		2000	500
JWH-250	JWH 250 N-pentanoic acid metabolite; JWH 250 N-(5-carboxypentyl) metabolite	>10000	>10000	1000	500