

Essential Reagents for Your Lab



Essential tools for every lab, including:

- *Protease Inhibitors*
- *Kinase Modulators*
- *Cell Cycle Modulators*
- *DNA Damaging Agents*
- *Transcription/Translation Inhibitors*

And more

For applications in

- *Cell Signaling*
- *Ubiquitin/Proteasome*
- *Neuroscience*
- *Cancer*
- *Immunology*
- *Epigenetics*



ONS BIO LLC
34194 Aurora Rd, #110
Solon, OH 44139 USA
Phone: 440.482.5005
Email: info@onsbio.com
Web: <https://onsbio.com>

Cat. #	Name	Activity
MC-1008	IBMX	Phosphodiesterase inhibitor; raises cellular cAMP and cGMP levels
MC-2073	Forskolin	Adenylate cyclase activator; raises cellular cAMP and cGMP levels
MC-2079	Ionomycin	Calcium ionophore; raises cytosolic Ca ²⁺
MC-1348	AEBSF	Protease inhibitor
MC-1346	Leupeptin	Trypsin-like and cysteine protease inhibitor
MC-1309	MG-132	Proteasome inhibitor
MC-2614	Z-VAD-FMK	Pan-caspase inhibitor
MC-1061	U-73122	Phospholipase C inhibitor
MC-2091	Okadaic acid	Protein phosphatase inhibitor
MC-2104	Staurosporine	Pan-kinase inhibitor; apoptosis inducer
MC-1104	Rapamycin	mTOR inhibitor; antiproliferative
MC-1182	FCCP	Mitochondrial oxidative phosphorylation inhibitor; apoptosis inducer
MC-2071	Cytochalasin D	Actin filament disruptor
MC-2387	Nocodazole	Microtubule disruptor
MC-2095	Taxol	Microtubule stabilizer
MC-1117	Cycloheximide	Protein biosynthesis inhibitor
MC-2100	Puromycin HCl	Protein biosynthesis inhibitor
MC-1134	Dexamethasone	Glucocorticoid; Transcriptional regulator;
MC-2165	PMA	PKC activator; Transcriptional regulator; Tumor promoter
MC-2054	Actinomycin D	Transcription inhibitor
MC-1123	Etoposide	Topoisomerase II inhibitor; DNA damaging agent
MC-2130	Doxorubicin HCl	DNA damaging agent
MC-1170	Mitomycin C	DNA damaging agent
MC-2794	5-EdU	DNA labeling agent

Please inquire for bulk or multi-vial discounts.

All products are of the highest quality and are in stock for immediate delivery.

Don't see your favorite biochemical tool here? Ask us.