



HEMO

PORTER®

SAMPLE TRANSPORT WITH AN EGDE.

HEMO PORTER® KEEPS SAMPLES COOL.

With the all new **Hemo Porter**® you can now eliminate the need of thermacole and ice pack to carry your samples from patient to lab.

The **Hemo Porter**® can be used to keep your samples cool as they are being transported from one place to another. Thanks to the inbuilt solid-solid phase change system.




HEMO PORTER®
**TEMPERATURE
CHECK.**

Receiving samples from patients at a higher temperature will affect test results.

Now monitor the temperature at which you receive the samples with the innovative **Hemo Porter**®

Extend your sample collection radius with temperature indicator.





HEMO PORTER®

SECURE YOUR SAMPLES.

With **Hemo Porter**® you can ensure complete safety of your samples while being transported.

When transporting the standard 13x75mm blood collection tube, the samples are held tight together in place with their caps.

Any jerks or fall will not let the tubes fall out or move out of the **Hemo Porter**®

HEMO PORTER®

VERTICAL INTEGRATION.

Hemo Porter® comes with ribs on the base as well as lid. They match neatly enabling vertical stacking of the porters both inside as well as outside the fridge.

The ribs on the base also ensure minimum contact with the surface hence reducing conductive heat transfer.





READY TO SHIP.

Hemo porter[®] is shipped in packs of two enough for a single collection person.

The box consists of QR code which can be scanned by any user to read how to use instructions.

The label on the base of the Hemo Porter also has a peelable label with instructions mentioned for use.