


WHAT ARE NOVATHREADS?


They are **hypodermic needles** preloaded with a **PDO absorbable suture.**

PDO ABSORBABLE SUTURE

It's a synthetic absorbable surgical suture composed of polydioxanone. It's one of the safest materials to be implanted in the body. It is fully absorbed by hydrolysis within 4 to 6 months and doesn't create any scar tissue. Being a monofilament, PDO has less chances of harboring bacteria thus reducing chances of infections.



P-DIOXANONE



POLYDIOXANONE (PDO)

Ph: +1 (570) 550-0700

Fax: +1 (305) 503-9602

www.novathreads.com


1444 Biscayne Blvd. Ste 116

Miami, FL 33132

Introducing NOVATHREADS



PRODUCT DESCRIPTION



NOVATHREADS BARBS


A thread is only as strong as its weakest point.

Novathreads barbs are designed to be more resilient than the competition and are arranged in a specific pattern to preserve maximum tensile strength.

**RESORBABLE
& EASY**

No more complications!





SHARP TIPS

NOVATHREADS NEEDLES

Novathreads needles have ultra thin walls that allow a thicker thread to be preloaded rather than the industry standard.

SMOOTH SMOOTH THREAD BAG OF 25PCS

31G x 1" SUTURE: 1¼" · 7-0 USP

29G x 1½" SUTURE: 2" · 6-0 USP

27G x 2½" SUTURE: 3" · 5-0 USP

25G x 3½" SUTURE: 5.5" · 4-0 USP



TWIST TWIST THREAD BAG OF 25PCS

29G x 1" SUTURE: 1¼" · 6-0 USP

29G x 1½" SUTURE: 2" · 6-0 USP

27G x 1½" SUTURE: 2" · 5-0 USP



BARB 2 BARB2 THREAD BAG OF 12PCS

21G x 3½" SUTURE: 6" · 2-0 USP



BARBS IN CANNULAS

Novathreads Barb1s and Barb3s are preloaded cannulas. The blunt tip of the cannula allows it to glide through the tissues preserving nerves and vessels.

BARB 1 BARB1 THREAD BAG OF 12PCS

23G x 2½" SUTURE: 3" · 4-0 USP

21G x 3½" SUTURE: 5.5" · 3-0 USP



BLUNT TIPS



BARB 3 BARB3 THREAD BAG OF 12PCS

23G x 2½" SUTURE: 3.5" · 4-0 USP

21G x 3½" SUTURE: 6.5" · 3-0 USP

19G x 4" SUTURE: 7" · 1-0 USP

