



As of July 1, 2019 all newly purchased or newly purposed vaccine storage units must be single stand-alone units.

The guide below will assist in determining the proper size needed for the estimated doses being stored in the unit at any single time.

	Volume at any	Maximum	Minimum Cubic	Type of Unit Needed
	one time	Doses	Feet	
	Very High	10,000+	Multiple units needed	Full size stand-alone unit
	High	2,000-10,000	May need multiple	
			units	Proces
Refrigerator	Medium High	1,000-2,000	40 cu. ft.	
	Medium	901-1,000	36 cu. ft.	
		801-900	21-23 cu. ft.	
		701-800	17-19.5 cu. ft.	
		400-700	16.7 cu. ft.	
	Low	100-399	4.9-6.7 cu. ft.	Under counter units
	Medium-High	501-6,000	7-14.8 cu. ft.	Full-size units
F	Low	201-500	5-5.6 cu. ft.	
Freezer		0-200	3.5-4.9 cu. ft.	Under counter units
				ABS overar to the to day. Mod Free 2201