



Description

Degreaser gel for sensitive skin preparation for EEG, EP, ECG, Stress Test, Holter, Telemetry recordings. Guarantees low impedance. Dermatologically tested for the most sensitive skins.

- **DYE-FREE**
- Single patient
- Pre-set quantity for single use
- 100% recyclable packaging

UMDNS 11425

Composition

Ingredients
Water, Aluminium Oxide, Propylene Glycol, Sodium Polyacrylate, Phenoxyethanol, Ethylhexylglycerin

Safety

Biocompatibility
(UNI EN ISO 10993) Yes

Patch test ✓ Non irritant

Chemical and physical properties

pH 5.5 – 6.5

Viscosity 60.000-80.000 cP

Specifications

Generic code NPREP010E-1

Package 1 box of 100 units (10 g each)

Environment

RoHS II ✓

REACH ✓