



Injectable Pump It™

The Genesis Instruments Injectable Pump It[™] is a low cost vaccination instrument with a fixed dosage volume ranging from 0 ml to 2 ml. Dosage / actuation is fixed with accuracy of less than 3% error. The unique 45-degree design balances weight for an accurate injection. This instrument is ideal when using 60 cc or 120 cc plastic serum bottles, as well as 100 cc glass serum bottles.



Oral Pump It™

The Genesis Instruments Oral Pump It[™] is an accurate and rapid dispensing instrument used for delivery of oral and topical medicants to companion and food producing animals. Capacity of the unit is from 0 ml to 10 ml per actuation on a fixed dose basis. This instrument affixes onto any threaded cylinder, Boston round or other container neck configuration whether plastic, glass or aluminum.

PRODUCT DESCRIPTION	2 ml Oral Pump It™	2 ml Oral Pump It™
2 ml Oral Pump It™, 1" Trigger & 28 mm Neck	1" Trigger &	1.45" Trigger &
2 ml Oral Pump It™, 1" Trigger & 33 mm Neck		28 mm Neck
2 ml Oral Pump It™, 1.45″ Trigger & 38 mm Neck		m
2 ml Oral Pump It™, 4″ Trigger & 28 mm Neck		
2 ml Oral Pump It™, 4″ Trigger & 38 mm Neck	5 ml Oral Pump It™	
3 ml Oral Pump It™, 1.45″ Trigger & 38 mm Neck	4" Trigger &	
3 ml Oral Pump It™, 4″ Trigger & 28 mm Neck	40 mm Neck	2 ml Oral Pump It™
5 ml Oral Pump It™, 1″ Trigger & 38 mm Neck		4'' Trigger &
5 ml Oral Pump It™, 1.45″ Trigger & 38 mm Neck		28 mm Neck 🔰 🔜 🥼
5 ml Oral Pump It™, 1.45″ Trigger & 40 mm Neck		
5 ml Oral Pump It™, 4″ Trigger & 33 mm Neck		
5 ml Oral Pump It™, 4" Trigger & 40 mm Neck	5 ml Oral Pump It™ 1.45″ Trigger & 38 mm Neck	
10 ml Oral Pump It™, 1″ Trigger & 28 mm Neck		
10 ml Oral Pump It™, 1.45″ Trigger & 28 mm Neck		10 ml Oral Pump It™ 4″ Trigger &
10 ml Oral Pump It™, 1.45″ Trigger & 38 mm Neck		
10 ml Oral Pump It™, 4″ Trigger & 33 mm Neck		33 mm Neck 🔰 🔤 🔤
10 ml Oral Pump It™, 4″ Trigger & 45 mm Neck		Attra L

Modular Pump It™

The Genesis Instruments Modular Pump It[™] is an accurate and rapid dispensing instrument used to deliver liquid nutritionals or vaccines by means of injection, oral, intranasal or topical administration. Manufacturers can assemble their specific pump method using the multiple component system. A variety of threaded or draw-off configurations and custom colors are available.

PRODUCT DESCRIPTION		
1 ml Draw-Off Modular Pump It™ with Luer Lock Trigger		
1 ml Draw-Off Modular Pump It™ with 1.45″ Trigger		
2 ml Draw-Off Modular Pump It™ with Luer Lock Trigger		
2 ml Draw-Off Modular Pump It™ with 4″ Trigger		
1 ml Threaded Modular Pump It™ with Luer Lock Trigger		
1 ml Threaded Modular Pump It™ with 1" Trigger		
1 ml Threaded Modular Pump It™ with 4″ Trigger		
2 ml Threaded Modular Pump It™ with Luer Lock Trigger		
2 ml Threaded Modular Pump It™ 1″ Trigger		
2 ml Threaded Modular Pump It™ with 1.45" Trigger		
2 ml Threaded Modular Pump It™ 4″ Trigger		

