



# MEDICAL CONTROL DIRECTIVE

## 2023-22

**DATE:** November 27, 2023

**TO:** Pinellas County EMS Agencies  
Pinellas County Emergency Communications  
Pinellas County Certified EMTs and Paramedics  
Pinellas County Certified Advanced Practice Paramedics, Nurses  
Pinellas County Online Medical Control Physicians  
Pinellas County Ambulance Billing and Financial Services  
ED Nurse Managers

**FROM:** Dr. Angus Jameson, EMS Medical Director



**RE:** December 2023 In-Service Training

**DISTRIBUTION/TRAINING PERIOD:** November 29 - December 31, 2023

**EFFECTIVE DATE:** December 1, 2023

- An “In-service Training” module separate from regularly scheduled CME has been developed and is assigned to all system clinicians through Vector Solutions for December 2023.
- The in-service includes:
  - Video message from Dr. Jameson
  - System updates
  - Substance Use Disorder (SUD) Survey
- All system clinicians are required to complete this training no later than December 31, 2023.
- Failure to complete this in-service by December 31, 2023, may result in clinical restriction or suspension.

**Attachments:**

- CS22.7 ALS Medical Response Bag
- CS22.15 Vehicle Supplemental Equipment and Supplies

**Distribution:**

- EMSChiefs e-mail distribution group
- Vector Solutions
- Pinellas County EMS Office of the Medical Director Webpage [www.pcemsomd.com](http://www.pcemsomd.com)

# CS22.7 ALS MEDICAL RESPONSE BAG

(This protocol reflects medical supplies, equipment and medications required in compliance with 64J-01 F.A.C.)

Bag			
StatPack G3 or G3+ Perfusion - Blue			
Top Exterior Pocket - Center - Glucometer Kit			
Item Name	PKG/UOM	Qty Rgd	Specific Notes
Glucometer (Bayer Contour)		1	
Glucometer Test Strips	BTL	1	Strips <b>MUST</b> be kept in original bottle - <b>DO NOT</b> mix strips from different bottles - <b>MUST</b> retain bottom of external packaging for initial and monthly quality control testing info
Lancets		10	Single patient use
Bandage, Adhesive, 1" (Band-Aid®)		10	
Alcohol Prep Pads		10	
Top Exterior Pocket - Interior Left Net			
Narcan Nasal Spray Device, 4 mg	Nasal Spray	2	
Top Exterior Pocket - Interior Right Net			
Naloxone, 1 mg/mL, 2 mL	PFS	3	
Storage Box, Opaque Plastic, Two-Part		1	Holds 3 - Naloxone Prefilled Syringes
Top Exterior Pocket - Lid Zippered Pocket			
Dextrose 10% in Water, 250 mL	Pre-mixed Bag	1	
IV Administration Set, 20 gtt (macro)		1	
Glucose Gel, 15g		2	
Storage Box, Clear Plastic, Two-Part		1	Holds the Glucose Gel Packs
Glucagon (Glucagen)		1	<b>One for One Replacement - Requires Incident Number of Use, Expired Item for Replacement, or Incident Report for all other situations</b>
Left Exterior Pocket - Center			
IV Start Kit		3	
Left Exterior Pocket - Interior Left Net			
IV Administration Set, 20 gtt (macro)		1	
Tourniquet (loose) - IV start		3	ORANGE - Non-trauma style tourniquet
Gauze, 4" x 4"	2 per pack	10	Sterile
Tape, Silk, 1"	RL	1	Single patient use
Tape, Self-Adherent, 1"	RL	1	Single patient use - color may vary
Left Exterior Pocket - Interior Right Net			
IV Catheter, 16g		2	Safety
IV Catheter, 18g		4	Safety
IV Catheter, 20g		4	Safety
IV Catheter, 22g		4	Safety
Gauze, Roll, 4"		1	Sterile
Left Exterior Pocket - Zippered Pocket			
0.9% Sodium Chloride, 1000 mL	Pre-mixed Bag	1	
0.9% Sodium Chloride, 10 mL	PFS	3	
Right Exterior Pocket - Center			
EZIO Driver, w/ Trigger Guard		1	Replace <b>ONLY</b> per manufacturer's instructions - <b>Daily Testing is NOT</b> recommended
Right Exterior Pocket - Interior Left Net			
IV Administration Set, 20 gtt (macro)		1	
EZIO Needle Set, 45 mm		2	<b>One for One Replacement - Requires Incident Number of Use, Expired Item for Replacement, or Incident Report for all other situations</b>
EZIO Stabilizer		1	

CS22.7 ALS MEDICAL RESPONSE BAG

# CS22.7 ALS MEDICAL RESPONSE BAG

CS22.7 ALS MEDICAL RESPONSE BAG

Right Exterior Pocket - Interior Right Net			
Item Name	PKG /UOM	Qty Rgd	Specific Notes
EZIO Needle Set, 25 mm		2	One for One Replacement - Requires Incident Number of Use, Expired Item for Replacement, or Incident Report for all other situations
EZIO Stabilizer		1	
Right Exterior Pocket - Zippered Pocket			
0.9% Sodium Chloride, 1000 mL	Pre-mixed Bag	1	
Pressure Infusion Bag, 1000 mL		1	
0.9% Sodium Chloride, 10 mL	PFS	3	
Top Center Interior Pocket			
Ketorolac, 30 mg/mL, 1 mL	PFS	2	
25g x 1" Hypodermic Needle		2	Safety Needle without syringe (Retractable)
Storage Box, Clear Plastic, Two-Part		1	Holds the Ketorolac and Hypodermic Needles
Acetaminophen, 10 mg/mL, 100 mL	Pre-mixed Bag	1	
Controlled Substance Box - Complete	Seahorse SE120 BLACK - Reference separate inventory		
Controlled Substance Box			
Seahorse 120 Black with Cyberlock (CL-C5N)			
Controlled Substance Content Shield (PCEMS)		1	
Etomidate, 2 mg/mL, 20 mL	Vial	2	
Midazolam, 5 mg/mL, 1 mL	Vial or PFS	4	
Fentanyl 50 mcg/mL, 1 mL, or 2 mL	2 mL Vial or 1 mL PFS	4	Max total amount per box - Vial format - 400 mcg or PFS format - 300 mcg
Top Center Interior Pocket - Lid Zippered Pocket			
Mucosal Atomization Device (MAD)		2	Utilized for administration of Midazolam and Fentanyl intranasally <b>ONLY</b>
Lower Center Interior Pocket - Upper Level			
Calcium Chloride, 100 mg/mL, 10 mL	PFS	2	
Atropine Sulfate, 0.1 mg/mL, 10 mL	PFS	2	
Sodium Bicarbonate, 1 mEq/mL, 50 mL	PFS	2	If unavailable reference SODIUM BICARBONATE KIT
Epinephrine, 0.1 mg/mL, 10 mL	PFS	6	If unavailable reference EPINEPHRINE 1 MG/ML - 1 ML VIAL KIT
Lidocaine, 20 mg/mL, 5 mL	PFS	2	If unavailable reference LIDOCAINE VIAL KIT
Adenosine Kit #1*	PFS	1	*1 - 6 mg/2 mL prefilled syringe or vial & 3 Way Stopcock
Adenosine Kit #2**	PFS	1	**2 - 6 mg/2 mL prefilled syringes or vials
Storage Box, Opaque Plastic, Two-Part		2	Used for Adenosine Kit #1 and Adenosine Kit #2
Storage Box, Opaque Plastic, Two-Part		12	Used for the protection of each medication in the prefilled syringe packaging format (PFS)
Lower Center Interior Pocket - Lower Level			
Medication Kit	Flambeau Multi-Compartment - See separate inventory		
Syringe Kit	Flambeau Multi-Compartment - See separate inventory		
Infusion Kit	Flambeau Open Core - See separate inventory		
Sodium Bicarbonate Kit ( <b>ONLY</b> if prefilled syringe format unavailable)	Flambeau Open Core - See separate inventory		

# CS22.7 ALS MEDICAL RESPONSE BAG

Lower Center Interior Pocket - Lid Zippered Pocket			
Item Name	PKG/UOM	Qty Rgd	Specific Notes
Shears, Trauma		1	
Sharps Container, Individual		2	Single patient use
Biohazard Waste Bag, Small <b>(RED)</b>		2	
Biohazard Waste Bag, Large <b>(RED)</b>		1	

Sodium Bicarbonate Kit (Flambeau 6734TE (T4000) Open Core Box)			
Sodium Bicarbonate 1 mEq/mL - 50 mL	Vial	2	
Syringe, 60 mL, Luer-Lock Tip		2	
Syringe, 20 mL, Luer-Lock Tip		2	
Syringe, 10 mL, Luer-Lock Tip		2	
Needle, 18g x 1.5" Blunt Fill with Filter		2	For drawing medications from vials <b>ONLY</b>

Epinephrine Vial Kit (Each kit in storage box, clear plastic, two-part)			
Epinephrine, 1 mg/mL, 1 mL	Vial	1	
0.9% Sodium Chloride, 10 mL	PFS	1	
Needle, 18g x 1.5" Blunt Fill with Filter		1	For drawing medications from vials <b>ONLY</b>

Lidocaine Vial Kit (Each kit in storage box, clear plastic, two-part)			
Lidocaine, 20 mg/mL, 5 mL	Vial	1	
Needle, 18g x 1.5" Blunt Fill with Filter		1	For drawing medications from vials <b>ONLY</b>
Syringe, 10 mL, Luer-Lock Tip		1	

Medication Kit (Flambeau 6747TE (T4007) Multi-Compartment Box)			
Ondansetron, 4 mg	Unit Dose ODT	2	
Ondansetron, 2 mg/mL, 2 mL	PFS	2	
Diphenhydramine, 50 mg/mL, 1 mL	PFS or vial	2	
Epinephrine, 1 mg/mL, 1 mL	Vial	2	
Amiodarone, 50 mg/mL, 3 mL	Vial	3	
Methylprednisolone Sodium Succinate, 125 mg/2 mL	Vial	2	
Nitroglycerin Aerosol Spray, 0.4 mg/spray	BTL	1	Replace per manufacturer's instructions
Spoon, Aspirin administration		3	Individually wrapped
Aspirin, Baby, 81 mg	BTL	1	
Ipratropium Bromide, 0.5 mg/2.5 mL	Unit Dose	2	
Albuterol Sulfate, 2.5 mg/3 mL	Unit Dose	4	
Diltiazem, 5 mg/mL, 5 mL	Vial	1	<i>Date deployed</i> _____ <b>MAX OF 30 DAYS OUT OF REFRIGERATION</b>
Norepinephrine, 1 mg/mL, 4 mL	Vial	1	

**CS22.7 ALS MEDICAL RESPONSE BAG**

# CS22.7 ALS MEDICAL RESPONSE BAG

Syringe Kit (Flambeau 6747TE (T4007) Multi-Compartment Box)			
Item Name	PKG/UOM	Qty Rgd	Specific Notes
Syringe, 1 mL with 25g x 1" Needle or 25g x 1" Hypodermic Needle		3	Safety syringe (Retractable) with needle or Safety Needle without syringe (Retractable)
Syringe, 3 mL with 25g x 1 1/2" Needle or 25g x 1 1/2" Hypodermic Needle		3	Safety syringe (Retractable) with needle or Safety Needle without syringe (Retractable)
Syringe, 3 mL, Luer-Lock Tip		2	
Syringe, 1 mL, Luer-Lock Tip		2	
Alcohol Prep Pads		10	
3-Way Stopcock		2	
Needle, 18 g x 1.5" Blunt Fill with Filter		2	For drawing medications from vials <b>ONLY</b>

Infusion Kit (Flambeau 6734TE (T4000) Open Core Box)			
Medication "ADD" Label, "PCEMS" Specific		4	
IV Administration Set with Flow Controller (20 gtt/mL)		1	
Dextrose 5% in Water, 100 mL	Pre-mixed Bag	1	
Magnesium Sulfate, 2 g/50 mL	Pre-mixed Bag	2	

CS22.7 ALS MEDICAL RESPONSE BAG

# CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

(This protocol reflects medical supplies, equipment and medications required in compliance with 64J-1 F.A.C.)

Equipment & Medical Supplies - <b>Patient Care Action Area</b>								
Item Name	PKG/ UOM	Ambulance			Fire			Specific Notes
		ALS	BLS - 911	BLS - VAN	ALS Transport Capable Rescue	ALS Medic Unit, Squad, Truck, Pumper, or Engine	BLS Engine, Squad, Truck, Pumper, Utility	
Pulse Oximeter, Fingertip		-	1	1	-	-	-	With lanyard, rubber boot, soft case, and Pelican 1010 case
Pelican 1010 Case		-	1	1	-	-	-	Storage of Fingertip Pulse Oximeter
Stethoscope, Adult/Pediatric		1	1	1	1	-	-	Sprague style
Blood Pressure Cuff, <b>Manual</b> , Infant		1	1	1	1	-	-	Cuff labeling reflects "CHILD"
Blood Pressure Cuff, <b>Manual</b> , Child		1	1	1	1	-	-	Cuff labeling reflects "SMALL ADULT"
Blood Pressure Cuff, <b>Manual</b> , Adult		1	1	1	1	-	-	
Blood Pressure Cuff, <b>Manual</b> , Large Adult		1	1	1	1	-	-	
Blood Pressure Cuff, <b>Non-invasive</b> , Child		1	1	-	1	-	-	For use with the Stryker LP15
Blood Pressure Cuff, <b>Non-invasive</b> , XLarge Adult		1	1	-	1	-	-	For use with the Stryker LP15
Glucometer (Bayer Contour)	EA	1	1	1	-	-	-	
Glucometer Test Strips	BTL	1	1	1	-	-	-	retain bottom of external packaging for quality control testing
Disinfectant Wipe, Alcohol	PK/30	2	2	2	1	-	-	70% isopropyl alcohol/30% DI Water

Equipment & Medical Supplies - <b>Reserve</b>								
Oxygen Cylinder, "PCEMS", M6		1	1	1	1	1	-	Spare - minimum 1000 psi
Oxygen Cylinder, "D"		1	1	1	1	1	-	Spare - minimum 1000 psi
Oxygen Cylinder, Onboard		1	1	1	1	-	-	Aluminum or Steel "M" - minimum 500 psi
Oxygen Regulator - Onboard Oxygen		1	1	1	1	-	-	
Oxygen Flowmeter with Hose Barb Adapter (Xmas Tree), Onboard Oxygen		2	2	1	2	-	-	min. 2, 4, 6, 8, 10, 15, 20, 25L flow settings and DISS Port

CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

# CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

Equipment & Medical Supplies - Reserve (cont.)								
Item Name	PKG/ UOM	Ambulance			Fire			Specific Notes
		ALS	BLS - 911	BLS - VAN	ALS Transport Capable Rescue	ALS Medic Unit, Squad, Truck, Pumper, or Engine	BLS Engine, Squad, Truck, Pumper, Utility	
Nasal Cannula, Adult		8	4	4	-	-	-	
Mask, Non-Rebreather, Adult		4	2	2	-	-	-	
Mask, Aerosol Mask, Adult		2	-	-	-	-	-	
Mask, Trach, Venturi with Diluters, Adult		2	2	-	-	-	-	
Nebulizer Setup		4	-	-	-	-	-	
King LTS-D Airway, Size 3		1	-	-	-	-	-	
King LTS-D Airway, Size 4		1	-	-	-	-	-	
King LTS-D Airway, Size 5		1	-	-	-	-	-	
Syringe, 60 mL, Luer-Lock Tip		2	-	-	-	-	-	
Tube Holder, Adult		2	-	-	-	-	-	
BVM Resuscitator, Adult or Small Adult		2	1	1	1	1	-	With adult mask and filter
BVM Resuscitator, Pediatric		1	1	1	1	1	-	With child, infant and neonate masks and filter
OPA 80mm, 90mm, 100mm, 110mm	Each Size	2	1	1	-	-	-	
EtCO2 Filterline Set, Adult/Pediatric		2	-	-	1	1	-	
EtCO2 Nasal Cannula, Adult		10	-	-	-	-	-	
CPAP Setup, Adult		1	-	-	1	1	-	
CPAP Setup, Small Adult/Pediatric		1	-	-	-	-	-	
Superset with Mask Elbow Adapter		1	-	-	-	-	-	
Tee Adapter		1	-	-	-	-	-	
Laryngoscope Handle, Medium		1	-	-	-	-	-	Single patient use, sterile & disposable
Suction Canister Set (canister, lid, suction tubing, vacuum tubing)		2	1	1	1	1	-	CHANGE ALL TUBING (SUCTION AND VACUUM) AND THE CANISTER AFTER EACH USE REGARDLESS OF ANY VISIBLE CONTENTS
Laryngoscope Blade, Mac "3"		1	-	-	-	-	-	Single patient use, sterile & disposable



# CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

Equipment & Medical Supplies - Reserve (cont.)								
Item Name	PKG/ UOM	Ambulance			Fire			Specific Notes
		ALS	BLS - 911	BLS - VAN	ALS Transport Capable Rescue	ALS Medic Unit, Squad, Truck, Pumper, or Engine	BLS Engine, Squad, Truck, Pumper, Utility	
Laryngoscope Blade, Miller "4"		1	-	-	-	-	-	Single patient use, sterile & disposable
Laryngoscope Blade, Mac "4"		1	-	-	-	-	-	Single patient use, sterile & disposable
ET tube, 6.0mm		1	-	-	-	-	-	Cuffed with stylet
ET tube, 7.0mm		1	-	-	-	-	-	Cuffed with stylet
ET tube, 7.5mm		1	-	-	-	-	-	Cuffed with stylet
ET tube, 8.0mm		1	-	-	-	-	-	Cuffed with stylet
ET tube, 8.5mm		1	-	-	-	-	-	Cuffed with stylet
Bougie, Pocket		2	-	-	-	-	-	Single use
Cold Pack		3	3	3	-	-	-	
Heat Pack		2	2	2	-	-	-	
1" Band-Aids		10	10	10	-	-	-	
2" Band-Aids		10	10	10	-	-	-	
1" Silk Tape		2	2	2	-	-	-	Single use
3" Silk Tape		2	2	2	-	-	-	Single use
1" Self-Adherent Tape		2	2	2	-	-	-	Single use
4" Roll Gauze, Sterile		2	2	2	-	-	-	Sterile
Multi-Trauma Dressing, 10" x 30"		2	2	2	-	-	-	Sterile
Splint, Moldable Padded Aluminum		2	2	2	2	2	-	
C-Collar, AMBU Perfit Ace		2	2	1	2	2	-	
C-Collar, AMBU Mini Perfit Ace		2	2	1	2	2	-	
IV Administration Set, 20 gtt (macro)		10	-	-	-	-	-	
IV Start Kit		12	-	-	-	-	-	
IV Catheter, 16g		2	-	-	-	-	-	
IV Catheter, 18g		10	-	-	-	-	-	
IV Catheter, 20g		12	-	-	-	-	-	
IV Catheter, 22g		4	-	-	-	-	-	
IV Administration Set with Flow Controller (20 gtt/mL		1	-	-	1	-	-	
25g x 1" Needle		3	-	-	-	-	-	Safety (Retractable) with or without syringe
25g x 1 ½" Needle		3	-	-	-	-	-	Safety (Retractable) with or without syringe

CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

# CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

Equipment & Medical Supplies - Reserve (cont.)								
Item Name	PKG/ UOM	Ambulance			Fire			
		ALS	BLS - 911	BLS - VAN	ALS Transport Capable Rescue	ALS Medic Unit, Squad, Truck, Pumper, or Engine	BLS Engine, Squad, Truck, Pumper, Utility	
Syringe, 20 mL, Luer-Lock Tip		2	-	-	-	-	-	
Syringe, 10 mL, Luer-Lock Tip		2	-	-	-	-	-	
Syringe, 3 mL, Luer-Lock Tip		5	-	-	-	-	-	
Syringe, 1 mL, Luer-Lock Tip		5	-	-	-	-	-	
3-way Stopcock		2	-	-	-	-	-	
Needle, 18g x 1.5” Blunt Fill with Filter		8	-	-	-	-	-	For drawing medications from vials <b>ONLY</b>
Naloxone, 1 mg/mL, 2 mL	PFS	6	-	-	-	-	-	
Narcan Nasal Spray Device, 4 mg	Nasal Spray	0	4	-	-	-	-	
Mucosal Atomization Device (MAD)		2	-	-	-	-	-	
Ketorolac Kit (3 - 30 mg/mL - 1 mL)	PFS	1	-	-	-	-	-	PFS
Acetaminophen 10 mg/mL - 100 mL	Pre- mixed Bag	2	-	-	-	-	-	
Dextrose 10% in Water 250 mL	Pre- mixed Bag	2	-	-	-	-	-	
0.9% Sodium Chloride, 1000 mL	Pre- mixed Bag	7	-	-	-	-	-	
0.9% Sodium Chloride, 10 mL	PFS	10	-	-	-	-	-	
Sodium Bicarbonate 1 mEq/mL 50 mL	PFS or Vial	2	-	-	-	-	-	
Epinephrine 0.1 mg/mL 10 mL <b>OR</b> Epinephrine 1 mg/mL - 1 mL vial kit if PFS unavailable	PFS or Vial Kit	5	-	-	-	-	-	
Ondansetron 4 mg	Unit Dose ODT	2	-	-	-	-	-	
Ondansetron 2 mg/mL - 2 mL	PFS	2	-	-	-	-	-	
Diphenhydramine 50 mg/mL - 1 mL	PFS or Vial	2	-	-	-	-	-	
Epinephrine 1 mg/mL - 1 mL	Vial	2	-	-	-	-	-	

# CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

Equipment & Medical Supplies - Reserve (cont.)								
Item Name	PKG/ UOM	Ambulance			Fire			Specific Notes
		ALS	BLS - 911	BLS - VAN	ALS Transport Capable Rescue	ALS Medic Unit, Squad, Truck, Pumper, or Engine	BLS Engine, Squad, Truck, Pumper, Utility	
Adenosine 3 mg/mL - 2 mL	PFS or Vial	2	-	-	-	-	-	
Storage Box, Clear Plastic, Two-Part		1	-	-				Used for protection of Adenosine PFS
Methylprednisolone Sodium Succinate 125 mg/2 mL	Vial	1	-	-	-	-	-	
Nitroglycerin Aerosol Spray 0.4 mg/spray	BTL	1	-	-	-	-	-	
Aspirin, Baby, 81 mg	BTL	1	-	-	-	-	-	Chewable tablet - unit dose
Spoon - Aspirin Administration		6	-	-				Individually Wrapped
Ipratropium Bromide 0.5 mg/2.5 mL	Unit Dose	2	-	-	-	-	-	
Albuterol Sulfate 2.5 mg/3 mL	Unit Dose	4	-	-	-	-	-	
Diltiazem 5 mg/mL - 5 mL	Vial	1	-	-	-	-	-	Good for 30 days out of refrigeration
Norepinephrine 1 mg/mL - 4 mL	Vial	1	-	-	-	-	-	
<b>Pelican 1015 Case</b>		1	-	Pending	-	-	-	For protection of the Nasal Naloxone
ECG Monitoring Electrodes		*	-	-	-	-	-	*100 total electrodes - packaging may vary
Alcohol Prep Pads	Ea	50	10	10	-	-	-	
Blood Specimen Draw Kit	Bx	2	-	-	2	2	-	
OB Birthing Kit	Ea	1	1	1	1	1	-	
Head Immobilizer	Set	2	1	1	1	1	-	
Patient Mover		2	2	1	1	1	-	
Restraint, Disposable	Pair	2	2	-	2	2	-	Single patient use
Restraint, Reusable, Poly Style, <b>Wrist</b>	Pair	2	-	-	2	-	-	<b>Reusable - NOT AN EXCHANGE ITEM</b>
Restraint, Reusable, Poly Style, <b>Ankle</b>	Pair	2	-	-	2	-	-	<b>Reusable - NOT AN EXCHANGE ITEM</b>

CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

# CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

Equipment & Medical Supplies - Reserve (cont.)								
Item Name	PKG/ UOM	Ambulance			Fire			Specific Notes
		ALS	BLS - 911	BLS - VAN	ALS Transport Capable Rescue	ALS Medic Unit, Squad, Truck, Pumper, or Engine	BLS Engine, Squad, Truck, Pumper, Utility	
Restraint Belt, Reusable, Poly Style	Individual	4	-	-	4	-	-	<b>Reusable - NOT AN EXCHANGE ITEM</b> - Used with wrist and ankle restraints
Restraint, Reusable, Protective Liner - <b>Wrist</b>	Pair	3	-	-	3	-	-	For use with Poly Style Restraints - Liner is single patient use
Restraint, Reusable, Protective Liner - <b>Ankle</b>	Pair	3	-	-	3	-	-	For use with Poly Style Restraints - Liner is single patient use
Triage Tag - FL Specific Version	Pack	1	1	1	1	1	-	Rev. 5/12 (50 tags/pack)
Triage Ribbon Dispenser System <b>(Fire ONLY!!)</b>		-	-	-	2	2	-	Complete with tape - green, red, yellow, black, magenta
Tamper Evident Security Bag		5	5	5	-	-	-	For securing patient valuables/medications
Patient Belongings Bag		5	5	5	-	-	-	Ambulance ONLY Specific Item
Bed Pan		2	2	2	2	-	-	
Urinal		2	2	2	2	-	-	
Infectious Linen Bag <b>(YELLOW)</b>		3	3	3	3	3	-	
Biohazard Waste Plastic Bag, Small <b>(RED)</b>		4	4	4	-	-	-	
Biohazard Waste Plastic Bag, Large <b>(RED)</b>		4	4	4	-	-	-	
Biohazard Waste Bag Impervious Container		1	1	1	1	1	-	
Sharps Container, Individual		2	2	2	3	3	-	Single Use
Sharps Disposal Container, Wall Mount with Key		1	1	1	1	1	-	Vehicle
Hand Sanitizing Wipe, Alcohol, Individual		50	50	50	25	25	25	Single use

# CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

Equipment & Medical Supplies - Reserve (cont.)								
Item Name	PKG/ UOM	Ambulance			Fire			Specific Notes
		ALS	BLS - 911	BLS - VAN	ALS Transport Capable Rescue	ALS Medic Unit, Squad, Truck, Pumper, or Engine	BLS Engine, Squad, Truck, Pumper, Utility	
Disinfectant Wipe, Alcohol	PK/30	2	2	2	2	2	2	70% isopropyl alcohol/30% DI Water
Hydrogen Peroxide Cleaner - Disinfectant	Spray Bottle	1	1	1	1	1	-	
Wipe, General Cleaning	BX	1	1	1	1	1	-	Disposable, single use
Nitrile Gloves, Non-sterile	PR	Multiple Pairs						Appropriate size
Primary Stretcher with appropriate patient belts per manufacturer	EA	1	1	1	1	-	-	<b>(NOT AN EXCHANGE ITEM)</b>
Sheet, Stretcher		10	10	10	5	-	-	
Pillow		2	2	2	2	-	-	Single use, disposable
Pillowcase		10	10	10	5	-	-	Single use, disposable
Blanket - Cot Quilt		1	1	1	-	-	-	<i>(Sunstar ONLY - for warmth)</i>
Blanket - Cotton for Warmth		4	4	4	4	4	-	Single Use, disposable
Blanket - <b>Yellow</b>		2	2	2	2	2	-	Single Use, disposable - Patient Rain Cover
Pedi-Mate PLUS Pediatric Restraint Device		1	1	1	1	-	-	<b>(NOT AN EXCHANGE ITEM)</b>
NeoMate Pediatric Restraint Device		1	1	1	1	-	-	<b>(NOT AN EXCHANGE ITEM)</b>
Vacuum Splint		1	1	1	1	1	-	Complete with RED Bag
Long Spine Board with Four Straps		2	1	1	1	1	-	
Scoop Stretcher		1	1	1	1	-	-	<b>(NOT AN EXCHANGE ITEM)</b>
Stair Chair		1	1	1	-	-	-	<b>(NOT AN EXCHANGE ITEM)</b>
Patient Slider		2	1	1	-	-	-	<b>(NOT AN EXCHANGE ITEM)</b>
Sager Splint		1	1	1	1	1	-	
Child Car Seat		1	1	1	1	-	-	<b>Check Expiration Date</b> <b>(NOT AN EXCHANGE ITEM)</b>
Cactus Pharmlock Controlled Substance Waste System		1	1	-	1	-	-	Includes bracket set complete and incorporated Cyberlock - <b>MUST</b> be mounted securely in the vehicle

CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

# CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

CS22.15 VEHICLE SUPPLEMENTAL EQUIPMENT & MEDICAL SUPPLIES

Equipment & Medical Supplies - Reserve (cont.)								
Item Name	PKG/ UOM	Ambulance			Fire			Specific Notes
		ALS	BLS - 911	BLS - VAN	ALS Transport Capable Rescue	ALS Medic Unit, Squad, Truck, Pumper, or Engine	BLS Engine, Squad, Truck, Pumper, Utility	
Cactus Pharmlock Controlled Substance Waste System Cartridge <b>(secured in the bracket)</b>		1	1	Pending	1	-	-	Holds a maximum of 3 liters of fluid - <b>MUST be replaced 90 days from the date the cartridge is unsealed for use</b>