



Indelible Voter Ink

QUALITY VOTER INK STAIN FOR ELECTIONS

Liquid ink with semi-permanent dye properties containing colour additives water and/or alcohol and silver nitrate. Used in elections to show a person has voted and to prevent people from voting more than once.

EMS supplier of choice for quality voter ink stain.

Best available prices.

Fastest production and shipping times.

Quality control for all shipments.

Custom volumes, bottle sizes & colour and SN concentrations.



Samples available upon request. For availability, quotes, prices, delivery and shipping times contact sales@ems-sge.com



Colours

Standard base colours of voter ink are violet-blue, red and green. All turn dark black when applied and exposed to sunlight.

Ink Specs

Concentrations 10% to 25% SN

Bottle volumes 5ml to 100ml standard
80ml ink in 100ml bottle

Die cut sponge inserts, heat sealed plastic aluminum sheet

80ml will mark from 500 up to 800 depending method of application and skin conditions,

Standard packaging for shipping, 100 bottles per carton, custom labels, individual packaging in separate plastic bags available

Custom pallet sizes, shrink wrapped, hard corners, custom labels and pallets.

Ink info

Voter ink is applied by two methods: dipping forefinger into a bottle with solution or applied with a felt tipped pen on the finger nail/cuticle.

Safe Over last plus 40 years it is estimated over 500 million voters have dipped their fingers in Voter Ink Stain.

Indelible stain will remain on skin for 72 hours and on the fingernail for up to three to four weeks.

Silver nitrate AgNO_3 is water soluble. When applied to skin it reacts to salt and forms silver chloride which clings to the skin. It is light sensitive and turns dark black when exposed to the sunlight.