

# CNM Testing UOL COVID-19 Test



## Loop-de-Loop™ RT-LAMP Test

The UOL COVID-19 Test is an inexpensive, user-friendly molecular assay using Loop-de-Loop™ RT-LAMP to amplify and detect SARS CoV-2 RNA with PCR-quality sensitivity within 40 minutes.

The UOL COVID-19 Test has been authorized by FDA under an EUA for use at sites operating under a CLIA Waiver, Certificate of Compliance, or Certificate of Accreditation.

## Test Overview


The UOL COVID-19 Test allows for convenient and streamlined testing via a reusable instrument and disposable, single-use test kits, with exportable results displayed on the companion app.



## Key Features

- + Positive results in as fast as 12 minutes<sup>1</sup> (40 minutes for negative results)
- 🎯 Limit of Detection: 400 cp/swab<sup>1</sup> (PCR: 1,800cp/mL<sup>2</sup> – 18,000cp/mL<sup>3</sup>)
- 📋 96.1% Accuracy compared to high-sensitivity PCR<sup>1</sup> (100% NPA)
- 👉 Easy sample collection via anterior nasal swab
- ☁️ Mobile app provides easy to interpret, exportable results

## Test Workflow



### NOTES:

1. UOL COVID-19 Test EUA (<https://www.fda.gov/media/156102/download>)
2. Roche cobas SARS-CoV-2 PCR assay EUA (<https://www.fda.gov/media/136049/download>)
3. CDC 2019-nCoV RT-PCR Diagnostic Panel EUA (<https://www.fda.gov/media/134922/download>)

For more information, contact us at [info@cnmtesting.com](mailto:info@cnmtesting.com)

[cnmtesting.com](http://cnmtesting.com) [LinkedIn](#) [f](#) [Instagram](#)