

# OSMOTOUCH 40

THE NEXT-GEN OSMOMETER

**Rapid and accurate analysis of osmolality in fluid samples.**

The OsmoTouch 40 from Astori Tecnica is a dedicated solution for the rapid analysis of the freezing point in biological fluid, pharmaceutical, and industrial solutions.



**FREEZING  
POINT**



**INCLUDED**




**OPTIONAL**

## Exciting new features:

- Real-time graph of the sample temperature
- Alphanumeric sample ID can be added easily
- Optional barcode reader for a quick sample identification
- 21 CFR Part 11 compliant software
- On-board memory for up to 4,000 results
- Up to 3 calibration ranges
- USB port for data export
- CryoSoft Touch software for PC monitoring & data management

## TECHNICAL SPECIFICATIONS

<b>Sample Capacity</b>	<b>Up to 40 samples</b>
<b>Reproducibility</b>	<b>±2.5 mOsm/Kg (0-600 mOsm/Kg) ±0.5% (600-3000 mOsm/Kg)</b>
<b>Resolution</b>	<b>1 mOsm/Kg</b>
<b>Sample volume</b>	<b>50 - 200 µL</b>
<b>Calibration</b>	<b>Automatic calibration</b>
<b>Display</b>	<b>Full colour 7" touchscreen</b>
<b>Printer/PC</b>	<b>2 ports for connecting an optional printer and PC</b>
<b>Barcode reader</b>	<b>Optional</b>
<b>Dimensions</b>	<b>13x25x15 in (wxdxh, head down)</b>
<b>Net weight</b>	<b>46 lbs</b>
<b>Operating Temperature</b>	<b>+5°C to 36°C</b>



CRYOSOFT TOUCH SOFTWARE

## ACCESSORIES

PART NUMBER	DESCRIPTION	PART NUMBER	DESCRIPTION
MAT-70100	Osmometer Calibration Standard 100 mOsm	AS2030N500	Sample Tubes (500/PK)
MAT-70323	Osmometer Reference Standard 323 mOsm	AS563	Tube Holder (80 position)
MAT-70900	Osmometer Calibration Standard 900 mOsm	AS67200	Cooling Liquid 250 mL
MAT-71500	Osmometer Calibration Standard 1500 mOsm	AS67210	Cooling Liquid 1L
MAT-72000	Osmometer Calibration Standard 2000 mOsm	AS67219	Thermal Printer Paper
MAT-72500	Osmometer Calibration Standard 2500 mOsm	AS08103	Filter, Liquid
MAT-73000	Osmometer Calibration Standard 3000 mOsm		