

Certificate of Analysis

Customer Name:	Confidential on behalf of RS COATINGS LTD	Date Received: 29.06.2020
Customer Contact:	CONFIDENTIAL	Date Analysed: 04.08.2020
Customer Email/Phone:	CONFIDENTIAL	Date Reported: 06.08.2020
Certificate Number:	BL057/2020	

 Tick box if customer has provided untreated control samples for inclusion in testing

Test Method: Determination of Antibacterial Activity using Test Based on MOD ISO 22196: 2011

BioLabTests Sample Reference Number	Sample	Test Bacteria	Contact Time*		Reduction against Initial	
			0 Hrs	24 Hrs	Log ₁₀	%
200/482	BioGrip – Non-Contact Touch Tool with Antimicrobial	MRSA	2.91 x10 ⁴	5.67 x10 ²	1.71	98.05%
		<i>E. coli</i>	1.74 x10 ⁴	≤100	≥2.24	≥99.43%

*Numbers represent Colony Forming Units at representative contact times.

Data show the changes in size of the initial MRSA (ATCC® 43300™) and *E. coli* (ATCC®8739™) populations due to their contact with the treated surface of the samples. The samples detailed in this CoA were tested for 24 hours (±1) at 35°C (±2) under a relative humidity of >90% as specified in the ISO 22196 method. The results detailed in this CoA relate only to the items tested. This certificate shall not be reproduced except in full, without written approval of BioLabTests.

Please note; sample was received in a ready to test state from the customer apart from sterilisation which is noted below.
 All testing is performed on site at the BioLabTests address above unless otherwise disclosed.



The material of the cover slips is polyethylene with a thickness of 0.045mm.

Please note – For all information regarding deviations from ISO 22196:2011, please request statement B-Lab-Stats-007

Tel: +44 (0) 333 240 8308

Fax: +44 (0) 2476 338081

info@biolabtests.com

www.biolabtests.com

Noted during testing	
	Note
No. of samples supplied per test:	6 of 6 supplied
Sample size:	Average size: 50 x 50mm
Sample sterilisation method:	DRY HEAT AT 65°C FOR ONE HOUR
Film size:	Average size: 40 x 40mm
Inoculation volume:	Volume: 0.05ml
Initial bacterial concentration:	MRSA: 5.83×10^5 CFU/ml <i>E. coli</i> : 3.48×10^5 CFU/ml



Ria Warren BSc (Hons), MRes
Quality Manager,
BioLabTests Ltd

END OF REPORT

