

BATRIK FIRST Wipes are Cleanroom wipes destined to do the first wipe-down of an endoscope after a procedure. They are made from hydroentangled non-woven material composed of 45% polyester and 55% cellulose. These wipes are manufactured and packaged in a cleanroom environment.

Features

- Compatible with virtually any solvent or solution for cleaning and wiping (such as **PrAmLi Power³ Enzymes™**)
- Double-bag packaging
- Steam autoclavable
- Production process guarantees extremely low particle generation and chemical extractables

Cleanroom Environment

- Can be used for a cleanroom environment of ISO Class 5 through 8



FIRST Wipes

Item #	Description	Size	Packaging
17-A0909	Nonwoven Cleanroom Wipes	9" x 9" (23 cm x 23 cm)	300 u. / bag 14 bags / case
17-A1212	Nonwoven Cleanroom Wipes	12" x 12" (31 cm x 31 cm)	150 u. / bag 10 bags / case