

## Product Description **QS FLOOR VACUUM CLEANER**

The **QS FLOOR VACUUM CLEANER** helps to quickly and easily remove fluids from the floor of the operating room during an operation. When the suction system is activated, all fluids are removed from the affected area.

Ensures a dry and stable floor!

Reduces cleaning intervals and setup costs!

Comfortable, safe and pain-free everyday work!

Fast, reliable and safe delivery!

With the help of the **QS FLOOR VACUUM CLEANER**, fluids from the floor of the operating room during an operation can be quickly and easily removed! The **QS FLOOR VACUUM CLEANER** is connected to the suction canister of a vacuum aspirator or to the central vacuum station using the pre-assembled hose. When placed on the affected area with the suction system activated, the **QS FLOOR VACUUM CLEANER** absorbs the fluids from the floor. Time-consuming and costly post-operative cleaning and storage of the **QS FLOOR VACUUM CLEANER** is not necessary, as it is a disposable product.

The **QS FLOOR VACUUM CLEANER** is made of pink-coloured transparent plastic, the transparent hose is made of a flexible and bendable material.

**QS FLOOR VACUUM CLEANER** dimensions: Ø 205 mm – Hose length 2,400 mm

Packing : 15 pcs

