



# DALMATIAN INVENTO Multifunction Standing Frame

Represented in Canada by: East West Medical Products Ltd.  
 www.eastwestmedical.ca  
 Email: info@eastwestmedical.ca









To order / Distributed in Canada by: Assistability  
 Unit G 7803 – 35 street SE  
 Calgary AB T2C 1V3  
 http://www.assistability.ca/  
 Email: info@assistability.ca  
 Tel: 1-800-403-9019/ 403-629-7697 / Fax: 403.800.9900







<b>Order #:</b>		<b>Date:</b>	
<b>Name:</b>	<b>Phone #:</b>	<b>Email:</b>	
<b>Bill to:</b>			
<b>Account #:</b>		<b>Company name:</b>	
<b>Street:</b>		<b>Street 2:</b>	
<b>City:</b>	<b>Province:</b>	<b>Postal:</b>	
<b>Ship to:</b>			
<b>Account #:</b>		<b>Company name:</b>	
<b>Street:</b>		<b>Street 2:</b>	
<b>City:</b>	<b>Province:</b>	<b>Postal:</b>	
<b>Name:</b>	<b>Phone #:</b>	<b>Email:</b>	
<b>Notes:</b>			








**STANDARD EQUIPMENT (Circle size and choose options below)**



	<ul style="list-style-type: none"> <li>• 3 in 1: Three utility functions for use in one device: seat, standing frame and lying function</li> <li>• Frame with castors and brakes</li> <li>• Pneumatic spring</li> <li>• Deep adjustable seat</li> <li>• Removable and adjustable abduction block</li> <li>• Knee support</li> <li>• Backrest</li> <li>• Pelvic belt ('hold &amp; pull' system)</li> <li>• Chest support (size 1, 2)</li> <li>• Independently adjustable chest and pelvic support (size 3)</li> <li>• 4 point safety vest ('hold &amp; pull' system)</li> <li>• Adjustable HexagonPro© Headrest</li> <li>• Forearm supports</li> <li>• Tray</li> <li>• Independently adjustable footrest</li> <li>• Innovative quadric mechanism</li> <li>• Breathable upholstery ATMO™</li> </ul> <p>SIZE 1 GRASS,            SIZE 2 SEA,            SIZE 3 LATE.</p>	<b>Size 1: \$5,443</b>
		<b>Size 2: \$5,539</b>
		<b>Size 3: \$5,869</b>

## ACCESSORIES

CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY	
SAFETY	DMI_002 	Belt for moving patients AM-P - Robust and ergonomic design as well as high durability facilitates daily care of the patient.	1	\$312	
	DMI_003 	Soft cover for tray - Soft cover protects the user from being hit by the tray. Indicated for patients with impaired limbs coordination.	1	\$327	
	TANGO Wheel: Noiseless, easy to maintain, on-marking, do not corrode.				
	DMI_015 	TANGO Wheel with brake (75 mm) (I)	1	\$126	
		TANGO Wheel with brake (75 mm) (A)	1	\$186	
	DMI_016 	TANGO Wheel with brake (125 mm) (I)	1	\$178	
		TANGO Wheel with brake (125 mm) (A)	1	\$219	
	DMI_020 	TANGO Wheel with brake (100 mm) (I)	1	\$152	
		TANGO Wheel with brake (100 mm) (A)	1	\$201	
	POSITIONING	DMI_101 	Thigh abduction belt - Can easily be attached to the chair by use of Velcro straps. Keeps legs in abduction and also prevents the patient from slipping off of the device.	1	\$256
DMI_102 		Head supporting belt - Attachable to the headrest by use of Velcro straps. Indicated for people with muscle weakness to prevent falling of the head.	1	\$163	
DMI_103 		3D footplates with sandals (U) (I) - Allow foot setting in all directions with the use of a hinge ball. Can adjust foot setting on footplates (front/back/left/right)	2	\$654	

<p>DMI_106</p> 	<p>Lower limb supporting cushion - Support lower legs</p>	<p>2</p>	<p>\$212</p>
<p>DMI_109</p> 	<p>Independently adjustable side supports(I) - Size 1,2</p>	<p>2</p>	<p>\$396</p>
<p>H Harness acts as a trunk stabilizer. HOLD&amp;PULL© system makes it easy to adjust the vest to the patient's body. You can quickly detach the vest if necessary.</p>			
<p>DMI_114</p> 	<p>H Harness ('hold &amp; pull' system) (I)</p>	<p>1</p>	<p>\$186</p>
<p>H Harness ('hold &amp; pull' system) (A)</p>			<p>\$238</p>
<p>DMI_123</p> 	<p>Adjustable extended headrest (I)</p>	<p>1</p>	<p>\$137</p>
<p>6 points safety vest provides a highly adjustable trunk stabilizer. HOLD&amp;PULL© system makes it easy to adjust the vest to the patient's body.</p>			
<p>DMI_125</p> 	<p>6 points safety vest (I)</p>	<p>1</p>	<p>\$283</p>
<p>6 points safety vest (A)</p>			<p>\$377</p>
<p>DMI_126</p> 	<p>Trunk belt - Universal solution applicable to many of our products. It secures the user from falling out of the device. Anatomical shape and Velcro straps allow for a comfortable fit.</p>	<p>1</p>	<p>\$282</p>

BODYMAP POSITIONING	New version of Hexagon Pro™ BodyMap headrest with Hexagon Pro™ and BodyMap cushion. Such a combination provides a perfect stable fit of the patient's head.			
	DMI_340 	HexagonPro™ BodyMap® Headrest DZ – SZ 1	1	\$1,196
		HexagonPro™ BodyMap® Headrest DZ – SZ 2	1	\$1,286
		HexagonPro™ BodyMap® Headrest DZ – SZ 3	1	\$1,406
COMFORT	DMI_401 	Tray for manual therapy (1) - Allows for active exercises assisting upper limbs with a particular focus on palm function.	1	\$826
	DMI_407 	Tray for nonverbal communication (1) - Set consists of 38 cards that may be put in random order into 24 places of the tray. It is covered with transparent and aesthetic plate made of Plexiglas®. Parents and child can create their own symbols they believe are significant for their communication.	1	\$1,247
	DMI_410 	Headrest cotton cover - Helps keep the headrest clean. Easy to assemble/disassemble. Can be washed.	1	\$30
			2	\$52
	DMI_412 	Elastico seat cushion - Precisely adapts to the patient's body. Made of memory foam covered with breathable and spatial fabric. It is indicated in anti-bedsore prophylaxis.	1	\$256
	DMI_413 	Elastico backrest cushion - Perfectly adapts to the patient's body. Made of dense shape memory foam covered with breathable and spatial fabric. It is indicated in anti-bedsore prophylaxis.	1	\$282
	DMI_419 	Shaped seat cushion (thigh shape)	1	\$175

STORAGE	<p>DMI_507</p> 	Oxygen tank holder	1	Price depends on tank size	
SPARE PARTS	DMI_700*	<p>Waterproof upholstery HYDRO™ (I)</p> <p>- Waterproof upholstery <b>Hydro™</b> is made of soft and comfortable material resistant to water and moisture. This upholstery is very durable and easy to clean and disinfect. This makes it perfect for indoor use when eating and drinking fluids, and outdoors during unfavorable weather conditions such as rain or snow. <b>Hydro™</b> is made of environmental friendly ecological materials.</p>	1	\$260	
	<p>The breathable upholstery Atmo™ is made of elastic, soft material. It is characterized by high durability and is resistant to stretching. Using three-dimensional material allows air to flow, which is crucial in case of pressure sores prevention. Also, it protects against bacteria and mould. Atmo™ is easy to clean with mild cleaning agents. Upholstery is made of environmental friendly materials. Atmo™ is a high level of comfort and safe for the skin.</p>				
	DMI_720*	Breathable upholstery ATMO™ - GRASS, SEA, LATTE (I)	1	\$0	
		Breathable upholstery ATMO™ - RED(I)	1	\$279	
	DMI_725*	<p>Medical upholstery MEDICUS™ (I)</p> <p>- Medical upholstery <b>Medicus™</b> is made of a waterproof material that ensures moisture circulation (vapor permeability). <b>Medicus™</b> is hypoallergenic (safe for allergy sufferers), prevents angers and bedsores. Upholstery is very durable and easy to clean and disinfect. This makes it perfect for using at home during daily activities such as eating and drinking. <b>Medicus™</b> meets the requirements for medical devices so it is recommended to use in health care centers.</p>	1	\$316	
	DMI_726	Standard feet sandal	1	\$234	
	DMI_727	Chest supports	2	\$389	
	<p>DMI_728</p> 	<p>4 points safety vest</p> <p>- Acts as a torso stabilizer. It has a large adjustment range, which facilitates its adjustment. Equipped with the <b>CHUCK AND COMBINE®</b> system.</p>	1	\$245	
DMI_739*	<p>Medical upholstery EXTREME-MED™ (I)</p> <p>- <b>Extreme-Med™</b> is high advanced medical upholstery. It is vinyl fabric with a lot of certificated, e.g. flammability test <b>CRIB7</b> and Moist-Heat Resistant <b>Tropical Test</b>. Furthermore, <b>Extreme-Med™</b> is also abrasion resistant, what is confirmed by <b>Martindale test</b> (500k cycles). The fabric is very durably, flame retardant and can be used in every kind of weather. Upholster is safe for patient: does not contain phthalates, is hypoallergenic, protects against the moulds, bacetrias and dermatophytosis. The fabric is very easy to clean and maintain. Using the special PCV layer, material does not absorb liquids and other stains. <b>Extreme-Med™</b> meets the requirements for medical devices so it is recommended to use in health care centers.</p>	1	\$334		

3 in 1 device:  
chair, standing  
frame and lying



Stander position






Positioning chair



Lying position

**AVAILABLE SIZES**



Symbol	Measurement	Unit	Size 1	Size 2	Size 3
W	Indicative user height*	[cm]	< 100	80 - 130	100 - 165
A	Seat width	[cm]	33	36	45
L	Seat depth	[cm]	16 - 25	20 - 31	27 - 42
C	Headrest height	[cm]	37 - 52	42 - 53	49 - 64
C1	Backrest height	[cm]	28	35	40
E	Backrest width	[cm]	13 - 23	18 - 26	18 - 32
F	Footrests length	[cm]	15 - 34	15 - 37	15 - 45
	Knee supports spacing	[cm]	18 - 26	21 - 29	25 - 37
	Knee supports height	[cm]	19 - 46	28 - 51	19 - 58
	Footplate length	[cm]	19	22	25
L1	Abduction block adjustment (from backrest)	[cm]	0 - 17	5 - 25	7 - 30
	Tray height	[cm]	67 - 80	74 - 97	95 - 112,5
	Stander tilt angle	[°]	-15 - +85	-15 - +85	-15 - +85
X	Backrest tilt angle	[°]	90 - 180	90 - 180	90 - 180
Z	Footrests tilt angle	[°]	90 - 180	90 - 180	90 - 180
	Maximum user weight	[kg]	25	45	60

**Note:** Measurement tolerance range is +/- 1 cm / ½ in.  
\*The value cannot be main factor for selection of product size.