

JUMBO positioning chair

Represented in Canada by: East West Medical Products Ltd. www.eastwestmedical.ca

Email: info@eastwestmedical.ca

To order / Distributed in Canada by: Assistability

Unit G 7803 – 35 street SE Calgary AB T2C 1V3

http://www.assistability.ca/

Email: info@assistability.ca

		Tel: 1-800-403-9	9019/ 403-6	29-7697 / Fax: 403.800.9900	
Order #:		Date:			
Name:	Phone #:		Email:		
	Bill to:				
Account #:		Company name:			
Street:		Street 2:			
City:	Province:		Postal:		
	Ship to:				
Account #:		Company name:			
Street:		Street 2:			
City:	Province:		Postal:		
Name:	Phone #:		Email:		
Notes:					
STANDARD EQI	JIPMENT (Circle size a	and choose options b	elow)		
	 Frame with castors and brakes Footrest with feet stabilizing belts Seat Seat cushion Abduction block Backrest Backrest cushion Independently adjustable chest and pelvic supports 4 points safety vest ('hold & pull' system) Breathable upholstery ATMO™ SIZE 1 GRASS, SIZE 2 SEA, SIZE 3 LATE 		Size 1: \$1,743		
			d pelvic	Size 2: \$1,843	
			Size 3: \$1,941		

		ACCESSORIES			
CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY	
	SLK_002	Belt for moving patients AM-P - Robust and ergonomic design as well as high durability facilitates daily care of the patient.	1	\$312	
SAFETY	SLK_003	Soft cover for tray - Soft cover protects the user from being hit by the tray. Indicated for patients with impaired limbs coordination.	1	\$327	
SAFETT	E	Tango wheels ae noise-free, maintenance-free, non-marking and do not corrode.			
	CLIV 045	TANGO Wheel with brake (75 mm) (I)	4	\$126	
	SLK_015	TANGO Wheel with brake (75 mm) (A)	1	\$186	
		TANGO Wheel with brake (125 mm) (I)		\$178	
	SLK_016	TANGO Wheel with brake (125 mm) (A)	1	\$219	
		TANGO Wheel with brake (100 mm) (I)	_	\$152	
	SLK_020	TANGO Wheel with brake (100 mm) (A)	1	\$201	
	SLK_101	Thigh abduction belts - Can easily be attached to the chair by use of Velcro straps. Helps position legs in abduction and prevents the patient from sliding.	1	\$256	
	SLK_102	Head supporting belt - Attachable to the headrest by use of Velcro straps. Indicated for people with muscle weakness to prevent falling of the head.	1	\$163	
POSITIONING	SLK_104	HexagonPro Headrest™ - Individual "ladyfingers" on the Hexagon Pro™ are adjusted to cushion and support any head position. It has been designed to create the most comfort for user by supporting and protecting their head. Available in many sizes and colours. It is a new version of Hexagon™ headrest equipped with hardware that allows precise adjustment of the headrest.		\$446	
	· ·	ovides trunk stability. HOLD&PULL© system makes it easier to adjust the name the vest if necessary.	vest to the pa	tient body. You	
	SLK_114	H Harness ('hold & pull' system) (I)	1	\$238	
		H Harness ('hold & pull' system) (A)		\$137	
	SLK_123	Adjustable extended headrest (I)	1	\$238	

	6 points safety ve to the patient's bo	st provides a highly adjustable trunk stabilizer. HOLD&PULL© system ma ody.	kes it easy to	adjust the vest
	SLK_125	6 points safety vest – S/M/L/XL (I)	1	\$283
		6 points safety vest – S/M/L/XL (A)	-	\$377
		exagon Pro™ BodyMap headrest with Hexagon Pro™ and BodyMap cushi t stable fit of the patient's head.	on. Such a cor	mbination
	SLK_340	HexagonPro™ BodyMap® Headrest DZ – Sz 1		\$1,196
BODYMAP POSITIONING	Manilynis	HexagonPro™ BodyMap® Headrest DZ – Sz 2	1	\$1,286
		HexagonPro™ BodyMap® Headrest DZ – Sz 3		\$1,406
	SLK_401	Tray for manual therapy (1) - Allows for active exercises assisting upper limbs with a particular focus on palm function.	1	\$826
	SLK_403	Tray - Tray made of ecological wood. Great for therapy, daily activities and meals.		\$396
	SLK_407	Tray for nonverbal communication (1) - Set consists of 38 cards that may be put in random order into 24 places of the tray. It is covered with transparent and aesthetic plate made of Plexiglas®. Parents and child can create their own symbols they believe are significant for their communication.	1	\$1,247
	SLK_410	Headrest cotton cover - Helps keep the headrest clean. Easy to assemble/disassemble. Can be washed.	1	\$30
	J _o .		2	\$52
S	SLK_412	Elastico seat cushion - Precisely adapts to the patient's body. Made of memory foam covered with breathable and spatial fabric. It is indicated in antibedsore prophylaxis.	1	\$256
	SLK_413	Elastico backrest cushion - Perfectly adapts to the patient's body. Made of dense shape memory foam covered with breathable and spatial fabric. It is indicated in anti-bedsore prophylaxis.	1	\$282
	SLK_421	Shaped seat cushion (pommel shape)	1	\$327
	SLK_419	Shaped seat cushion (thigh shape)	1	\$212

	SLK_429	Soft footrest cover - Comfortable pad that protects the feet from direct contact with the footplate		\$230
STORAGE	SLK_507	Oxygen tank holder - Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank.	1	Price depends on tank size
	upholstery is very fluids, and outdoo friendly ecological		use when ea is made of er	ting and drinking nvironmental
	- Medical upholstery Medicus™ is made of a waterproof material that ensures moisture circulation (vapor permeability). Medicus™ is hypoallergenic (safe for allergy sufferers), prevents angers and bedsores. Upholstery is very durable and easy to clean and disinfect. This makes it perfect for using at home during daily activities such as eating and drinking. Medicus™ meets the requirements for medical devices so it is recommended to use in health care centers.			
SPARE PARTS	- Extreme-Med™ is high advanced medical upholstery. It is vinyl fabric with a lot of certificated, e.g. flammability test CRIB7 and Moist-Heat Resistant Tropical Test. Furthermore, Extreme-Med™ is also abrasion resistant, what is confirmed by Martindale test (500k cycles). The fabric is very durably, flame retardant and can be used in every kind of weather. Upholster is safe for patient: does not contain phthalates, is hypoallergenic, protects against the moulds, bacetrias and dermatophytosis. The fabric is very easy to clean and maintain. Using the special PCV layer, material does not absorb liquids and other stains. Extreme-Med™ meets the requirements for medical devices so it is recommended to use in health care centers.			
• · · · · · · · · · · · · · · · · · · ·	SLK_700	Waterproof upholstery HYDRO™(I)	1	\$223
	CLV 720	Breathable upholstery ATMO ™ - GRASS, SEA, LATTE (I)	1	\$0
	SLK_720	Breathable upholstery ATMO ™ - RED(I)	1	\$242
	SLK_725	Medical upholstery MEDICUS™ (I)	1	\$279
	SLK_728	4 points safety vest - Acts as a torso stabilizer. It has a large adjustment range, which facilitates its adjustment. Equipped with the CHUCK AND COMBINE ® system.		\$245
	SLK 739	Medical upholstery EXTREME MED™(I)	1	\$297

COLOURS

More about our upholsteries you can find <u>here</u>

BREATHABLE UPHOLSTERY ATMO™
















Sea Size 2



Latte Size 3



*Optional

WATERPROOF UPHOLSTERY HYDRO™











Green Grass *Optional



Blue Sea *Optional



Brown Sugar *Optional



*Optional

MEDICAL UPHOLSTERY - EXTREME MED™

















Citrus Green *Optional



Candy *Optional



*Optional

MEDICAL UPHOLSTERY - MEDICUS™

















Med Skin *Optional



*Optional

Real upholstery colors may differ from the presented images.

^{*}Optional colors available for special request. Production time ca. 14 days.



Symbol	Measurement	Unit	Size 1	Size 2	Size 3
W	Indicative user height*	[cm]	70 - 130	100 - 150	110 - 180
Α	Seat width	[cm]	35	39	46
L	Seat depth	[cm]	13 - 31	19 - 38	29 - 54
С	headrest height (optional)	[cm]	54	58	75
C1	Backrest height	[cm]	34	38	55
E	Backrest width	[cm]	17 - 27	22 - 32	22 - 35
F	Foot length adjustment	[cm]	18 - 35	24 - 45	28 - 45
L1 Seat depth to the abduction block Footer length	Seat depth to the abduction block	[cm]	0 - 24	3 – 29	8 – 45
	Footer length	[cm]	19	25	26
X	Seat tilt angle	[°]	10	10	10
Z	Backrest tilt angle	[°]	90 – 110	90 – 110	90 – 110
	Max. user weight	[kg]	30	45	75

Note: Measurement tolerance range is +/- 1 cm / $\frac{1}{2}$ in. *The value cannot be main factor for selection of product size