

| | |
|-----------------|--------------|
| Order #: | Date: |
|-----------------|--------------|

| | | |
|--------------|-----------------|---------------|
| Name: | Phone #: | Email: |
|--------------|-----------------|---------------|

Bill to:

| | |
|-------------------|----------------------|
| Account #: | Company name: |
|-------------------|----------------------|

| | |
|----------------|------------------|
| Street: | Street 2: |
|----------------|------------------|

| | | |
|--------------|------------------|----------------|
| City: | Province: | Postal: |
|--------------|------------------|----------------|

Ship to:

| | |
|-------------------|----------------------|
| Account #: | Company name: |
|-------------------|----------------------|

| | |
|----------------|------------------|
| Street: | Street 2: |
|----------------|------------------|

| | | |
|--------------|------------------|----------------|
| City: | Province: | Postal: |
|--------------|------------------|----------------|

| | | |
|--------------|-----------------|---------------|
| Name: | Phone #: | Email: |
|--------------|-----------------|---------------|

Notes:

STANDARD EQUIPMENT (Circle size and choose options below)



- Base
 - Seat
 - Removable abductor
 - Chest belt
 - Adjustable chest support
 - Pelvic belt
 - Adjustable pelvic support
 - Castors with brakes (size 3 only)
 - Breathable upholstery ATMO™
- SIZE 1 GRASS,
SIZE 2 SEA,
SIZE 3 LATE

Size 1: \$1,448


Size 2: \$1,579

ACCESSORIES

| CATEGORY | SYMBOL | EQUIPMENT | PIECES | QTY |
|--|--|---|--------|---------|
| SAFETY | NKK_003  | Soft cover for tray - Soft cover protects the user from being hit by the tray. Indicated for patients with impaired limbs coordination. | 1 | \$327 |
| | NKK_102  | Head supporting belt - Attachable to the headrest by use of Velcro straps. Indicated for people with muscle weakness to prevent falling of the head. | 1 | \$163 |
| POSITIONING | NKK_104  | HexagonPro Headrest™ -Individual “ladyfingers” on the Hexagon Pro™ are adjusted to cushion and support any head position. It has been designed to create the most comfort for user | 1 | \$483 |
| | NKK_123  | Adjustable extended headrest | 1 | \$520 |
| BODYMAP POSITIONING | New version of Hexagon Pro™ BodyMap headrest with Hexagon Pro™ and BodyMap cushion provides comfort and stabilization of the patient’s head. | | | |
| | NKK_340  | HexagonPro© Bodymap Headrest™ - Sz 1 | 1 | \$1,196 |
| | HexagonPro© Bodymap Headrest™ - Sz 2 | \$1,286 | | |
| COMFORT | NKK_403  | Tray - Made of ecological friendly wood. It lets patient play and supports therapy. Ideal for daily activities and meals. | 1 | \$593 |
| | NKK_410  | Headrest cotton cover - Helps to keep the headrest clean and can be washed. Easy to remove/put on. | 1 | \$30 |
| | | | 2 | \$52 |
| NKK_423  | Castors with brakes (75mm) | 5 | \$316 | |

| | | | | |
|-------------|--|--|-------|-------|
| SPARE PARTS | NKK_700 | Waterproof upholstery HYDRO™ (I)* - Waterproof upholstery Hydro™ is made of soft and comfortable material resistant to water and moisture. This upholstery is very durable and easy to clean and disinfect. This makes it perfect for indoor use when eating and drinking fluids, and outdoors during unfavorable weather conditions such as rain or snow. Hydro™ is made of environmental friendly ecological materials. | 1 | \$223 |
| | The breathable upholstery Atmo™ is made of elastic, soft material. It is characterized by high durability and is resistant to stretching. Using three-dimensional material allows air to flow, which is crucial in case of pressure sores prevention. Also, it protects against bacteria and mould. Atmo™ is easy to clean with mild cleaning agents. Upholstery is made of environmental friendly materials. Atmo™ has a high level of comfort and is safe for the skin. | | | |
| | NKK_720 | Breathable upholstery ATMO™ - GRASS, SEA, LATTE (I)* | 1 | \$0 |
| | | Breathable upholstery ATMO™ - RED(I)* | 1 | \$242 |
| | Medical upholstery Medicus™ is made of a waterproof material that ensures moisture circulation (vapor permeability). Medicus™ is hypoallergenic (safe for allergy sufferers) and helps prevent bedsores. Upholstery is very durable and easy to clean and disinfect. This makes it perfect for using at home during daily activities such as eating and drinking. Medicus™ meets the requirements for medical devices so it is recommended to use in health care centers. | | | |
| | NKK_725 | Medical upholstery MEDICUS™ (I)* | 1 | \$260 |
| NKK_739 | Medical upholstery EXTREME MED™ (I)* | 1 | \$279 | |

AVAILABLE SIZES



| Symbol | Measurement | Unit | Size 1 | Size 2 |
|----------|-----------------------------------|-------------|-----------------|-----------------|
| W | Indicative user height* | [cm] | < 120 | < 160 |
| A | Seat width | [cm] | 50 | 56 |
| L | Seat depth | [cm] | 29 | 32 |
| C | Backrest height | [cm] | 33 - 41 | 47 - 58 |
| E | Pelvic supports width | [cm] | 26 | 30 |
| E1 | szerokość peloty piersiowej | [cm] | 26 | 30 |
| L1 | Seat depth to the abduction block | [cm] | 12 - 22 | 14 - 24 |
| | Maximum user weight | [kg] | 25 | 35 |

Note: Measurement tolerance range is +/- 1 cm / ½ in.

*The value cannot be main factor for selection of product size.