AKCESMED Our passion - your independence!

CAT II INVENTO Supine/Prone Standing Frame

Represented in Canada by: East West Medical Products Ltd. www.eastwestmedical.ca

Email: info@eastwestmedical.ca

To order / Distributed in Canada by: Assistability
Unit G 7803 – 35 street SE
Calgary AB T2C 1V3

http://www.assistability.ca/

Size 2: \$4,688

AADL Order Form

			Email: orders@assistability.ca		
		Tel: 1-800	0-403-9019/ 403-629-7697 / Fax: 403.800.9900		
Order #:		Date:			
Name:	Phone #:		Email:		
		Bill to:			
Account #:		Company name:			
Street:		Street 2:	Street 2:		
City:	Province:		Postal:		
		Ship to:			
Account #:		Company name:	Company name:		
Street:	treet:				
City:	Province:		Postal:		
Name:	Phone #:		Email:		
	STANDARD EQUIPMEN	NT (Circle size and option	ns below)		
	 Foot plate w Gas spring Adjustable k Tray Chest and position (adjustable) Chest and position (adjustable) 	castors and brakes ith foot stabilization nee support elvic pads with safety elvic supports (adjusteadrest for supine ports)	stable)		

· Safety lock for changing the frame angle

Innovative quadric mechanism
Breathable upholstery ATMO™


Tilt angle indicator


SIZE 1 GRASSSIZE 2 SEA

		ACCESSORIES				
CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY		
	Tango wheels are noiseless, easy to maintain, and do not leave marks or corrode.					
SAFETY	KTI_015	Tango Wheel (75 mm) (I)	1	\$126		
		Tango Wheel (75 mm) (A)	1	\$186		
	KTI_016	Tango Wheel (125 mm) (I)	1	\$178		
		Tango Wheel (125 mm) (A)		\$219		
	VTI 020	Tango Wheel (100 mm) (I)	1	\$178		
	KTI_020	Tango Wheel (100 mm) (A)	1	\$219		
POSITIONING	KTI_102	Head supporting belt - Attached to the headrest by use of Velcro straps. Indicated for people with muscle weakness to prevent falling of the head.	1	\$163		
	KTI_103	Sandals with 3D Adjustment(U)(I) - Allow for a smooth foot setting in all directions with the use of the hinge ball. You can correct your foot setting on the footplates (front/back, left/right) and also set foot in space (plantar/dorsal flexion, pronation/supination of the foot).	1	\$736		
	Soft knee supports stabilize the lower limbs. Made of soft foam.					
	KTI_110 Soft knee supports (I)			\$124		
	as the second	Soft knee supports (A)	2	\$156		
COMFORT	KTI_403	Tray - Easy to attach and remove. It is also child proof, which means that your child cannot release the tray. Light and easy to clean. Perfect for playing and feeding. Provides extra trunk support as well.	1	\$259		
	KTI_425	Tray with bowl i (I) - allows he storage of toys or other things needed for therapy.		\$455		
SPARE PARTS	durable and eas	olstery Hydro™ is made of soft and comfortable material resistant to water and y to clean and disinfect. This makes it perfect for indoor use when eating and dr other conditions such as rain or snow. Hydro™ is made of environmental friendl	rinking fluids, a	and outdoors dur		
		Waterproof upholstery HYDRO™ (I)*	1	\$242		
	KTI_700	Waterproof upholstery HYDRO™ (A)*	1	\$297		
	stretching. Using protects against	upholstery Atmo™ is made of elastic, soft material. It is characterized by high du three-dimensional material allows air to flow, which is crucial in case of pressu bacteria and mould. Atmo™ is easy to clean with mild cleaning agents. Upholst s. Atmo™ is a high level of comfort and safe for the skin.	ure sores preve	ention. Also, it		

		Breathable upholstery ATMO ™GRASS, SEA (I)*	1	\$0
	KTI_720	Breathable upholstery ATMO ™LATTE, RED (I)*		\$253
		Breathable upholstery ATMO ™GRASS, SEA, LATTE, RED (A)*	1	\$331
	Medical upholstery Medicus™ is made of a waterproof material that ensures moisture circulation (vapor permeability). Medicus™ is hypoallergenic (safe for allergy sufferers), prevents angers and bedsores. Upholstery is very durable and easy to clean and disinfect. This makes it perfect for using at home during daily activities such as eating and drinking. Medicus™ meets the requirements for medical devices so it is recommended to use in health care centers.			
	KTI_725	Medical upholstery MEDICUS™ (I)*		\$346
		Medical upholstery MEDICUS™ (A)*	1	\$457
	Extreme-Med™ is high advanced medical upholstery. It is vinyl fabric with a lot of certificated, e.g. flammability tes Moist-Heat Resistant Tropical Test. Furthermore, Extreme-Med™ is also abrasion resistant, what is confirmed by N (500k cycles). The fabric is very durably, flame retardant and can be used in every kind of weather. Upholster is saf does not contain phthalates, is hypoallergenic, protects against the moulds, bacteria and dermatophytosis. The falt to clean and maintain. Using the special PCV layer, material does not absorb liquids and other stains. Extreme-Med requirements for medical devices so it is recommended to use in health care centers.			by Martindale test is safe for patient: ne fabric is very easy
	VTI 720	Medical upholstery EXTREME-MED ™ (I)*	1	\$353
	KTI_739	Medical upholstery EXTREME-MED ™ (A)*		\$464

AVAILABLE SIZES





Symbol	Measurement	Unit	Size I		Size II	
	Measurement	Onic	Supine	Prone	Supine	Prone
w	Indicative user height*	[cm]	70 - 100 13 - 27		80 - 130 18 - 31	
A/A1	Pelvic/chest supports width	[cm]				
c	Pelvic supports height adjustment	[cm]	30 - 59	34 - 59	27 - 78	40 - 78
C1	Chest supports height adjustment	[cm]	46 - 75	50 - 75	50 - 97	65 - 97
B/B1	Pelvic/chest supports depth	[cm]	14		14	
D	Headrest height	[cm]	72 - 106	=	80 - 135	-
G	Knee supports width	[cm]	14 – 23	14 - 23	14 - 23	14 - 23
Н	Knee supports height	[cm]	15 – 31	15 - 54	15 – 54	15 - 54
	Tilt angle	[°]	0 – 90		0 – 90	
	Max user weight	[kg]	35		45	

Note: Measurement tolerance range is +/- 1 cm / ½ in.
*The value cannot be main factor for selection of product size.